



Jul 22, 2020

© Marmoset NHP Free Floating Anti-GFP Antibody Stain After HCR FISH

Guoping Feng Lab¹

¹Massachusetts Institute of Technology

In Development This protocol is published without a DOI.

F Guoping Feng Lab

ABSTRACT

A protocol for staining NHP tissue for anti-GFP antibody after staining for markers of intrest with HCR FISH.

PROTOCOL CITATION

Guoping Feng Lab 2020. Marmoset NHP Free Floating Anti-GFP Antibody Stain After HCR FISH. **protocols.io**

https://protocols.io/view/marmoset-nhp-free-floating-anti-gfp-antibody-stain-bij9kcr6

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 14, 2020

LAST MODIFIED

Jul 22, 2020

PROTOCOL INTEGER ID

39265

MATERIALS

NAME	CATALOG #	VENDOR
Goat serum	G9023	Sigma Aldrich
GFP Polyclonal Antibody	A-11122	Invitrogen - Thermo Fisher
Goat anti-Rabbit IgG (H L) Highly Cross-Adsorbed Secondary Antibody Alexa Fluor 488	A-11034	Invitrogen - Thermo Fisher

Primary Antibody 3d

- Permeabilized marmoset brain sections (→ 100 μm , fixed previously in [M]4 % volume PFA) in [M]1.2 % volume Triton X-100 in PBS for ⓒ 00:30:00 .
- Block sections with [M]10 % volume NS (natural serum depending on Secondary Antibody) and PBT [M]0.1 % volume Triton X-100, 1X PBS] for © 01:00:00.
- 3 Incubate primary antibody (1:500 whole mount [M].002 % volume) with [M]3 % volume NS and PBT-Azide [M]0.2 % volume NaN3, 1X PBT] at § 4 °C for © 72:00:00.

Secondary Antibody 2d

- 4 Wash sections in 1X PBS for **© 00:10:00** 6X.
- 5 Incubate with secondary antibody (1:500 whole mount [M].002 % volume) in PBT at 8 4 °C for © 24:00:00.
- 6 Wash sections in wash buffer for © 00:10:00 6X.
- 7 Stain for DAPI (1:10000 in PBS [M].0001 % volume).
- 8 Wash sections in 1X PBS for © 00:05:00 2X.