



AUG 29, 2023

## 🌐 Qiagen RNEasy RNA extraction protocol

Michael Dan  
Jessica Pardy<sup>1</sup>, Siemon<sup>2</sup>,  
Justin  
Donovan<sup>1</sup>, Dilan  
Joseph<sup>1</sup>,  
christopher.degroot<sup>1</sup>, Richard Gibson<sup>1</sup>

<sup>1</sup>University of Western Ontario; <sup>2</sup>Western University



Michael Dan Siemon  
Western University

OPEN  ACCESS



### ABSTRACT

This protocol was employed by Western University for RNA extraction of wastewater samples for wastewater-based epimology in London, Ontario, Canada and surrounding area. The protocol was adapted from QIAamp® Viral RNA Mini Handbook for use with the Qiagen RNeasy extraction kit.

### MATERIALS

QIAamp RNeasy viral minikit - Qiagen

#### DOI:

[dx.doi.org/10.17504/protocols.io.dm6gp337dvzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gp337dvzp/v1)

#### External link:

<https://www.qiagen.com/us/resources/download.aspx?id=c80685c0-4103-49ea-aa72-8989420e3018&lang=en>

**Protocol Citation:** Jessica Pardy, Michael Dan Siemon, Dilan Joseph, Justin Donovan, christopher.degroot, Richard Gibson 2023. Qiagen RNEasy RNA extraction protocol.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.dm6gp337dvzp/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

This protocol was used by Western University for RNA extraction of wastewater samples for wastewater-based epidemiology. It was used from 2020 to April 2023.

**Created:** Aug 24, 2023

**Last Modified:** Aug 29, 2023

**PROTOCOL integer ID:**  
86937

## Protocol modified from QIAamp® Viral RNA Mini Handbook

- 1 Thoroughly mix wastewater sample then aliquot 40 mL into 50 mL Falcon tube. Centrifuge at 12000 RPM for 90 min. Decant supernatant, assume 280 µl pellet.
- 2 Pipet 1120 µl prepared Buffer AVL into Falcon tube with sample
- 3 Pulse-vortex samples for 15s every 2 minutes for 8 minutes
- 4 Add 1120ul ethanol then centrifuge at 4000RPM for 5 minutes
- 5 Process supernatant through QiAmp Mini column by adding 750 µl and centrifuging at 8000 rpm for 1 min. Remove flow-through and repeat until all supernatant has been processed
- 6 Add 500 µl Buffer AW1 to QiAmp Mini column and centrifuge at 6000 x g (8000 rpm) for 1 min. Replace the collection tube with a clean 2 ml collection tube

- 7** Add 500 µl Buffer AW2 to QiAmp Mini column and centrifuge at 6000 x g (8000 rpm) for 1 min.  
Replace the collection tube with a clean 2 ml collection tube
- 8** Dry membrane by centrifuging at 13000 rpm for 3 min
- 9** Place the QIAamp Mini column in a clean 1.5 ml eppendorf tube.  
Add 60 µl Buffer AVE equilibrated to room temperature.
- 10** Centrifuge at 6000 x g (8000 rpm) for 1 min