



Jul 21, 2021

# Chemical Fixation and embedding of cultured cells for Transmission Electron Microscopy

Marta Orlando<sup>1</sup><sup>1</sup>Charité Universitätsmedizin Berlin

1 Works for me



Share

This protocol is published without a DOI.

Marta Orlando

## DISCLAIMER

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## ABSTRACT

This protocol can be followed to fix primary cultured neurons (or cells in general) and embed them in epoxy resin. Later ultrathin sections can be obtained with an ultramicrotome for Transmission Electron Microscopy imaging. The reagents used in this protocol are toxic!

## PROTOCOL CITATION

Marta Orlando 2021. Chemical Fixation and embedding of cultured cells for Transmission Electron Microscopy. **protocols.io**  
<https://protocols.io/view/chemical-fixation-and-embedding-of-cultured-cells-bwsbpean>

## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 21, 2021

## LAST MODIFIED

Jul 21, 2021

## PROTOCOL INTEGER ID

51747

## GUIDELINES

work under chemical hood with safety protections

## SAFETY WARNINGS

Please work under a chemical hood with gloves and lab coat for all steps of this protocol!

glutaraldehyde, osmium tetroxide are highly toxic  
also sodium cacodylate buffer contains arsenic and is toxic

uranyl acetate is radioactive

do not inhale/ingest/touch without gloves and read the safety information for each reagent before starting to work.

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

1 

Fix for 1 hr at RT in 1.25% glutaraldehyde in 66mM sodium cacodylate buffer

2 

Wash 3 times with 0.1 M sodium cacodylate buffer

3 

Postfix 1 hr in 1% OsO<sub>4</sub> (optional: for better contrast add 1.5% K<sub>4</sub>Fe(CN)<sub>6</sub> to the osmium solution) in 0.1M sodium cacodylate buffer

4 

Wash 3x10' with 0.1 M sodium cacodylate buffer

5 

Wash 3x10' with distilled H<sub>2</sub>O

6 

En bloc stain for 45' in 1% uranyl acetate solution in distilled H<sub>2</sub>O

7 

Wash 3x10' with distilled H<sub>2</sub>O

8 

Dehydrate by incubating in increasing EtOH concentrations: 10 minutes in 70%, 10 minutes in 80%, 10 minutes in 96%, 2x15 minutes in 100%

9 

Infiltrate in 30% Epoxy resin in EtOH

10 

10. Infiltrate in 70% Epoxy resin in EtOH

11 

11. Infiltrate in 90% Epoxy resin in EtOH

12 

12. Incubate in 100% Epoxy resin

13 

13. Flat embed in Epoxy resin

14 

14. Bake for 48 hrs at 65°C

15 15. Remove the glass coverslip from the resin block by thermal shock