



Oct 13, 2020

Protocol 3: Antibiotic Resistance

Forked from [Protocol 3: Antibiotic Resistance](#)

1

¹UCSC**1** Works for me

This protocol is published without a DOI.

UCSC BME 22L



PROTOCOL CITATION

2020. Protocol 3: Antibiotic Resistance. **protocols.io**
<https://protocols.io/view/protocol-3-antibiotic-resistance-bndbma2n>

FORK FROM

Forked from Protocol 3: Antibiotic Resistance, Alyssa Ayala

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 13, 2020

LAST MODIFIED

Oct 13, 2020

PROTOCOL INTEGER ID

43139

MATERIALS TEXT

Perishables

1 - Sterile Water Tube

Non-perishables

1 - 10-100 µL Variable Volume Adjustable Pipette

1 - Box 1-200 µL Pipette Tips

14 - Petri Plates

1 - Microcentrifuge Tube Rack

5 - Inoculation Loops

1 - Bag of Microcentrifuge Tubes


DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

BEFORE STARTING

The final step is testing the success of the bacterial transformation. There are multiple ways bacterial cells can become antibiotic resistant. As you may recognize in the video, the E Coli in this experiment is made resistant to streptomycin through modifying the antibiotic binding target (rpsL gene) through transformation.



- 1 **Plate the CRISPR solution.** Take LB Strep/Kan/Arabinose agar plate out of the fridge and let sit to adjust to room temperature. Shake the microcentrifuge tube containing the CRISPR experiment and then pour or pipette contents onto the agar plate. Use an inoculation loop to spread contents onto the plate. Let dry for a few minutes. Only let the plate dry in a sanitary space or other contaminants may join the plate.
- 2 **Incubate contents on LB plate.** Put the lid back on, flip the plate over, make sure to label it, and leave the bacteria to culture in room temperature. It will take 1 to 2 days to grow.
- 3 **Congratulations!** You have now completed your first CRISPR experiment! Go you.
- 4 **Take a picture of your end culture** and input it on your Lab Notebook. If you did not see any growth, identify factors that may have contributed to this result.