

Feb 25, 2020

Preparation of leukocytes by differential lysis of erythrocytes

Mesut Karataş¹, Şenol Doğan², Emrulla Spahiu³, Adna Ašić¹, Larisa Bešić¹, Yusuf Turan¹

¹International Burch University, ²Universität Leipzig, ³Middle East Technical University

1

Works for me

[dx.doi.org/10.17504/protocols.io.bcv2iw8e](https://doi.org/10.17504/protocols.io.bcv2iw8e)


Emrulla Spahiu

Middle East Technical University



ABSTRACT

Leukocytes are isolated by centrifugation after specific lysis of erythrocytes

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

<https://doi.org/10.1101/2020.02.23.961599>

MATERIALS TEXT

REAGENTS






1. Lysis Buffer (155 mmol/L ammonium chloride; 10 mmol/L sodium bicarbonate; 0.1 mmol/L EDTA).

Dissolve 8.3 g ammonium chloride, 0.84 g sodium bicarbonate and 29.3 mg EDTA in about 900 mL reagent-grade water. Titrate to pH 7.4 with HCl, then make the volume to 1 litre. Store at 4°C.

2. Isotonic Saline (0.9%, w/v)

Dissolve 9 g NaCl in 1 litre reagent-grade water. Store at 4 degrees Celsius.

- 1 Centrifuge **10 ml** EDTA blood to pellet all cells **1500 rpm, 4°C 00:10:00** 10m
- 2 Remove plasma into a clean container and freeze.
- 3 Restore original blood volume with 0.9% saline and transfer the blood suspension into a **50 ml** conical centrifuge tube.
- 4 Add **40 ml** cold lysis buffer
- 5 Stand **On ice**, mixing occasionally, until erythrocytes are lysed (the red cell suspension remains red in colour, but becomes transparent). This should take only about **00:05:00** to **00:10:00** 10m
- 6 Centrifuge **1500 rpm, 4°C 00:05:00** 5m
- 7 Discard supernatant and resuspend leucocyte pellet in **5 ml** cold lysis buffer.
- 8 Stand **On ice** for **00:10:00** 10m
- 9 Dilute cell suspension to **50 ml** with cold 0.9% saline. Mix and centrifuge **1500 rpm, 4°C 00:05:00** 5m

- 10 Discard supernatant, then resuspend leucocyte pellet in  **10 ml** 0.9% saline. Care should be taken to obtain an even cell suspension without being too vigorous and causing cell disruption.
- 11 Divide the cell suspension into two equal aliquots, into two  **10 ml** conical centrifuge tubes.
- 12 Centrifuge  **1500 rpm, 4°C 00:05:00** 5m
- 13 Remove all supernatant and dry walls of centrifuge tube with a tissue.
- 14 Store leucocyte pellets and plasma at  **-20 °C** (or  **-80 °C**)



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited