



Mar 01, 2022

Immunofluorescent Staining

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5s5q6g6 **Haley Geertsma**
University of Ottawa

This protocol is used to stain cryosectioned mouse brain tissue.

DOI

dx.doi.org/10.17504/protocols.io.b5s5q6g6

Haley Geertsma 2022. Immunofluorescent Staining. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5s5q6g6>



_____ protocol ,

Mar 01, 2022

Mar 01, 2022

58941

- 1 To cryo-sectioned brain tissue, wash with 1X phosphate buffered saline (PBS) for 3x 5-minute^{15m} washes.
- 2 Incubate in blocking buffer for 1 hour at room temperature. 1h
Blocking buffer: 10% serum + 0.5% Triton X-100 in 1X PBS
- 3 Wash tissue with 1X PBS. 5m
- 4 Incubate in primary antibody diluted in blocking buffer overnight at 4°C. 1d

- 5 Wash tissue with 1X PBS for 5x 5-minute washes. 30m
- 6 Incubate in secondary antibody diluted in blocking buffer for 1 hour at room temperature. 1h
- 7 Wash tissue with 1X PBS for 5x 5-minute washes. 30m
- 8 If tissue wasn't previously mounted on a slide, mount on a superfrost plus slide and let dry at room temperature for at least 10 minutes. 15m
- 9 Coverslip with fluorescent mounting medium and a #1.5 coverslip. Outline the coverslip with clear nailpolish. 1m