



JUL 16, 2023

Immunofluorescence staining of myenteric and submucosal plexuses

Connor Monahan¹

¹Columbia University



Connor Monahan

Columbia University

ABSTRACT

This protocol details myenteric and submucosal plexus immunostaining.

ATTACHMENTS

[786-2005.pdf](#)

MATERIALS

Materials



Normal Donkey Serum Jackson ImmunoResearch Laboratories, Inc. Catalog #017-000-121



Rabbit anti-tyrosine hydroxylase antibody; AB_390204 Merck Millipore (EMD Millipore) Catalog #AB152



CD3 antibody | 145-2C11 Bio-Rad Laboratories Catalog #MCA2690



IBA Ab Wako Catalog #019-19741



Recombinant Anti-CD68 antibody [FA-11] Abcam Catalog #ab53444

- 1% Triton-X
- PBS
- DAPI

OPEN ACCESS

DOI:

dx.doi.org/10.17504/protocols.io.q26g7p85kgwz/v1

Protocol Citation: Connor Monahan 2023. Immunofluorescence staining of myenteric and submucosal plexuses. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7p85kgwz/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jul 12, 2023

Last Modified: Jul 16, 2023

PROTOCOL integer ID:
84907

Keywords: myenteric plexus, submucosal plexus, immunostain, immunofluorescence, enteric neuron, ileum, macrophage, T cell


Procedure

5h

1

Note

Note: The plexuses are fragile. At each wash, use a micropipette and a dissecting microscope to carefully remove the solution from the well not suck up the plexus.



Cut out two small regions of each plexus sheet in PBS at  4 °C .

2

Add each tissue section to separate wells in a 96 well dish.

3





Add 10% normal donkey serum (Jackson ImmunoResearch, Cat #017-000-121; West Grove, PA) and 1% Triton-X in PBS. Leave in blocking solution for  01:00:00 at  Room temperature on a rotator.

1h

4



Incubate the tissue  Overnight at  Room temperature on a rotator with primary antibodies in 10% normal donkey serum, 1% Triton-X in PBS.

1h

4.1

Stain sections with antibodies against tyrosine hydroxylase (TH) (1:500, Millipore-Sigma, Cat # AB152; Burlington, MA) and CD3 (1:400, Bio-Rad Laboratories, Cat #MCA2690; Hercules, CA), ANNA1 (1:32,000; kind gift by the Gershon laboratory (Margolis et al., 2016)), Iba1 (1:500, WAKO, Cat # 019-19741; Richmond, VA), CD68 (1:1000, Abcam, Cat # ab53444; Waltham, MA).

5




On Day 2, remove the primary antibody solution from each well and wash.

5.1 Wash with PBS-Tween (0.1%) for  00:10:00 . (1/3) 10m

5.2 Wash with PBS-Tween (0.1%) for  00:10:00 . (2/3) 10m

5.3 Wash with PBS-Tween (0.1%) for  00:10:00 . (3/3) 10m

6 Incubate in secondary antibody (1:1000) in blocking solution for  02:00:00 at 2h



 Room temperature , covered on a rotator.

7 Wash again.



7.1 Wash for  00:10:00 each in 0.1% PBS-Tween. (1/3) 10m

7.2 Wash for  00:10:00 each in 0.1% PBS-Tween. (2/3) 10m

7.3 Wash for  00:10:00 each in 0.1% PBS-Tween. (3/3) 10m

8 Mount on slide with vectashield medium with DAPI

9 Imaging:



- 9.1 For imaging enteric neurons, for each plexus collect 2-3, 2x2 tile z-stack 640.17x640.17 μm confocal images at 20x magnification.
- 9.2 For imaging macrophages, for each plexus collect 2-3, 2x2 tile z-stack 390.09 x 390.09 μm confocal images at 20x magnification.

10 Analysis:



- 10.1 Count the number of ANNA1⁺, TH⁺ cells, and IBA1⁺ cells for each stacked image using Fiji.
- 10.2 Within the SP, threshold the TH⁺ signal, then analyze mean fluorescent intensity (MFI) and the area of the TH signal. Keep the thresholding consistent across each image, animal, and condition within each experiment.
- 10.3 Within each animal, sum the the number of ANNA1⁺, TH⁺ cells, and IBA1⁺ cells across all images separately then divide by the acquisition area.
- 10.4 For each experiment, normalize to the CFA only condition.

