




May 23, 2022

# In Vitro B-galactosidase and B-hexosaminidase Activity Assay (Total Cell Lysate)

[Laura Smith](#)<sup>1</sup>

<sup>1</sup>Department of Clinical and Movement Neurosciences, Queen Square Institute of Neurology, University College London (UCL)

1

[dx.doi.org/10.17504/protocols.io.kqdg3p8r7l25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3p8r7l25/v1) Laura Smith

Synthetic substrate can be used at acidic pH (pH 4.1) to assess activity of lysosomal hydrolases,  $\beta$ -Galactosidase and  $\beta$ -Hexosaminidase. The substrate is cleaved by  $\beta$ -Galactosidase and  $\beta$ -Hexosaminidase and produces a fluorescent product proportional to activity.

DOI

[dx.doi.org/10.17504/protocols.io.kqdg3p8r7l25/v1](https://dx.doi.org/10.17504/protocols.io.kqdg3p8r7l25/v1)

Laura Smith 2022. In Vitro B-galactosidase and B-hexosaminidase Activity Assay (Total Cell Lysate). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kqdg3p8r7l25/v1>

 protocol ,

Apr 28, 2022

May 23, 2022

61620

## Reagents

1

Buffers: McIlvaine citrate-phosphate (MV), pH 5.4

Substrates:

- $\beta$ -Galactosidase: 4-Methylumbelliferyl  $\beta$ -D-galactopyranoside (Sigma M1633)
- $\beta$ -hexosaminidase: 4-Methylumbelliferyl N-acetyl- $\beta$ -D-glucosaminide (Sigma M2133):

Standard: 4-methylumbelliferone (MWt. 176)

Stopping solution: 0.25 M glycine buffer pH 10.4 reagent.

Inhibitors: CBE (Sigma 5424) and DNJ (Enzo BML-SL230-005).

#### Preparation of reagents

### 2 McIlvaine citrate-phosphate (MV):

A	B
	pH 4.1
0.1 M citric acid	30 ml
0.2 M Na <sub>2</sub> HPO <sub>4</sub>	20 ml, then 1:1

0.1M citric acid monohydrate (Mwt = 210.14 g/mol) – 5.2535 g in 250 mL dH<sub>2</sub>O

0.2M Na<sub>2</sub>HPO<sub>4</sub> (Mwt = 141.9 g/mol) – 7.098 g in 250 mL dH<sub>2</sub>O

#### Substrates:

### 3 $\beta$ -Galactosidase: 4-Methylumbelliferyl $\beta$ -D-galactopyranoside (Sigma M1633):

Take 3.4 mg in 10 mL dH<sub>2</sub>O and heat at 80°C until dissolved (0.34mg/ml)

Make up in two bijoux tubes; vortex to dissolve or sonicate and store at -20°C.

For each experiment, heat at 60-80°C in oven to ensure all powder is solubilised.

### 4 $\beta$ -hexosaminidase: 4-Methylumbelliferyl N-acetyl- $\beta$ -D-glucosaminide (Sigma M2133):

Take 27 mg in 10 mL dH<sub>2</sub>O and heat at 80°C until dissolved (2.71mg/ml)

Make up in two bijoux tubes of 5 mL and store at -20°C.

For each experiment, heat at 60-80°C in oven to ensure all powder is solubilised.

### 5 Standard: 4-methylumbelliferone (MWt. 176) (Sigma M1381):

Desired concentration is 1 nmol in 200  $\mu$ L aliquots.

Take 2 mg in 1.5 mL dH<sub>2</sub>O, then do a 1:1000 dilution in dH<sub>2</sub>O to give a 0.2 $\mu$ g/200 $\mu$ L concentration.

Aliquot into 200  $\mu$ Ls in eppendorfs.

When using in experiment, add 1 mL glycine buffer to Eppendorf and load 200  $\mu$ L into each well.

## 6 Stopping buffer: 0.25 M glycine buffer pH 10.4 reagent:

Make up 64g NaOH in 200 mL dH<sub>2</sub>O.

Make up 150g glycine in 1600 mL dH<sub>2</sub>O.

When all mixed, add both together.

Ensure pH is 10.4 and make up to 2 L with dH<sub>2</sub>O.

### Apparatus / Instrumentation

## 7 37°C water bath

Plate reader, excitation 360 nm, emission 460 nm, sensitivity=50

### Sample preparation

## 8 Enzyme:

Samples are resuspended in water or 1% (v/v) TX-100 in PBS. For TX-100 lysis, cells incubated on ice for 15 mins and debris/nuclei removed by centrifugation at 17,000 x g, 10 min, 4 °C.

Supernatant containing GBA enzyme placed in fresh tube. All samples are sonicated in water bath for 1 minute. Protein concentration measured with BCA protein assay.

### Method Protocol

## 9 Dilute in water a portion of the sonicate to give a protein concentration of 0.25 - 4 mg / ml.

## 10 Set up mix in eppendorf tubes for each well as follows:

Make a master mix e.g multiply by number of wells.

B-gal:

20µl B-gal/HEX buffer pH 4.1

10µl B-gal substrate solution

B-hex:

20µl B-gal/HEX buffer pH 4.1

10µl B-gal substrate solution

## 11 Add 10 µl diluted enzyme sample in to duplicate wells of a 96 well plate. (Note: may have to optimise volume of sample loaded depending on sensitivity of fluorescence machine used).

## 12 Add 10 µl of lysis buffer used (water of TX-100 in PBS) to duplicate wells to serve as substrate

blanks.

13 Add 30 µl of reaction mixture to each well

14 Incubate at 37°C for 30 minutes.

15 Add 200 µL stopping solution to each well

16 Standard:

To 1 nmol 4-methylumbelliferone standard in 200 µl H<sub>2</sub>O add  
1.0 ml stopping reagent. Mix and 200 ml to empty wells to serve as a fluorescence standard.

17 Read fluorescence, excitation 360 nm, emission 460 nm.

18 Calculate activity in nmol / hr / mg protein