



FEB 06, 2023

## 🌐 DNA extraction - Zooplankton - 96 wells

coline.royaux<sup>1,2,3</sup>, Nicolas Rabet<sup>1,2,3</sup>, Céline Bonillo<sup>1,2,3</sup>

<sup>1</sup>Sorbonne Université; <sup>2</sup>Muséum National d'Histoire Naturelle; <sup>3</sup>UMR BOREA



Coline Royaux

Université Pierre et Marie Curie (Paris VI), Muséum National...

### ABSTRACT

This protocol was used to extract DNA from whole or parts of zooplanktonic freshwater crustaceans (Copepoda, Branchiopoda, ...) from New Caledonia.

OPEN  ACCESS

#### DOI:

[dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1](https://dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1)

**Protocol Citation:** coline.royaux, Nicolas Rabet, Céline Bonillo 2023. DNA extraction - Zooplankton - 96 wells.

**protocols.io**

<https://dx.doi.org/10.17504/protocols.io.j8nlk4x7wg5r/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Oct 15, 2021

**Last Modified:** Feb 06, 2023

**PROTOCOL integer ID:**  
54125

- 1 Prepare your 96-well extraction plate with one individual per well. Alternate genus in the wells to detect eventual contamination between wells.

1d

**1.1** Collect one individual from a sample 1m

**1.2** Note its genus and determine its sex with a binocular microscope 1m

**1.3** For big individuals (more than 5 mm), dissect a few legs and put it in the well. Be careful not to damage the rest of the body and put it in a tagged Eppendorf tube.  
For little individuals (less than 5 mm), put the whole body. 3m

**Note**

If necessary, use alcohol to get the biological material to fall at the bottom of the well

**1.4** When all 96 wells are filled, the biological material has to dry to go to lysis 12h

**Note**


If necessary, use a micropipette to empty an excess of alcohol in the well

**Safety information**

Make sure the plate is closed when you want to transport it elsewhere

**2** Prepare the lysis 15m

**2.1** Mix  18 mL T1 buffer and  2.5 mL K proteinase in a Multi-Channel Reservoir and 10m

distribute  200 µL of the mix with a multimitropipette in each well

#### Note

 180 µL T1 buffer and  25 µL K proteinase in each well

**2.2** Close your extraction plate with a heated aluminium foil and an adhesive plastic film

3m


**3** Put your plate in a proofer at  56 °C  Overnight (6h or more) to lyse the tissues

6h

**4** Perform the DNA extraction with a DNA extraction robot








**4.1** Remove the adhesive film and aluminium foil from the plate and put it in the robot

**4.2** It will deposit  200 µL BQ1 buffer and  200 µL ethanol in each wells and mix it

**4.3** Then,  600 µL of the wells content (lysate, BQ1, ethanol) are transfered on the tissue binding plate. Reagents excess are emptied in a waste container.

**4.4** The tissue binding plate is then dried by a  00:05:00 aspiration to bind DNA to the silica membrane of the binding plate

5m

- 4.5** The silica membrane is then washed with  600 µL BW buffer and twice with  900 µL B5 buffer per well. Each wash is intercalated by a  00:05:00 aspiration dry. 5m
- 4.6** The waste container is then removed from under the binding plate which is dried again by a  00:10:00 aspiration 10m
- 4.7** An empty extraction plate is placed under the binding plate to retrieve the genomic DNA from it
- 4.8** DNA is eluted from the tissue binding plate with  100 µL BE buffer in each well and is collected in the new extraction plate underneath
- 4.9** After a  00:03:00 rest, the binding plate is dried for  00:02:00 and the elution is repeated with  100 µL BE buffer 5m
- 4.10** Retrieve the new extraction plate containing the genomic DNA and discard the rest