



NOV 07, 2022



IN DEVELOPMENT

 Vector Linearization**This protocol is published without a DOI.**[Felipe FE Edaes](#)¹¹Harvard Felipe FE Edaes

COMMENTS 0

ABSTRACT

Protocol for the linearization of a previously obtained vector.

PROTOCOL CITATION

Felipe FE Edaes 2022. Vector Linearization . **protocols.io**
<https://protocols.io/view/vector-linearization-civ8ue9w>



LICENSE

_____ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 06, 2022

LAST MODIFIED

Nov 07, 2022

PROTOCOL INTEGER ID

72352

GUIDELINES

Add phosphatase to remove phosphate, which remains from the end of linear vectors, thus preventing cells from relinearizing the vectors.

Equipment

1

Equipment

NanoDrop™ One UV-Vis Spectrophotometer

NAME

spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU



<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

LINK


Sample Volume (Metric): Minimum 1 µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

SPECIFICATIONS

Measure DNA Concentration of the Plasmid

2 Use  1.0 µL to  1.5 µL of the elution buffer (used to elute the DNA) as blank.

3 Measure blank once to confirm the accuracy.

4 To measure DNA sample, add  1-1.5 µL (equivalent to blank)

Calculate the amount of Plasmid to be used

5 If $1\ \mu\text{g}$ is $1000\ \text{ng}$, this value should be divided by the amount of DNA measured.

5.1 If the required amount is $3\ \mu\text{g}$ and the measured sample was $450\ \text{ng}/\mu\text{L}$, the proper calculation should be: $3000/450 = \sim 6.67 = 7.00\ \mu\text{L}$

3h

Vector Linearization

6 $3\ \mu\text{g}$ Plasmid DNA

7 $1\ \mu\text{L}$ Enzyme A

8 $1\ \mu\text{L}$ Enzyme B

9 $2\ \mu\text{L}$ 10X Buffer

10 Remaining H₂O $20\ \mu\text{L}$ (up to)

3h

11 Incubate the reaction for $03:00:00$ at $37\ ^\circ\text{C}$, then place it On ice.

11.1 Place it On ice and add $0.5\ \mu\text{L}$ **phosphatase** to the linearized vector.

