



Nov 23, 2020

Toxicity Assay

Elizabeth Fozo¹¹In-house protocol

1

Works for me

This protocol is published without a DOI.

Eadewunm

ABSTRACT

Toxicity Assay (based on growth curves protocol)

PROTOCOL CITATION

Elizabeth Fozo 2020. Toxicity Assay. **protocols.io**
<https://protocols.io/view/toxicity-assay-bpwhmpb6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 19, 2020

LAST MODIFIED

Nov 23, 2020

PROTOCOL INTEGER ID

44713

DISCLAIMER:

DISCLAIMER: THIS IS A WORK IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Toxicity Assay (based on growth curves protocol)

Rules of thumb

- 1 Never remove > 1/3 of a culture for a growth curve
- 2 Read OD600 (optical density at 600 nm) of 200 microliters of culture added to 800 microliters of fresh media -> 1,000 microliters of fluid in cuvette; = 1:5 dilution

Remember to set the factor equal to 5 when reading from the spectrometer

- 3 Easier to measure 800 microliters of fresh media into cuvettes ahead of time; simply store in cuvette box
- 4 Easier to measure media into a 50 mL conical for filling fresh cuvettes; less likely to get a contaminated bottle
- 5 Easier to use same blank multiple times; suggest every time you refill conical with media for cuvettes, make a new blank with that and use with cuvettes from that mixture.
- 6 Place cuvettes into 10% bleach solution when finished.
*This is a must for BSL 2 organisms such as *E. faecalis* and EHEC.

Making an overnight culture

- 7 Measure 5mL appropriate media with appropriate antibiotics into fresh conical (make sure to label conical with date, strain!)
- 8 Take 1 colony from a plate with appropriate strain (use pipette to poke that colony)
- 9 Inoculate media
- 10 Grow overnight (usually start culture ~3-5 p.m. night before culture is needed) In the morning:

In the morning

- 11 Measure OD600
- 12 Use $V_1C_1 = V_2C_2$ to make a 0.01 OD600 dilution in media for growth curve
- 13 Place culture in flasks, mix, place in the appropriate incubator, start count up timer.

Toxicity Assay:

- 14 The first time point is at 30 min on the count-up timer.
- 15 Take time point every 30 min, write down data and graph as you go

- 16 Use count up timer to measure time. When it reaches 30 min, 1 hour, 1.5 hours, etc. remove cultures from the incubator, pipette into cuvettes (remember 1/5 dilution, e.g. 200 microliters culture in 800 microliters media)
- 17 Remember to mix culture gently before adding to a cuvette (the culture can be clumpy at the bottom, and you want a homogenous solution)
- 18 Once the culture reaches OD 0.3 (at about 2-3 hours for E. coli in LB), split cell culture in half.

Try to reach as close to OD 0.3 as possible. May require you to take samples off the 30-minute cycle.
- 19 Induce $\frac{1}{2}$ cells with the appropriate inducer (for example, arabinose) at the appropriate final concentration (for ara, 0.2% is saturating).
- 20 Continue incubation for all cultures, taking OD600 every 30 minutes after induction.