



Version 1

Nov 30, 2020

Mollusk pedal mucus effects on epilithic biofilms V.1

Clara María Arboleda-Baena¹, Rodrigo De La Iglesia¹, Claudia Pareja¹, Isadora Pla¹¹Pontificia Universidad Catolica de Chile*In Development*

This protocol is published without a DOI.

The microscopic world, trophic and non trophic pathways in the intertidal rocky shore.Clara María Arboleda-Baena
Pontificia Universidad Catolica de Chile

PROTOCOL CITATION

Clara María Arboleda-Baena, Rodrigo De La Iglesia, Claudia Pareja, Isadora Pla 2020. Mollusk pedal mucus effects on epilithic biofilms. **protocols.io**<https://protocols.io/view/mollusk-pedal-mucus-effects-on-epilithic-biofilms-bp8kmruw>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 30, 2020

LAST MODIFIED

Nov 30, 2020

PROTOCOL INTEGER ID

45036

MATERIALS TEXT

Laminar flow cabinet

Polycarbonate cell culture plate of 6-Wells

Sterile glass slide (20 x 7 cm)

1 L Sterile seawater

Sterile scalpel

Cryovials

0.01% acridine orange

500 mL PBS

Fluorescent Microscope with an excitation filter (FS38) 470/40 nm and emission filter 525/50 nm.

Image J

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

- 1 We cultured epilithic biofilms on 24 cover glass slides inside a Polycarbonate cell culture plate of 6-Wells, for one week in K medium (Keller et al. 1987).
- 2 We collected, during nocturnal low tides, animals of each mollusk species from wave-exposed platforms and brought them to the laboratory in coolers and then placed to acclimatize for a week in separate aquariums by circulating seawater and constant aeration.
- 3 Collection of pedal mucus was done under a laminar flow cabinet (Connor 1986). Animals were carefully removed from their containers, washed in sterile-filtered seawater (0.2 µm pore-size filters), and then, each of them, was placed on an individual inclined sterile glass slide (20 x 7 cm).
- 4 Sterile seawater was added to stimulate individual movement and mucus production as they moved through the glass. The animals moved in the glass slides for a maximum of five minutes and then were removed.
- 5 The mucus on the replicate glass slides was removed with a sterile scalpel and put it in individual cryovials with sterile seawater.
- 6 The experiment consisted of placing the pedal mucus collected from the individuals, of each of the species, in separate replicated wells that had been cultured with biofilm. The control wells received no mucus and treatments were randomly assigned to the cell culture plates.
- 7 After one week, we stained the cover glass slides with 200 µl 0.01% acridine orange for 5 minutes, and then we removed the biofilm.

The acridine orange, AO, is a fluorescent compound that emit red fluorescence when attached to single-stranded templates (RNA) and green fluorescence if the nucleic acids are double-stranded (DNA) (Rigler 1966) and can be used to get gross estimates of the total amount of bacterial cells in the biofilm.

After 3 min incubation in the dark, the staining solution was removed, and the plate was washed twice with 500 µl of PBS solution.
- 8 We took five photographs of each cover glass slide under a Fluorescent Carl Zeiss AXIO Scope A1 Microscope using excitation filter (FS38) 470/40 nm and emission filter 525/50 nm.
- 9 Then, from the photographs, we measured the cover of the photosynthetic epilithic biofilm using the software Image J and the workflow:
- 10 Covers under different treatments were analyzed with one-way ANOVA with grazer species as a fixed factor. A Tukey post-hoc test was performed to determine the pattern of differences.

