



JAN 22, 2024

Ubiquitin immunoprecipitation using an anti-ubiquitin nanobody

Cole S Sitron¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry

ASAP Collaborative Research Network



Cole Sitron

ABSTRACT

This protocol describes a method to detect ubiquitination on a protein of interest. This technique relies on immunoprecipitation (IP) of ubiquitinated proteins from a cell lysate through the use of an anti-ubiquitin nanobody coupled to agarose beads. The eluate can be run on SDS-PAGE to determine whether the protein of interest was recovered in the IP and, therefore, ubiquitinated.

ATTACHMENTS

[io_Ubiquitin immunoprecipitation using an anti.docx](#)

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.dm6gp36k1vzp/v1

Protocol Citation: Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Ubiquitin immunoprecipitation using an anti-ubiquitin nanobody.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.dm6gp36k1vzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: Jan 16, 2024

Last Modified: Jan 22, 2024

PROTOCOL integer ID: 93785

Funders Acknowledgement:
Aligning Science Across
Parkinson's
Grant ID: ASAP-000282

MATERIALS

Materials

Cell culture

2-4 50-80% confluent wells of cells growing in a 6-well plate.

Equipment

Equipment	
Bioruptor® Plus sonication device	NAME
Bioruptor®	BRAND
B01020001	SKU
https://www.diagenode.com/en/p/bioruptor-plus-sonication-device	

- CLARIOstar Plus Plate Reader (BMG Labtech) (or equivalent)
- Microcentrifuge

Specialized Reagents

- Ubiquitin Selector NanoTag Biotechnologies Catalog #N2510
- MiniSpin Columns (50 units) NanoTag Biotechnologies Catalog #A1001-L
- TrypLE® Express Enzyme (1X), phenol red Thermo Fisher Catalog #12605036
- Pierce® Rapid Gold BCA Protein Assay Kit Thermo Fisher Catalog #A53225

Buffers

- Triton Lysis Buffer:

A	B
Tris pH 8	20 mM
NaCl	150 mM
Triton X-100	1%
cOmplete, Mini EDTA-free Protease Inhibitor Cocktail (Roche cat. no. 11873580001)	1X
Freshly-added 20 mM N-ethylmaleimide (2.5 mg/ml; Sigma Aldrich E3876-25G)	2.5 mg/ml
50 uM PR-619 (from 100 mM DMSO stock; Sigma-Aldrich 662141-25MG)	50 uM
Benzonase (Max Planck Institute of Biochemistry Core Facility)	7.5 U/ml

⊗ cOmplete™ EDTA-free Protease Inhibitor Cocktail **Roche Catalog #11873580001**

⊗ N-Ethylmaleimide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E3876**

⊗ DUB Inhibitor V, PR-619 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #662141**

- Triton Wash Buffer:

A	B
Tris pH 8	20 mM
NaCl	150 mM
EDTA	1 mM
Triton X-100	1%
SDS	0.1%
PR-619	50 uM

- 4X

⊗ NuPAGE™ LDS Sample Buffer (4X) **Thermo Fisher Scientific Catalog #NP0007** +

5% beta-mercaptoethanol (

⊗ 2-Mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**


-100ml)

- 10% FBS (⊗ Fetal Bovine Serum **Gibco - Thermo Fischer Catalog #10270106**) in 1X PBS (diluted from

⊗ PBS (10X), pH 7.4 **Thermo Fisher Scientific Catalog #70011051**)

- 1X PBS

Cell collection and lysis

1 Remove the medium from the wells and trypsinize with  500 μ L TrypLE Express.



2 Quench the TrypLE Express with  500 μ L 10% FBS and move the cells into a centrifuge tube.






3 Pellet cells at  1500 x g for  00:03:00 at  4 $^{\circ}$ C .



3m

4 Wash cells with  1 mL PBS and pellet again.



5 Resuspend cells in  150 μ L Triton Lysis Buffer and incubate  On ice for  00:05:00 .



5m



6 Sonicate in the BioRuptor for 5 cycles of  00:00:30 on and  00:00:30 off at  4 $^{\circ}$ C .

1m

7 Clarify the lysate by centrifugation  18000 x g for  00:10:00  4 $^{\circ}$ C .



10m

8 During the centrifugation, prepare BSA standards according to the Rapid Gold BCA Protein Assay Kit manual and begin equilibrating  50 μ L of Ubiquitin pan-selector resin (pipette with a cut p200 tip) in  500 μ L



of lysis buffer.

9

Collect supernatant and place into a new tube On ice .



10

Prepare a small 1:10 dilution of samples and perform Rapid Gold BCA Protein Assay according to manufacturer's instructions.



Immunoprecipitation (IP)

1h 7m

11

Dilute samples into Triton Lysis Buffer in two tubes: 1 mg in 500 μ L (for IP) and 100 μ g in 25 μ L (for input).



12

Pellet the beads at 1000 x g for 00:02:00 at 4 $^{\circ}$ C .

2m



13

Pull off the supernatant, add the 1 mg of lysate to the beads, and incubate with rotation at 4 $^{\circ}$ C for 01:00:00 .



14

During the incubation denature the input sample with 25 μ L 4X NuPAGE LDS Sample Buffer at 95 $^{\circ}$ C 00:05:00 .

5m



15

Pellet the beads at 1000 x g for 00:02:00 at 4 $^{\circ}$ C .

2m



16 Wash the beads with either  1 mL Triton Wash Buffer.



17 Pellet the beads and repeat for a total of 3 washes.



18 Remove the supernatant and add  100 μ L of 2X NuPAGE LDS Sample Buffer.



19 Elute by boiling at  95 $^{\circ}$ C for  00:05:00 .

5m

20 Transfer the beads to a MiniSpin column placed in a tube.



21 Centrifuge  3000 x g for  00:03:00 .

3m



22 The eluate has now been collected in the tube and is ready for SDS-PAGE analysis.

