



MAR 19, 2024

DNA Barcoding Protocol for Arthropods

Valerie Warhol¹

¹Science Club



Valerie Warhol
Science Club

ABSTRACT

This is a basic protocol for doing extraction and amplification of DNA barcodes (COI) for arthropods (typically small insects or spiders).

PROTOCOL MATERIALS



Guanidine Hydrochloride 6M **Carolina Biological Supply Catalog #C33427**

Step 5

OPEN  ACCESS



Protocol Citation: Valerie Warhol 2024. DNA Barcoding Protocol for Arthropods.

protocols.io

<https://protocols.io/view/dna-barcoding-protocol-for-arthropods-dawi2fce>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited




Protocol status: Working
We use this protocol and it's working

Created: Mar 19, 2024

Last Modified: Mar 19, 2024




PROTOCOL integer ID: 96938

Prepare sample and equipment

- 1 Make sure all instruments, such as forceps and pestle, are clean and sterile.
- 2 Prepare water bath at  65 °C .
- 3 Dissect sample from specimen (typically 1 leg). Return specimen to freezer.
- 4 Let sample dry for 5–10 minutes to remove ethanol.
- 5 Prepare a clean 1.5 mL hinged tube by writing sample ID on it and filling with  250 µL of  Guanidine Hydrochloride 6M **Carolina Biological Supply Catalog #C33427**

Lyse cells

11m


- 6 Put sample in tube. Grind sample with pestle until broken up into tiny pieces.
- 7 Incubate sample tube in  65 °C water bath for  00:10:00 .
- 8 Remove tube and lower temperature of water bath to  57 °C .

10m

9 Centrifuge tube for 00:01:00 at maximum speed to pellet debris.


1m

10 Remove  Silica Resin **Carolina Biological Supply Catalog #C33426** from refrigerator.

11 Label a clean 1.5 μ L tube with sample number. Transfer  150 μ L of the supernatant to the clean tube. Discard old tube containing debris.

Bind DNA

5m 30s

12 Add  3 μ L of silica resin to tube. Mix well by pipetting up and down several times.

13 Close tube and incubate for 00:05:00 in  57 °C water bath.

5m

14 Centrifuge for 00:00:30 at maximum speed to pellet the resin.


30s

15 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet.

Wash


1m

- 16 Remove molecular grade water from refrigerator and

 Wash Buffer **Carolina Biological Supply Catalog #C33428**


 from freezer.

- 17 Add

 500 µL

 of ice-cold wash buffer to the pellet. Mix well by pipetting up and down several times to resuspend the silica resin.


- 18 Close the tube and centrifuge for

 00:00:30

 at maximum speed to pellet the resin. 30s


- 19 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet.

- 20 Again, add

 500 µL

 of ice-cold wash buffer to the pellet. Mix well by pipetting up and down to resuspend the silica resin.

- 21 Close the tube and centrifuge for

 00:00:30





 at maximum speed to pellet the resin. 30s

- 22 Return wash buffer to freezer.

- 23 Use a pipette with a fresh tip to remove the supernatant, being careful not to disrupt the pellet. Spin the tube briefly to collect any remaining drops of supernatant, and then remove these with a pipette.

Elute DNA

30s

- 24 Add  100 μ L of molecule grade water to the silica resin and mix by pipetting up and down several times.
- 25 Incubate the mixture at  57 $^{\circ}$ C for 5 minutes.
- 26 Centrifuge for  00:00:30 at maximum speed to pellet the resin. 30s
- 27 Label a clean 1.5 μ L tube with sample number. Transfer  90 μ L of the supernatant to the clean tube, being careful not to disturb the pellet. Discard old tube containing the resin.
- 28 Store sample in freezer until ready to PCR. If going directly to PCR, put sample in refrigerator.

Amplify COI (PCR)

10m

- 29 For each DNA sample, label a PCR microtube with sample ID.
- 30 Turn on PCR thermal cycler and connect to computer.

31 Remove molecular grade water from refrigerator. Remove template DNA, PCR master mix, and primers from freezer. Let thaw for 00:10:00 . 10m

A	B	C	D
Primer Name	Direction	Sequence	Concentration
LCO1490	Forward	GGTCAACAAATCATAAAGATATTGG	10 µM
HCO2198	Reverse	TAAACTTCAGGGTGACCAAAAAATCA	10 µM

32 Add 32 µL molecular grade water to each microtube.

33 Mix forward primer by flicking. Add 1.5 µL forward primer to each microtube.

34 Mix reverse primer by flicking. Add 1.5 µL reverse primer to each microtube.

35 Mix template DNA by flicking. Add 5 µL template DNA to each microtube.

36 Mix EZ PCR Master Mix 5X miniPCR by inverting. Add 10 µL PCR master mix to each microtube.

37 Put microtubes in PCR thermal cycler and start touchdown PCR.

1. Denature: 95° C 30 sec; Anneal: 60° C 30 sec; Extend: 72° C 45 sec.
2. Denature: 95° C 30 sec; Anneal: 59° C 30 sec; Extend: 72° C 45 sec.
3. Denature: 95° C 30 sec; Anneal: 58° C 30 sec; Extend: 72° C 45 sec.
4. Denature: 95° C 30 sec; Anneal: 57° C 30 sec; Extend: 72° C 45 sec.

5. Denature: 95° C 30 sec; Anneal: 56° C 30 sec; Extend: 72° C 45 sec.
6. Denature: 95° C 30 sec; Anneal: 54° C 30 sec; Extend: 72° C 45 sec.
7. Denature: 95° C 30 sec; Anneal: 53° C 30 sec; Extend: 72° C 45 sec.
8. Denature: 95° C 30 sec; Anneal: 52° C 30 sec; Extend: 72° C 45 sec. Repeat this step for 28 cycles.
9. Final extension: 72° C 5 min.

(An initial long denature step is not needed for this protocol.)

38 Remove microtubes from PCR thermal cycler. Verify PCR using gel electrophoresis.

39 Store PCR products in freezer ( -20 °C).