




NOV 16, 2023

🌐 Cell isolation from dorsal mouse skin for single-cell RNA-seq

Stella Yue Ting
Cenk Celik¹, Lee¹, Guillaume
Thibault¹

¹Nanyang Technological University



Cenk Celik
Nanyang Technological University

ABSTRACT

We identified heterogeneity of the skin cell populations upon *Enterococcus faecalis* infection.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1

Protocol Citation: Cenk Celik, Stella Yue Ting Lee, Guillaume Thibault 2023. Cell isolation from dorsal mouse skin for single-cell RNA-seq.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.yxmvmn8m9g3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 08, 2022

Last Modified: Nov 16, 2023

PROTOCOL integer ID:
68336

Funders
Acknowledgement:
NMRC
Grant ID: MOH-000566

MATERIALS

Consumables

- Hank's Balanced Salt Solution $\text{Ca}^{2+}/\text{Mg}^{2+}$ -free
- Dulbecco's Phosphate Buffered Saline $\text{Ca}^{2+}/\text{Mg}^{2+}$ -free
- Dispase
- Collagenase type I
- Liberase TM Research Grade
- 0.25% Trypsin/EDTA
- Bovine Serum Albumin
- Trypan Blue
- Ice

Instruments

- CO_2 supplemented 37°C incubator
- Chromium Instrument/X, 10X Genomics
- Countess 3 Automated Cell Counter, Invitrogen
- 15ml and 50ml tubes-compatible centrifuge
- 1.5ml tube-compatible centrifuge

Surgical tools

- 6-mm punch
- Scissors
- Forceps
- Bistouri
- Blades G23
- 10mL syringes
- 21G needles

Other labware

- 50ml centrifuge tubes
- 15ml centrifuge tubes
- 1.5ml Protein LoBind tubes
- Disposable haemocytometer (4 per sample)
- 10mm dishes
- 6-well plates
- Pasteur pipettes

PROTOCOL MATERIALS

✕ Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Thermo Fisher Scientific Catalog #88284

In 6 steps

✕ Liberase TM Merck MilliporeSigma (Sigma-Aldrich) Catalog #000000005401119001

Step 2

✕ Bovine Serum Albumin Merck MilliporeSigma (Sigma-Aldrich) Catalog #A7030

Step 3

✕ Standard Biopsy Punches, Disposable Standard biopsy punch; 6mm Thermo Fisher Catalog #12460412

Step 7

✕ scalpel blades

Step 8

✕ Falcon® 40 µm Cell Strainer Corning Catalog #352340

Step 21

✕ Dispase II, powder Thermo Fisher Catalog #17105041

Step 0.1

✕ 1X Dulbecco's Phosphate Buffered Saline (DPBS) Thermo Fisher Scientific Catalog #14190094

In 2 steps

✕ Collagenase Type I powder Thermo Fisher Scientific Catalog #17100017

Step 1

✕ Trypsin EDTA Gibco - Thermo Fischer Catalog #25-051-Cl.

Step 6

✕ 10 mL syringes Becton Dickinson (BD) Catalog #BD 309695

In 2 steps

✕ Falcon® 70 µm cell strainer Corning Catalog #352350

In 2 steps

Enzyme solutions

20m

1 Dispase Solution

2m

Dissolve 0.5g of ✕ Dispase II, powder Becton Dickinson (BD) Catalog #17105041 in a 50ml of

✕ 1X Dulbecco's Phosphate Buffered Saline (DPBS) Becton Dickinson (BD) Catalog #14190094 to

make 10 mg/mL of Dispase II solution

0.2 Filter the solution by using a 0.2µm filter

0.3 Aliquot 1.25mL in tubes and store at -20 °C

1 Collagenase type I Solution

Dissolve 1g of Collagenase Type I powder Becton Dickinson (BD) Catalog #17100017 in 4mL of

Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson (BD) Catalog #88284 to

make a 250 mg/mL collagenase type I solution

1.1 Make a stock solution

i.e. $67500\text{U/mL} \times 4\text{mL} = 50\text{U}/\mu\text{L} \times V_{\text{max}}$

$V_{\text{max}} = 5.4\text{mL}$


$5.4 - 4.0 = 1.4\text{mL}$ to add the stock to make the stock solution

1.2 Filter the solution by using a 0.2µm filter

1.3 Aliquot 20µL in tubes and store at -20 °C

2 Liberase TM Solution


Punch through the cap with a syringe to add 2mL of

 Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson (BD) Catalog #88284

into 5mg  Liberase TM Becton Dickinson (BD) Catalog #000000005401119001 bottle to make

[M] 2.5 mg/mL liberase solution

2.1 Dilute into [M] 0.2 mg/mL in 2m

 Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson (BD) Catalog #88284

(1.04 Wünsch unit/mL)


2.2 Filter the solution by using a 0.2µm filter 2m


2.3 Aliquot 750µL in tubes and store at  -20 °C 5m

Working Solutions 2m

3 Working Solution 1 (WS1) 2m

Make a [M] 0.5 Mass / % volume of


 Bovine Serum Albumin Becton Dickinson (BD) Catalog #A7030 in 50mL of


 Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson (BD) Catalog #88284

3.1 Filter the solution by using a 0.2µm filter and keep it at  4 °C 2m

4 Working Solution 2 (WS2) 2m

Make a [M] 0.04 Mass / % volume of Bovine serum albumin solution in 50mL of

 1X Dulbecco's Phosphate Buffered Saline (DPBS) Becton Dickinson
(BD) Catalog #14190094

4.1 Filter the solution by using a 0.2µm filter and keep it at  4 °C

2m

5 Enzyme Cocktail 1 (EC1)–4 samples

3m


From the enzyme stocks that have been prepared before, mix:

- 1.25mL of Dispase
- 0.75mL of Liberase
- 20µL of type I Collagenase

in a 15mL tube

5.1 Add 8.0 mL of

1m


 Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson
(BD) Catalog #88284


5.2 Warm up in a  37 °C water bath before using

30m

6 Enzyme Cocktail 2 (EC2)–4 samples

2m

Dilute the  Trypsin EDTA Becton Dickinson
(BD) Catalog #25-051-CI. in

 Hanks Balanced Salt Solution (HBSS) without Ca²⁺ Mg²⁺ Becton Dickinson
(BD) Catalog #88284 to











make a  0.05 Mass / % volume of trypsin/EDTA solution










6.1 Warm up in a  37 °C water bath before using

30m


Tissue Dissociation and Cell Isolation Procedure

7m

- 7 Punch out the skins from mice using a sterile  Standard Biopsy Punches, Disposable Standard biopsy punch; 6mm Becton Dickinson (BD) Catalog #12460412 and immediately float on WS1 on ice 5m
- 8 Peel off any fat tissue beneath the skin using a blade  scalpel blades Becton Dickinson (BD) 2m
- 9 Mince the fat tissue-cleaned skin samples into smaller pieces using sterile scissors while holding the tissue with sterile forceps on ice 2m
- 10 Incubate the minced tissues in 2 ml of EC1 in a 6-well plate at a  37 °C incubator with a 5% CO₂ supplementation for two hours 2h
- 
 02:00:00
- 10.1 Shake the plate orbitally in the incubator every  00:15:00 for better digestion 15m
- 
- 11 Remove a piston from a sterile,  10 mL syringes Becton Dickinson (BD) Catalog #BD 309695 to smash the partly-digested tissues in the plate 1m
- 12 Place a  Falcon® 70 µm cell strainer Becton Dickinson (BD) Catalog #352350 into a 50mL tube and sift through the cell/tissue mixture 1m
- 13 Wash the EC1-digested tissue on the cell strainer with WS1 thrice 3m
- 

- 14 Save the flow-through on ice
- 15 Invert the cell strainer with tissue remnants and place it on a sterile 6-well plate 1m
- 16 From the sieve-through side of the cell strainer, add 2mL of EC2 aiming at tissue remnants 1m
- 17 Incubate the undigested tissue at a  37 °C incubator with a 5% CO₂ supplementation for  00:15:00 15m
- 

- 18 Remove a piston from a sterile,  10 mL syringes Becton Dickinson (BD) Catalog #BD 309695 to smash the partly-digested tissues in the plate 1m
- 19 Place a sterile  Falcon® 70 µm cell strainer Becton Dickinson (BD) Catalog #352350 into the 50mL tube obtained at Step 15 to pool and sift through the cell/tissue mixture 1m
- 
- 20 Wash the trypsin-digested tissue on the cell strainer with WS1 thrice 1m
- 
- 21 Place a sterile  Falcon® 40 µm Cell Strainer Becton Dickinson (BD) Catalog #352340 into a clean 50mL tube and sift through the cell suspension 1m

22

Centrifuge the pooled cell suspensions at  300 x g, 20°C, 00:10:00

10m

23

Carefully aspirate the supernatant and resuspend the pellet in 1mL of WS2

2m



24

Centrifuge the pooled cell suspensions again at  300 x g, 20°C, 00:05:00

5m

25

Carefully aspirate the supernatant and resuspend the pellet in 500µL of WS2

1m



26

Determine the cell viability by mixing 10µl of the cell suspensions with 10µl of trypan blue

2m



Equipment	
Countess 3 FL Automated Cell Counter	NAME
Automated Cell Counter	TYPE
Thermofisher scientific	BRAND
AMQAF2000	SKU
https://www.thermofisher.com/th/en/home/life-science/cell-analysis/cell-analysis-instruments/automated-cell-counters/models/countess-3-fl.html	LINK

26.1

Take at least four counts per sample

8m



26.2 Remember to tick "Trypan Blue correction" on the Countess cell counter

1m



27 Average cell numbers for each sample and dilute each cell suspension with WS2 to adjust the cell number to 700-1200 cells/ μ L as indicated in [10X Genomics Chromium Next GEM Single Cell 3' Reagent Kits v3.1 Protocol \(CG000315\)](#)

5m

27.1 Aim for cell viability above 70% for each sample



Expected result

Cell viability: 70% per count

27.2 If the cell viability is below 70%, you may choose to remove dead cells using a dead cell removal kit



28 Immediately proceed with [Chromium Next GEM Single Cell 3' Reagent Kits v3.1 Protocol \(CG000315\)](#)