



SEP 23, 2023

Mitochondrial isolation protocol

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna (AT)

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes how to isolate crude mitochondrial fractions from HeLa cells.

ATTACHMENTS

[752-1920.pdf](#)

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.kqdg3x4zzg25/v1

Protocol Citation:

Elias Adriaenssens 2023. Mitochondrial isolation protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kqdg3x4zzg25/v1>

MANUSCRIPT CITATION:

The methodology is based on Wieckowski et al. 2009 Nat Protocols. Wieckowski, M., Giorgi, C., Lebiedzinska, M. et al. Isolation of mitochondria associated membranes and mitochondria from animal tissues and cells. Nat Protoc 4, 1582–1590 (2009).

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status:

Working
We use this protocol and it's working

Created: Jun 26, 2023

Last Modified: Sep 23, 2023

PROTOCOL integer ID:
84020

Keywords: Mitochondrial isolation, HeLa cells

MATERIALS

MATERIALS

Mitochondrial isolation buffer:

A	B
HEPES (pH 7.4)	5 mM
mannitol	250 mM
EGTA	0.5 mM

make fresh on the day itself!

Oligomycin/Antimycin A cocktail:

A	B
Oligomycin	10 μ M
Antimycin A	4 μ M


In case cells were treated for more than 8 hours, we also added 10 μ M Q-VD-OPh (A1901, ApexBio) to suppress apoptosis.

RIPA buffer:


A	B
Tris-HCl pH 8.0	50 mM
NaCl	150 mM
sodium deoxycholate	0.50%
SDS	0.10%
NP-40	1%


supplemented by cOmplete EDTA-free protease inhibitors (11836170001, Roche) and phosphatase inhibitors (PhosphoSTOP, 4906837001, Roche).

 Oligomycin A ApexBio Technology Catalog #A5588

 Antimycin A from Streptomyces sp. Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8674

 Q-VD-OPh ApexBio Technology Catalog #A1901

 cOmplete™, Mini, EDTA-free (Protease Inhibitor) Roche Catalog ##11836170001

 Roche PhosSTOP™ Merck MilliporeSigma (Sigma-Aldrich) Catalog #4906837001

Mitochondrial isolation



1h 15m

- 1 Seed HeLa cells in 15 cm dishes and grow until confluence. Treat the cells with DMSO or the mitophagy-inducing cocktail Oligomycin/Antimycin A (O/A) if needed.

- 2 Collect HeLa cells by trypsinization and resuspension in DMEM medium. Centrifuge the cells at  300 x g, 4°C, 00:10:00 .

10m



- 2.1 Resuspend the cell pellet in  1 mL of ice-cold PBS (1x) and spin again for  300 x g, 4°C, 00:05:00 .

5m

- 3 Remove the PBS and resuspend the cell pellet in  1 mL mitochondrial isolation buffer, which was cooled to  4 °C .


Note

Note that the mitochondrial isolation buffer is prepared fresh on the day itself.

- 4 To lyse cells without damaging the mitochondria, pipet the  1 mL cell suspension up and down (15-20x) with a 26.5G needle.

Note

This should lyse the plasma membrane but leave organelles intact.

- 5 Spin the suspension down at  600 x g, 4°C, 00:10:00 .

10m





- 5.1 Keep the supernatant and centrifugation  600 x g, 4°C, 00:10:00 .

10m


Note

The mitochondria are located in the supernatant at this speed, the pellet contains intact cells, cell nuclei, and other large cell debris.

- 6 After two spins at  600 x g , subject the supernatant to  7000 x g, 4°C, 00:10:00 .

10m



- 6.1 The pellet contains the mitochondria, so remove the supernatant and resuspend the pellet in  1 mL mitochondrial isolation buffer.

- 6.2 Then, centrifugation  7000 x g, 4°C, 00:10:00 .

10m

Note

The supernatant from the first spin at $7000 \times g$ can be stored as a cytosolic fraction.

- 7 After two spins at $7000 \times g$, subject the resuspended pellet to $10000 \times g$, 4°C , 00:10:00.

10m



- 7.1 The pellet contains the mitochondria, so remove the supernatant and resuspend the pellet in 1 mL mitochondrial isolation buffer.

- 7.2 Then, centrifugation $10000 \times g$, 4°C , 00:10:00.

10m

- 8 After two spins at $10000 \times g$, the pellet can be resuspended in RIPA lysis buffer or used for further procedures such as mitochondrial import assays, proteinase K protection assays, etc.