



SEP 23, 2023

OPEN ACCESS



DOI:
[dx.doi.org/10.17504/protocol
s.io.8epv5x11ng1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5x11ng1b/v1)

Protocol Citation: Elias
Adriaenssens 2023. Plasmid
Construction and Gibson
cloning. **protocols.io**
[https://dx.doi.org/10.17504/p
rotocols.io.8epv5x11ng1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5x11ng1b/v1)

License: This is an open
access protocol distributed
under the terms of
the [Creative Commons
Attribution License](#), which
permits unrestricted use,
distribution, and reproduction
in any medium, provided the
original author and source
are credited

Protocol status: Working
We use this protocol and it's
working

Created: Jun 28, 2023

Last Modified: Sep 23,
2023

PROTOCOL integer ID:
84153

Keywords: Plasmid
Construction

Plasmid Construction and Gibson cloning

Elias Adriaenssens¹

¹Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna

ASAP Collaborative Research Network

Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

ABSTRACT

This protocol describes plasmid construction.











ATTACHMENTS

[759-1927.pdf](#)

Plasmid Construction

5h

- 1 Obtain the sequences of all cDNAs by amplifying existing plasmids, HAP1 cDNA, or through gene synthesis (Genscript).

- 2 For insect cell expressions, the sequences are codon optimized and the gene is synthesized (Genscript).
- 3 With the exception of the NAP1-6xAla mutant, which was obtained through gene synthesis (Genscript), all other plasmids were generated by Gibson cloning.
- 4 For Gibson cloning, generate inserts and vector backbones by PCR amplification or excise from agarose gels after restriction enzyme digestion at  37 °C for  02:00:00 . 2h
- 5 Purify the inserts and plasmid backbones with Promega Wizard SV gel and PCR Cleanup System (Promega).
- 6 Mix the purified inserts and backbones in a molar 3:1 ratio, respectively, supplemented by a 2x NEBuilder HiFi DNA assembly enzyme mix (New England Biolabs).

- 7 Incubate Gibson reactions for  01:00:00 at  50 °C and then transform into DH5-alpha competent *E. coli* cells. 1h

- 8 Grow transformed Gibson reactions  Overnight on agar plates containing the appropriate selection marker (ampicillin, kanamycin, or chloramphenicol). 1h

- 9 Pick single colonies, grown  Overnight in liquid cultures, and pellet for DNA plasmid extraction using the GeneJet Plasmid Miniprep kit (Thermo Fisher). 1h

- 10 Submit the purified plasmid DNA for DNA Sanger sequencing (MicroSynth AG) to verify insert sequences by Sanger sequencing.
- 11 Further analyze positive clones by whole plasmid sequencing (Plasmidsaurus).