

FEB 12, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io. 5qpvo3zzzv4o/v1

Protocol Citation: Benoit

Quinquis , Alexandre Famechon , nathalie.galleron, Victoria Meslier, mathieu.almeida 2024. DNA yield estimation and custom DNA quality evaluation before shotgun metagenomic sequencing.

protocols.io

https://dx.doi.org/10.17504/protocols.io.5qpvo3zzzv4o/v1

ONA yield estimation and custom DNA quality evaluation before shotgun metagenomic sequencing

Benoit Quinquis¹, Alexandre Famechon¹, nathalie.galleron¹, Victoria Meslier¹, mathieu.almeida¹

¹Université Paris Saclay, INRAE MetaGenoPolis, 78350 Jouy-en-Josas, France



Victoria Meslier

Université Paris Saclay, INRAE MetaGenoPolis, 78350 Jouy-en-...

DISCLAIMER

NA

ABSTRACT

This protocol describes the procedures used for DNA quantification and a custom DNA quality evaluation before proceeding to shotgun metagenomic sequencing.

ATTACHMENTS

Quant-iT-DNA-BRprotocol.pdf DNA_yield_&_custom_qual ity_evaluation_for_shotgun _metagenomic_sequencin g_MetaGenoPolis.pdf

IMAGE ATTRIBUTION

NA

GUIDELINES

NA

MATERIALS

NA

protocols.io

MANUSCRIPT CITATION:

Meslier, V.; Menozzi, E.; David, A.; Morabito, C.; Lucas Del Pozo, S.; Famechon, A.; North, J.; Quinquis, B.; Koletsi, S.; Macnaughtan, J.; et al. Evaluation of an Adapted Semi-Automated DNA Extraction for Human Salivary Shotgun Metagenomics. Biomolecules 2023, 13, 1505.

https://doi.org/10.3390/biom1310

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's

working

1505

Created: Feb 09, 2024

Last Modified: Feb 12, 2024

PROTOCOL integer ID: 94964

Keywords: ASAPCRN, DNA quantification, DNA quality, shotgun metagenomics

Funders Acknowledgement:

Agence Nationale de la Recherche

Grant ID: ANR-11-DPBS-0001

Aligning Science Across

Parkinson's

Grant ID: ASAP-000420

SAFETY WARNINGS



See safety recommendations from the kits used in this protocol

ETHICS STATEMENT

NA

BEFORE START INSTRUCTIONS

NA

DNA quantification

1 using the DNA-binding fluorescent Quant-it dsDNA BR assay



2 Genomic 50 kb kit by capillary electrophoresis using the Fragment Analyzer, followed by a custom DNA quality analysis