



Feb 08, 2022

🌐 Identification of the different phylogroups from the *Klebsiella pneumoniae* species complex by real-time PCR

Elodie Barbier¹, Carla Rodrigues², Sylvain Brisse², Pascal Piveteau³

¹Université de Bourgogne Franche-Comté - INRAE Dijon - France;

²Institut Pasteur, Université de Paris, Biodiversity and Epidemiology of Bacterial Pathogens, Paris, France;

³INRAE, UR OPAALE, F-35044, Rennes, France

1



dx.doi.org/10.17504/protocols.io.b4s5qwg6

Klebsiella Research and Surveillance

Elodie Barbier

Université de Bourgogne Franche-Comté - INRAE Dijon - France

KpSC phylogroup ID

We developed and optimized a real-time PCR for the identification of the different phylogroups (Kp1 to Kp6) belonging to *Klebsiella pneumoniae* species complex (KpSC) by amplification of five specific gene targets.

DOI

dx.doi.org/10.17504/protocols.io.b4s5qwg6

Elodie Barbier, Carla Rodrigues, Sylvain Brisse, Pascal Piveteau 2022. Identification of the different phylogroups from the *Klebsiella pneumoniae* species complex by real-time PCR. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4s5qwg6>



Klebsiella ; phylogroups; PCR ; identification

protocol ,

Feb 08, 2022

Feb 08, 2022

57917

Primers and Probes

1

The **KpSC phylogroup ID** : Primers, Probes and Amplicon length

A	B	C	D	E	F	G
PhG -Target	Forward (5' – 3')	Reverse (5' – 3')	Probe	Quencher	Reporter	Amplicon length (bp)
Kp1 - group_8970	CCTCAAACACGCCAATATGC	TACCGCGACGAGTAAAGTGG	GATCCATTGATTCCATTGGAACCGG	BHQ-1	JOE	85
Kp2 - group_10724	TTGTTGATTGGCAGGCCTTC	AAATGCTGTGACCACCGTTG	CGCCTGGCGTCAGTGGCCCCGAGC	BHQ-1	JOE	70
Kp3/Kp5 - <i>tetR</i>	GCGGGCCGGCTTTTC	CGCATCCCAGGGTATATTCG	AGGCCGAAGCCGCAATGATATTACTCA	BHQ-1	FAM	71
Kp4 - <i>nanK</i>	AAACAGCGGATGCTCTGAAC	ATGAGTCCGCCAAGATTCTG	ATAGCCTCGACCGGGATTATTCGTGA	BHQ-2	TAMRA	145
Kp6 - <i>mdtG_1</i>	GCGTTTCAATGATCCTGTCC	ATCTGCGTCTGGAAAACAG	GTGCGCTGATTGGTATTGGATATGGC	BHQ-2	TAMRA	76

PCR mix preparation

2 PCR mix preparation for each set of primers and probe

A	B
Mix Reagents	Volume per well (μL)
Takyon ROX Probe	12.5
MasterMix dTTP blue	
Forward (500 nM)	1
Reverse (500 nM)	1
Probe (250 nM)	0.5
DNA template	2.5
PCR grade water	2.5

Master Mix: Takyon™ ROX Probe MasterMix dTTP blue(Ref.UF-RPMT-B0701, Eurogentec, Belgium)

Amplification program

3

A	B	C
	Temperature (°C)	Time (minutes:seconds)
Holding stage	95.0	03:00
Cycling stage (x40)	95.0	00:15
60.0	01:00	

Thermocycler: ABI StepOne Real-Time PCR system (Thermo Fischer Scientific)

Positive and negative controls

4

Negative control: PCR grade water

Positive controls:

- Kp1 DNA (*Klebsiella pneumoniae* ATCC 13883^T) at concentration 2 ng/μl ; 2.5 μl used *i.e.*, 5 ng.
- Kp2 DNA (*Klebsiella quasipneumoniae* subsp. *quasipneumoniae* 01A030^T)
- Kp3 or Kp5 DNA (*Klebsiella variicola* subsp. *variicola* F2R9^T; *Klebsiella variicola* subsp. *tropica* 1266^T)
- Kp4 DNA (*Klebsiella pneumoniae* 07A044^T)
- Kp6 DNA (*Klebsiella quasivariicola* KPN1705^T)