



FEB 14, 2024

## 🌐 Open field test in rats

eduard.bentea<sup>1</sup>, María Sanchiz Calvo<sup>1</sup>, Veerle Baekelandt<sup>1</sup>

<sup>1</sup>Katholieke Universiteit Leuven



eduard.bentea

### ABSTRACT

Protocol for performing the open field test in rats. The open field test measures spontaneous locomotion during 5-min free exploration of an open field arena (1m x 1m), and can be used to measure hypolocomotion in animal models of Parkinson's disease. Rats with decreased nigrostriatal dopaminergic neurotransmission show decreased mobility during the open field test.

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.kxygx3dozg8j/v1](https://dx.doi.org/10.17504/protocols.io.kxygx3dozg8j/v1)

**Protocol Citation:** eduard.bentea, María Sanchiz Calvo, Veerle Baekelandt 2024. Open field test in rats. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kxygx3dozg8j/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working



**Created:** Feb 06, 2024

Last Modified: Feb 14, 2024

PROTOCOL integer ID: 94756

**Keywords:** ASAPCRN, motor function, behavioral test, rats

**Funders Acknowledgement:**  
ASAP (Aligning Science Across Parkinson's)

Test		1h 5m
1	Bring cages to the behavioral room for at least  01:00:00 prior to the test for habituation	1h
2	Place rat individually in a clean open field arena (1m x 1m)	
3	Video tape from above for  00:05:00 and record parameters of spontaneous activity (distance traveled, movement time, immobility) using an activity tracking software.	5m
4	Clean arena with mild soap solution	