



VERSION 2
DEC 05, 2022

WORKS FOR ME

1

Detection of phospho IRF3 and p65 in PBMCs by flow cytometry V.2

COMMENTS 0

DOI

dx.doi.org/10.17504/protocols.io.e6nvwdwzlmk/v2

[Andrea Bausch](#)¹

¹In Situ Therapy



Andrea Bausch

ABSTRACT

This protocol describes a method to detect phosphorylated IRF3 and phosphorylated p65 (subunit of NFκB) in PBMCs by flow cytometry. An extracellular staining to identify different immune subsets is followed by an intracellular staining steps for the detection of phospho proteins

DOI

dx.doi.org/10.17504/protocols.io.e6nvwdwzlmk/v2

PROTOCOL CITATION

Andrea Bausch 2022. Detection of phospho IRF3 and p65 in PBMCs by flow cytometry. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwdwzlmk/v2>
Version created by [Andrea Bausch](#)



KEYWORDS

pIRF3, phosphoFlow, pP65, IRF3, p65

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 05, 2022

LAST MODIFIED

Dec 05, 2022

PROTOCOL INTEGER ID

73541

Before start

1 Prepare DC Medium and pre-warm to 37 °C :

⊗ RPMI 1640 Medium, GlutaMAX™ Supplement **Thermo Fisher Catalog #61870036**

+ 5% Human Serum

+ 1% ⊗ MEM Non-Essential Amino Acids Solution (100X) **Thermo Fisher Catalog #11140050**

+ 1% ⊗ Sodium Pyruvate (100 mM) **Thermo Fisher Scientific Catalog #11360070**

Prepare Flow Buffer and cool to 4 °C :

⊗ DPBS no calcium no magnesium **Gibco - Thermo Fisher Catalog #14190169**

+2% heat inactivated FBS

+ 2mM

⊗ Ethylenediaminetetraacetic acid disodium salt solution BioUltra for molecular biology pH 8.0 ~0. **Merck Millipore Sigma Catalog #03690**

Put an aliquot of Methanol in the freezer and an aliquot of

⊗ ROTI®Histofix 4 % 500ml **Carl Roth Catalog #P087.4** and PBS in the fridge

Cool centrifuge to 4 °C

2h

Preparation and stimulation of cells

2 Thaw PBMCs and adjust concentration to 1.11×10^7 cells/ml in DC Medium

3 Seed 90 µL cell suspension in 96-well round bottom plate (1×10^6 cells per well).

4 Prepare a 10x stock solution of your stimulus.

Note

A 1 µM ⊗ diABZI STING agonist-1 **MedChemExpress Catalog #HY-112921A** solution (final concentration in well) can be used as positive control

5 Add 10µl of your prepared stimulus to your cells

- 6 Incubate at 37°C, 5%CO₂ for the desired time



Note

Maximum phosphorylation was detected after 1h for pIRF3 and after 30min for pP65 when stimulated with STING agonists

1h

Extracellular Staining

- 7 All following steps have to be performed On ice using ice cold buffers. All centrifugation steps need to be performed at 4 °C



Put an Aliquot of PBS, Flow Buffer and HistoFix on ice

- 8 Spin cells down 300 x g, 4°C, 00:05:00

Discard supernatant

Resuspend in 200µl Flow Buffer

Spin cells down 300 x g, 4°C, 00:05:00

Discard supernatant

10m

- 9 Prepare Human BD Fc Block™ **Becton-Dickinson Catalog #564220** and staining mixes. Keep in the fridge until use.

- 9.1 Fc Block: Dilute 1:20 in Flow Buffer. Volume needed: 50 µL per sample

- 9.2 Master Mix:
Prepare in BrilliantStainBuffer:




Master Mix for 1 sample:

A	B	C	D	E
Fluoroch	Markers	Dilut	µl/test	µl Ab
BV421	CD3	50	1	1
BV510	CD14	20	2.5	2.5
BV605	HLA-DR	50	1	1



A	B	C	D	E
BV711	CD56	50	1	1
BV785	CD123	50	1	1
FITC	CD19	5	10	10
AF647	CD11c	100	0.5	0.5
eFluor78	Via eF78	500	0.1	0.1
			µl Brilliant	32.9

9.3 Prepare intracellular staining mix
Dilute pIRF3 or pP65 antibody 1:50 in Flow Buffer
Volume needed: 50µl per sample

9.4 Prepare staining mixes for controls
e.g. Viability Dye Single Staining, FMO controls










10 Resuspend in  50 µL FcBlock (1:20 pre-dilution in Flow Buffer) 15m
Incubate  00:15:00  On ice

11  300 x g, 4°C, 00:05:00 5m
Discard supernatant




12 Resuspend cells in 50µl MasterMix 20m
Incubate for  00:20:00  On ice in the dark (cover with aluminum foil to protect from light)

13 Add 150µl of Flow Buffer per well 10m
Spin cells down  300 x g, 4°C, 00:05:00
Discard supernatant
Resuspend in 200µl Flow Buffer
Spin cells down  300 x g, 4°C, 00:05:00
Discard supernatant

Fixation and permeabilization

- 14 Resuspend cells in 100µl of FlowBuffer
Add 100µl of ice cold 4% HistoFix
Incubate  On ice for  00:10:00 in the dark 10m
- 15 Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant
Resuspend in 200µl ice cold PBS
Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant 10m
- 16 Resuspend cells slowly in 180µl ice cold Methanol
Incubate for  00:30:00 - 60min at  4 °C in the fridge 30m
- 17 Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant
Resuspend in 200µl ice cold FACS Buffer
Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant
Resuspend in 200µl ice cold FACS Buffer
Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant 15m

Intracellular Staining

- 18 Resuspend cells in 50µl intracellular staining mix
Incubate  00:30:00 at 4C 40m
- 19 Add 150µl of Flow Buffer per well
Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant
Resuspend in 200µl Flow Buffer
Spin cells down  500 x g, 4°C, 00:05:00
Discard supernatant 10m
- 20 Resuspend in 100µl of Flow Buffer and keep at 4C protected from light until acquisition