



Mar 10, 2022

Environmental DNA (eDNA) metabarcoding protocol for fish species

Omneya Ahmed¹, Tomas Larsson¹, Mats Töpel¹, Alexander Eiler¹¹eDNA solutions AB

1

dx.doi.org/10.17504/protocols.io.b4awqsfe**eDNAsolutions**Tech. support phone: +0700264843 email: omneya@ednasolutions.se

Omneya Ahmed Osman

Environmental DNA metabarcoding universal primers targeting the hypervariable region of the 12S rRNA gene

DOI

dx.doi.org/10.17504/protocols.io.b4awqsfe

Omneya Ahmed, Tomas Larsson, Mats Töpel, Alexander Eiler 2022. Environmental DNA (eDNA) metabarcoding protocol for fish species. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4awqsfe>



protocol ,

Jan 25, 2022

Mar 10, 2022

57398

Serial dilutions of mock community was prepared as a positive control

MATERIALS

 Agencourt Ampure XP **Beckman**

Coulter Catalog #A63AA0

 UltraPure™ DNase/RNase-Free Distilled Water **Thermo Fisher**

Scientific Catalog #10977015

 10 mM dNTPs **Life**

Technologies Catalog #10297-018

 Q5 High-Fidelity DNA Polymerase - 500 units **New England**

Biolabs Catalog #M0491L

The 1st part of the protocol is performed in the pre-PCR room.

The 2nd part in the post-PCR room.

Never bring back PCR products to the pre-PCR room.

Always add a negative control samples in each PCR run

Laboratory work space and equipment were sterilized by UV-light and DNase solution and 70% ethanol. Filter pipet tips were used in all steps of the laboratory work.

- 1 DNA extraction can be performed using Qiagen DNeasy power water sterivex kit. The quality of the extracted DNA was estimated using Nanodrop.

Qiagen DNeasy power water sterivex kit: <https://www.qiagen.com/se/resources/resourcedetail?id=c5fe7d5f-070a-4ebe-ac04-4bbf05a13e91&lang=en>

- 2 Perform the first PCR (triplicates/duplicates of each sample) using Illumina adaptor attached primers that target the gene of your choice.

2.1 MiFish primers

A modified version of the universal primers targeting the hypervariable region of the 12S rRNA gene (163-185 bp) (Miya et al., 2015) was used. The sequence of the primer set is

MiFish-UF: 5'ACA CTC TTT CCC TAC ACG ACG CTC TTC CGA TCT NNN NNN
GTC GGT AAA ACT CGT GCC AGC

MiFish-UmR: 5'AGA CGT GTG CTC TTC CGA TCT NNN NNN CAT AGT GGG GTA
TCT AAT CCC AGT TTG.

2.2 Mock community

DNA extract of 10 fish species were pooled and used as a positive control

-The fish species are *Clupea harengus*, *Glyptocephalus cynoglossus*, *Scomber scombus*, *Thunnus alalunga*, *Pleuroneates platessa*, *Pollachius virens*, *Salmo salar*, *Gadus morhua*, *Reinhardtius hippoglossoides* and *Melanogrammus aeglefinus*.

First PCR reaction

3

First PCR reaction for fish amplification

3.1

Components	Working conc.	Final conc.	1 reaction (μl)
5xQ5 Reaction Buffer	5X	1X	5
MiFish_F	10 μM	0,3 μM	0,75
MiFish_R	10 μM	0,3 μM	0,75
dNTPs	2 mM	0,2 mM	2,5
Q5 HF DNA polymerase	2 U/μl	0.02 U/μl	0,25
Template DNA			5
Nuclease-Free water			10,75
Total			25

For environmental sample add 5 μ l and for mock community add 1 μ l as a template.

STEP	TEMP.	TIME
Initial Denaturation	98 C	30 sec
	98 C	10 sec
15 cycles	66 C	30 sec
	72 C	30 sec
Final Extension	72 C	2 min
Hold	6 C	∞

7 Check second PCR products with Agarose gel electrophoresis (1%)

8 Perform purification with magnetic beads (Agencourt AMPure)
https://research.fhcrc.org/content/dam/stripe/hahn/methods/mol_biol/Agencourt%20AMPure%20XP.pdf

9 Quantification of the concentration of second PCR product before pooling using PicoGreen assay
<http://tools.thermofisher.com/content/sfs/manuals/PicoGreen-dsDNA-protocol.pdf>

Calculate PCR samples concentration and volume before pooling

10 Pool the PCR samples in equal DNA amount (ng) or for unequal length amplicons, in equal molecule amount (mol). You will get one tube including a mix of all the samples.

To calculate the volume of each sample to be pooled (DNA amount mixing):

- Use the lowest concentration sample to define the minimum amount of DNA (ng) that you have available from a single sample:

DNA concentration (ng/μL) of the lowest concentration sample multiplied with its volume (μL).



TH (b)