



DEC 01, 2023

Flow cytometry


Pranay Srivastava¹

¹Massachusetts General Hospital, ASAP

ASAP Collaborative Research Network

Team Chen

[1 more workspace](#) ↓



Fang Zhou

Massachusetts General Hospital, Aligning Science Across Park...

ABSTRACT

This protocol is to assess immune cell profile in spleen (Chauhan et al. 2018).

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.ewov1qm47gr2/v1

Protocol Citation: Pranay Srivastava 2023. Flow cytometry. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1qm47gr2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 30, 2023

Last Modified: Dec 01, 2023

Keywords: ASAPCRN

- 1** Spleen was removed in a 35 mm petri plate with 5 ml RPMI 1640 and digested mechanically and passed through 70 µm filter screen.
- 2** The cell suspension was centrifuged, and the pellet was incubated in RBC lysis buffer. The resulting cell suspension was washed in 1xPBS and blocked with Fc Block (Biolegend, Cat# 101302, 1 µl/50 µl).
- 3** The cells were incubated with MC1R (Invitrogen, Cat# PIPA521911, 1.39 µg) antibody followed by fluorophore conjugated primary antibodies for extracellular markers (Biolegend Cat# 101235, CD11b-BV421 (0.25 µg); Cat# 127641, Ly6G-BV-650 (0.25 µg); Cat# 128041, Ly6C-BV785 (0.125 µg); Cat# 100516, CD4-APC (0.25 µg); Cat# 100751, CD8a-BV510 (0.5 µg); Cat# 152405, CD19-PerCP-Cy5.5 (0.25 µg); Cat# 102036, CD25-BV605 (0.3 µg)) and AF488 (Invitrogen, Cat# A11034, 1:200).
- 4** Zombie dye (Biolegend, Cat# 423101, 1µl/sample) was used to differentiate between live and dead cells. Helper T cells and cytotoxic T cells were identified by CD4+ and CD8+, respectively. CD4+CD25+ cells were used to mark Tregs. CD19+ cells were used as marker for B cells. Monocytes were identified as CD11b+Ly6G-Ly6C^{high} cells and neutrophils were marked by CD11b+Ly6C-Ly6G+. For monocytes and neutrophils, CD11b-positive cells were first extracted from the live cell subset by expansion with SSC-A, followed by Ly6C. After expansion with Ly6C and Ly6G, we excluded Ly6G-positive cells. Monocytes were identified as CD11b+Ly6G-Ly6C high cells, and neutrophils were marked by CD11b+Ly6C-Ly6G+.
- 5** SORP 5 Laser BD Fortessa X-20 (BD Bioscience) and FlowJo v10.7.1 (Becton Dickson & Company) software was used for data acquisition.