



NOV 12, 2023

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1

Protocol Citation: Leonardo A Parra-Rivas 2023. Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn3nrml5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Nov 12, 2023

Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS)

Leonardo A Parra-Rivas¹

¹University of California, San Diego

Subhojit Roy Lab



Leonardo A Parra-Rivas

University of California, San Diego

ABSTRACT

Liquid Chromatography with Tandem Mass Spectrometry (LC-MS/MS)

- 1 The data acquisition by LC-MS/MS was adapted from a published procedure

CITATION

Seyfried NT, Dammer EB, Swarup V, Nandakumar D, Duong DM, Yin L, Deng Q, Nguyen T, Hales CM, Wingo T, Glass J, Gearing M, Thambisetty M, Troncoso JC, Geschwind DH, Lah JJ, Levey AI (2017). A Multi-network Approach Identifies Protein-Specific Co-expression in Asymptomatic and Symptomatic Alzheimer's Disease..

LINK

<https://doi.org/10.1016/j.cels.2016.11.006>

- 2 Derived peptides were resuspended in a loading buffer (0.1% trifluoroacetic acid, TFA) and separated on a Waters' Charged Surface Hybrid (CSH) column (150 µm internal diameter x 15 cm; particle size: 1.7 µm).
- 3 The samples were run on an EVOSEP liquid chromatography system using the 15 samples per day preset gradient (88 min) and were monitored on a Q-Exactive Plus Hybrid Quadrupole-Orbitrap Mass Spectrometer (Thermo Scientific).
- 4 The mass spectrometer cycle was programmed to collect one full MS scan followed by 20 data-dependent MS/MS scans. The MS scans (400-1600 m/z range, 3×10^6 AGC target, 100 ms maximum ion time) were collected at a resolution of 70,000 at m/z 200 in profile mode.
- 5 The HCD MS/MS spectra (1.6 m/z isolation width, 28% collision energy, 1×10^5 AGC target, 100 ms maximum ion time) were acquired at a resolution of 17,500 at m/z 200. Dynamic exclusion was set to exclude previously sequenced precursor ions for 30 seconds. Precursor ions with

+1, and +7, +8, or higher charge states were excluded from sequencing.