



Jan 28, 2022

# RNA extraction from colonial tunicates

Marta Wawrzyniak<sup>1</sup>, Simon Blanchoud<sup>1</sup><sup>1</sup>University of Fribourg

1

[dx.doi.org/10.17504/protocols.io.b33nqqme](https://dx.doi.org/10.17504/protocols.io.b33nqqme)

Blanchoud lab, UNIFR

Marta Wawrzyniak  
University of Fribourg

This protocol has been successfully used with *Botrylloides diegensis* and has been adapted from the following publication:

[An efficient low-cost laboratory workflow for the study of blood cells and RNA extractions in marine invertebrates](#)

DOI

[dx.doi.org/10.17504/protocols.io.b33nqqme](https://dx.doi.org/10.17504/protocols.io.b33nqqme)

Marta Wawrzyniak, Simon Blanchoud 2022. RNA extraction from colonial tunicates. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.b33nqqme>



ascidians, colonial tunicates, RNA extraction

protocol ,

Jan 20, 2022


Jan 28, 2022

57166

Change gloves frequently, particularly as the protocol progresses from crude extracts to more purified materials. Use sterile tubes. Perform all steps on ice and use RNase-free water unless otherwise stated.

Liquid nitrogen  
Sterile tubes and plastic pestles  
Extraction buffer : 0.2M Tris-HCl pH 7.5 , 0.1M LiCl, 5mM EDTA, 1/10 of the total volume of SDS 10%  
Phenol pH 4 (4 C)  
Chloroform  
LiCl (for 50mL: 12.6g 6M LiCl ; 6.3g 3M LiCl)  
SC-EtOH: Sodium acetate + 100% Ethanol (1/3 : 2/3)  
70% and 100% Ethanol  
RNase-free water

- 1 Clean the slide from which you will take the colony of your interest. See [Cleaning colonial ascidians](#).
- 2 Isolate a cleaned colony composed of approx. 20 zooids.

2.1 Transfer to a tube and spin at maximum speed for  **00:02:00** .

2m

2.2 Remove the excess water and shock-freeze the tube in liquid nitrogen.

3 Add  **400 µL** of extraction buffer to the frozen sample and macerate with a plastic pestle.


4 Add  **100 µL** more of extraction buffer and  **500 µL** of 1:1 phenol:chloroform.

5 Mix the tube by inversion a couple of times until it gets cloudy.

6 Centrifuge the homogenate at 1400 g for  **00:05:00** at  **4 °C** .

5m


7 Carefully collect  **400 µL** of the upper phase into a new tube.

7.1 Note: if desired this sample could be used for DNA extraction - carefully transfer  **200 µL** of the interphase into a new tube (See [DNA extraction from colonial tunicates](#)).

8 Add  **500 µL** of **[M]6 Molarity (M)** LiCl to the supernatant.

9 Incubate the mixture at  **-80 °C** for  **01:00:00** . 1h

10 Centrifuge at 1400 g for  **00:10:00** at  **4 °C** . 10m

11 Discard the supernatant and resuspend the pellet in  **1 mL** of **[M]3 Molarity (M)** LiCl.









12 Shake slowly for  **00:15:00** at  **Room temperature** on a linear shaker. 15m

13 Centrifuge at 1400 g for  **00:10:00** at  **4 °C** . 10m

14 Discard the supernatant and resuspend the pellet in  **1 mL** of SC-EtOH solution.

15 Incubate at  **-80 °C** for  **00:15:00** . 15m

16 Centrifuge at 1400 g for  **00:15:00** at  **4 °C** . 15m

- 17 Discard the supernatant and wash the pellet with  **1 mL** of **70 % volume** Ethanol.
- 18 Centrifuge at 1400 g for  **00:05:00** at  **4 °C** . 5m
- 19 Discard the supernatant and place the tubes up-side-down on a paper towel for  **00:05:00**<sup>15m</sup> to  **00:10:00** .
- 20 Resuspend the pellet in RNase-free water (  **20 µL** to  **100 µL** depends on the amount of pellet).
- 21 Quantify the RNA concentration and quality using the NanoDrop, the capillary electrophoresis and/or the Bioanalyzer.
- 22 Store at  **-80 °C** .