



APR 19, 2023

## The extraction of cross-kingdom biofilm DNA using Zymo Reasearch Quick-DNA Fecal/Soil Microbe Kits.

Lindsay Kate Newbold<sup>1</sup>, Daniel S Read<sup>2</sup>, Joetay<sup>2</sup>,  
Jonathan Warren<sup>3</sup>, kerry.walsh<sup>3</sup>

<sup>1</sup>UK Centre for Ecology & Hydrology;

<sup>2</sup>UK Centre for Ecology & Hydrology (UKCEH); <sup>3</sup>Environment Agency



Lindsay Kate Newbold

UK Centre for Ecology & Hydrology

OPEN  ACCESS

**Protocol Citation:** Lindsay Kate Newbold, Daniel S Read, Joetay, Jonathan Warren, kerry.walsh 2023. The extraction of cross-kingdom biofilm DNA using Zymo Reasearch Quick-DNA Fecal/Soil Microbe Kits.. **protocols.io** <https://protocols.io/view/the-extraction-of-cross-kingdom-biofilm-dna-using-crzm746>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Mar 29, 2023

**Last Modified:** Apr 19, 2023

**PROTOCOL integer ID:**  
79629

### DISCLAIMER

#### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](https://protocols.io) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](https://protocols.io), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## ABSTRACT

The extraction of DNA from environmental samples, including soils, sediments, water, and biofilms, is a crucial stage in the analysis of environmental microbial communities and wider environmental communities via environmental DNA (eDNA). The choice of DNA extraction methodologies significantly influences not only the yield and purity of the extracted DNA, but also the composition of the community as determined by DNA analysis techniques such as amplicon and metagenomic sequencing. Numerous studies have investigated the influence of different kit-based DNA extraction protocols on metabarcoding of various taxa, including bacteria, fungi, diatoms, microeukaryotes, and invertebrates (Dopheide et al. 2019; Giangacomo et al. 2021; Matsuoka et al. 2022; Shaffer et al. 2022; Vasselon et al. 2017; Pearman et al. 2020; Majaneva et al. 2018; Deiner et al. 2015; Djurhuus et al. 2017; Kawato et al. 2021). Here we present an amended version of the manufacturers (Zymo research) recommended protocol for Quick-DNA Fecal/Soil Microbe Kits optimised for the extraction of cross-kingdom biofilm community DNA suitable for downstream applications such as metabarcoding.

## GUIDELINES

This protocol is based upon the recommended manufacturers protocol for Zymo Quick-DNA Fecal/Soil Microbe Miniprep Kit (D6012) and Quick-DNA Fecal/Soil Microbe 96 Kit (D6011). When using the 96 Kit volumes should be adjusted as per manufacturers recommendations.

## SAFETY WARNINGS




Please read material SDS sheet before starting protocol.







\_d6011\_quick-dna\_fecalsoil\_microbe\_96\_kit.pdf

## BEFORE START INSTRUCTIONS

Add  500  $\mu$ L beta-mercaptoethanol per 100ml Genomic Lysis Buffer.

## Sample Preparation

15m

- 1 In field samples are taken by the Environment Agency, and  5 mL biofilm suspension preserved in equal amounts of preservation buffer:  3.5 Molarity (M) ammonium sulphate,  17 millimolar (mM) sodium citrate,  13 millimolar (mM) EDTA. Samples are frozen in

## Equipment

### Centrifuge tube

15 ml Centrifuge tube

Corning

CLS430791-500EA

<https://www.sigmaaldrich.com/GB/en/product/sigma/cls430791>

NAME

TYPE

BRAND

SKU


LINK

prior to laboratory preparation.

- 2 Once back in the laboratory field preserved samples are first thawed and then vortexed to create a homogenous mixture.

- 3 Samples are pelleted at  3000 x g, 5°C for  00:15:00 minutes

15m

- 4 Ensure samples are fully pelleted, then remove  500 µL of buffer to separate microcentrifuge tube.

## Equipment

### 2 ml Eppendorf Safe-Lock tube

microcentrifuge tube

Eppendorf

0030120094

<https://www.eppendorf.com/gb-en/eShop-Products/Laboratory-Consumables/Tubes/Eppendorf-Safe-Lock-Tubes-p-0030120094>

NAME



TYPE

BRAND

SKU












LINK

Remaining buffer is gently aspirated to avoid disrupting pellet.



- 5 Pellet is resuspended in  500  $\mu\text{L}$  of reserved preservation buffer, and stored in the same fully labelled microcentrifuge tube. Preserved sample is then frozen at   $-80^{\circ}\text{C}$  prior to analysis.

## Sample Lysis

1h 11m


- 6 Thaw sample to room temperature, vortex to resuspend pellet and add  100  $\mu\text{L}$  to ZR BashingBead Lysis tube. Add  900  $\mu\text{L}$  DNA/RNA shield, and lyse for  00:40:00 on a horizontal vortex mixer (Vortex Genie). 40m
- 7 Briefly spin for  10000 x g,  $18^{\circ}\text{C}$ , 00:01:00 . Add  40  $\mu\text{L}$  proteinase K, resuspend and mix through pipetting and incubate at   $65^{\circ}\text{C}$  for  00:20:00 . 21m
- 8 Centrifuge ZR BashingBead tube at  10000 x g,  $18^{\circ}\text{C}$ , 00:05:00 to pellet lysate. 5m
- 9 Transfer up to  400  $\mu\text{L}$  supernatant to a Zymo-Spin III-F Filter in a collection tube and centrifuge at  8000 x g,  $18^{\circ}\text{C}$ , 00:01:00 . 1m
- 10 Add  1200  $\mu\text{L}$  of Genomic Lysis buffer in a fume hood to the filtrate in the collection tube from step 4. Mix well.

### Note

Genomic Lysis buffer should have beta-mercaptoethanol added to a final dilution of 0.5% (v/v) ie  500  $\mu\text{L}$  per  100 mL buffer. Beta-mercaptoethanol is a strong reducing agent added to clean tannins and other polyphenolics found in plant extracts. Please refer to material SDS and use appropriate PPE.

## DNA purification

9m 30s

- 11** Transfer  800  $\mu\text{L}$  of the mixture from step 10 to a Zymo-Spin IICR column in a collection tube and centrifuge at  10000 x g, 18°C, 00:01:00 . Discard flow through from the collection tube and repeat with remaining mixture. 1m
- 12** Add  200  $\mu\text{L}$  DNA Pre-Wash buffer to the Zymo-Spin IICR Column in a **new** collection tube and centrifuge  10000 x g, 18°C, 00:01:00 . 1m
- 13** Add  500  $\mu\text{L}$  g-DNA Wash Buffer to the Zymo-Spin IICR Column and centrifuge  10000 x g, 18°C, 00:01:00 . 1m
- 14** Transfer the Zymo-Spin IICR column to a clean  1.5 mL centrifuge tube and add  100  $\mu\text{L}$  \* DNA Elution buffer directly to the column matrix. Centrifuge  10000 x g, 18°C, 00:00:30 to elute the DNA. 30s
- ### Note
- \* If samples exhibit a low biomass it's recommended that this volume be reduced to  50  $\mu\text{L}$  in order to concentrate DNA recovered.
- 15** Place a Zymo-Spin III-HRC Filter in a clean collection tube and add  600  $\mu\text{L}$  prep solution. Centrifuge  8000 x g, 18°C, 00:03:00 . Discard flow through and place in clean labelled 1.5 ml microcentrifuge tube. 3m
- 16** Transfer Eluted DNA from step 14 to a prepared Zymo-Spin III-HRC Filter and centrifuge at exactly  16000 x g, 18°C, 00:03:00 . 3m
- 17** The filtered DNA is now suitable for PCR and other downstream applications.

