



Aug 05, 2021

Nanovesicles extraction

Veerle Baekelandt¹¹KU Leuven

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.bw5wpg7e

Joris Van Asselberghs

ABSTRACT

Nanovesicles extraction (Exosomal isolation, SHSY-5Y)

DOI

dx.doi.org/10.17504/protocols.io.bw5wpg7e

PROTOCOL CITATION

Veerle Baekelandt 2021. Nanovesicles extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bw5wpg7e>



LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 05, 2021

LAST MODIFIED

Aug 05, 2021

PROTOCOL INTEGER ID

52118

- 1 seed cells for 24:00:00 in 3 x 15 cm dishes at a density of 20 x 10⁶ cells per plate 1d
 - 2 wash cells of all the plates with 2 times 10 mL of 1x PBS
 - 3 Add 10 mL DMEM with 1% exosome depleted FBS to each plate for 24:00:00 1d
- Preparation of exosome depleted FBS:
Commercial FBS was filtrated through 0.22 µm filter and 140000 x g , overnight
- 4 collect medium.
purify nanovesicles with differential centrifugation (see substeps)

4.1  **300 x g, 00:15:00** (remove cells) 15m

4.2  **15000 x g, 00:30:00** (remove cell debris) 30m

4.3 filter through 0.22 µm filter (filtration of apoptotic bodies)

4.4  **140000 x g, 4°C, 03:00:00** 3h

5 decant supernatant and collect nanovesicles in  **50 µl** cell lysis buffer (for western blotting) or 50 µl PBS+ PI (for intact cell cross-linking)




Intact cell cross-linking
by Joris Van Asselberghs

PREVIEW

RUN






5.1 seed cells (SHSY-5Y) in 10 cm dishes at a density of 6×10^6 cells per plate.
incubate for  **24:00:00** 1d

5.2 wash cells 2 times with 1X PBS

5.3 harvest cells with 1x PBS + PI cocktail (PBS + PI buffer)






5.4  **1000 x g, 00:15:00** and resuspend in  **400 µl** 15m
 **1X PBS (Phosphate-buffered saline)** **Contributed by users** + PI buffer

5.5 divide  **400 µl** resuspended cells over 2 vials 30m
tube 1: add 1 µM disuccinimidyl glutarate (DSG; Thermo Fisher Scientific)
tube 2: add 1 µM dithiobis succinimidyl propionate (DSP; Thermo Fisher Scientific)
incubate  **00:30:00** at  **37 °C**

Preparation of 50X DSP, 50X DSG

50X DSP: 2.02mg into 100µL DMSO
50X DSG: 1.632mg into 100µL DMSO

Note: Prepare the crosslinkers freshly before use

- 5.6 quench reaction for  **00:15:00** on  **Room temperature** in both tubes with  **20 Milimolar (mM)** Tris (^{15m}
 **7.4**)
- 5.7 sonicate cells 2 times at 30 Hz with 15 ON-OFF intervals of  **00:00:10** (1s pulses) ^{10s}