



VERSION 2

FEB 09, 2023

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.n92ldp2m8l5b/v2

Protocol Citation: Boyang Liu, Liangyu Cui, Zhangwen Deng, Yue Ma, Diancheng Yang, Yanan Gong, Yanchun Xu, Shuhui Yang, Song Huang 2023. DNA extraction protocol for snake genomic sequencing by using AxyPrep™ Multisource Genomic DNA Miniprep Kit. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n92ldp2m8l5b/v2> Version created by Boyang Liu

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Feb 06, 2023

Last Modified: Feb 09, 2023

DNA extraction protocol for snake genomic sequencing by using AxyPrep™ Multisource Genomic DNA Miniprep Kit V.2

In 2 collections

Boyang Liu¹, Liangyu Cui¹, Zhangwen Deng², Yue Ma¹, Diancheng Yang^{3,4}, Yanan Gong^{3,4}, Yanchun Xu¹, Shuhui Yang¹, Song Huang^{3,4}

¹College of Wildlife and Protected Area, Northeast Forestry University, Harbin 150040, China;

²Guangxi Forest Inventory and Planning Institute, Nanning 530011, China;

³Anhui Province Key Laboratory of the Conservation and Exploitation of Biological Resource, College of Life Sciences, Anhui Normal University, Wuhu 241000, China;

⁴Huangshan Noah Biodiversity Institute, Huangshan 245000, China















Boyang Liu

ABSTRACT

Genomic DNA was extracted from the muscle of a snake using the AxyPrep™ Multisource Genomic DNA Miniprep Kit (Axygen, China), with some modifications from the kit instructions as outlined in this protocol. The extracted DNA was then suitable for WGS sequencing.

Keywords: DNA extraction,
AxyPrep™ Multisource
Genomic DNA Miniprep Kit,
nucleic acid isolation

- 1 The following protocol is a modified DNA extraction protocol using the AxyPrep™ Multisource Genomic DNA Miniprep Kit.

Approximately 10mg of snake tissue was minced into small pieces using sterilized scissors on ice.
 On ice
- 2 Add 350 µl Buffer PBS and 15 µl SDS Lysis Buffer and shake for 30s. 
 00:00:30
- 3 Add 150 µl Buffer C-L and 20 µl Proteinase K. Immediately vortex for 1 min to mix well. After brief centrifugation, incubate at 56°C for 10 min. 
 00:01:00
 56 °C  00:10:00
- 4 Add 350 µl Buffer P-D, vortex for 30s to mix well, and centrifuge at 12,000×g for 10 min. 
 00:00:30
 12000 rpm, 00:10:00
- 5 Place the DNA preparation tube in a 2 ml centrifuge tube, transfer the mixture in step 4 to the preparation tube, and centrifuge at 12,000×g for 1 min.
 12000 rpm, 00:01:00
- 6 Discard the filtrate, put the preparation tube back into the original 2 ml centrifuge tube, add 500µl Buffer W1, and centrifuge at 12,000×g for 1 min.
 12000 rpm, 00:01:00
- 7 Discard the filtrate, put the preparation tube back into the original 2 ml centrifuge tube, add 700µl Buffer W2, centrifuge at 12,000×g for 1 min.

🔄 12000 rpm, 00:01:00

7.1 In the same way, cleanse again with 700µl Buffer W2.

Note

1. Confirm that absolute ethanol has been added to Buffer W2 concentrate according to the volume specified on the reagent bottle.
2. Cleaning again with Buffer W2 can ensure that the salt is completely removed and eliminates the impact on the enzyme digestion reaction.

8 Discard the filtrate, put the preparation tube back into the original 2ml centrifuge tube, and centrifuge at 12,000×g for 1 min.

🔄 12000 rpm, 00:01:00

9 Place the DNA preparation tube in another clean 1.5 ml centrifuge tube, add 100-200µl Eluent or deionized water to the center of the preparation tube membrane, stand at room temperature for 1 min, and centrifuge at 12,000×g for 1 min to elute the DNA.

1m

🕒 00:01:00 | 🌡 Room temperature

🔄 12000 rpm, 00:01:00