



APR 24, 2023

# Promega Wizard DNA extraction - Drosophila whole body protocol

finley.thomas.<sup>1</sup>

<sup>1</sup>UCL



finley.thomas.

## ABSTRACT

## OPEN ACCESS

**Protocol Citation:** finley.thomas. 2023. Promega Wizard DNA extraction - Drosophila whole body protocol. **protocols.io** <https://protocols.io/view/promega-wizard-dna-extraction-drosophila-whole-body-cs8wwhxe>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Apr 21, 2023

**Last Modified:** Apr 24, 2023

**PROTOCOL integer ID:**  
80886

**Keywords:** Promega, DNA extraction, Wizard, Drosophila, Insect

Edited version of the 'Promega Wizard Genomic DNA Purification' protocol optimised to work on pooled whole-body fruit fly samples.

## Lysis & Incubation

1 Add 10 flies into a 1.5ml Eppendorf tube and cover with 300 µL Nuclei Lysis Solution .

2m

- 2 Use a **micro-tube pestle to crush** the flies until they are fully ground.


5m



- 2.1 For the initial crushing, use your hands and the pestle only.


*Once the flies are partially disintegrated an electronic homogeniser can be used to increase efficiency.*

Take care not to froth the solution up too much.


- 3 Once the flies are thoroughly homogenised, add another  300 µL of Nuclei Lysis Solution .


2m

*By this point, the solution should be opaque and red-brown in colour.*

- 4 **Vortex** the homogenised solution then **incubate** at  65 °C for 20 minutes .


20m

- 5 Remove the samples from the heat and lower the set temperature to  55 °C .

- 6 Allow for the **samples to cool** for a few moments, then add  17.5 µL Proteinase K .

2m

(This is not included in the Promega Wizard kit).





- 7 **Vortex the samples** thoroughly then **incubate** at  55 °C for 3 hours .

3h

Vortex the samples regularly throughout this time.

## Protein Precipitation

4m

- 8 Remove the samples from the heat block and **allow to cool** to room temperature.
- 9 Add  200 µL Protein Precipitation Solution to each sample & **vortex** for **20 seconds** (until thoroughly mixed).
- 10 Leave the samples  On ice for 5 minutes . 5m
- 11 **Centrifuge** the samples at high speed:  13200 rpm , for 4 minutes . 4m
- 12 **Carefully transfer the supernatant** (liquid) to a new (labelled) micro-centrifuge/Eppendorf tube. 2m
- Take care not to disturb the pellet.**
- 13 *If any tissue or protein precipitate (white mass) remains in the supernatant, **repeat steps 11 & 12.***
- 14 Once the supernatant is clear of tissue mass and protein precipitate, add  600 µL isopropanol (at room temperature). 2m

## DNA elution and cleansing

31m

- 15 **Gently invert** the samples to mix the isopropanol and supernatant. 5m
- White threads of DNA may or may-not form; if no threads are visible after ~5 minutes, continue on with the protocol.*

- 16 **Centrifuge** the samples at  13200 rpm , for 1 minute . 1m
- 17 Taking care not to disturb the pellet (of DNA), **remove the supernatant**. 4m  
*If the supernatant does not contain the pellet, it can be discarded.*
- 18 Add  600 µL 70% ethanol , then **gently invert** the tube several times (to wash the DNA). 3m
- 19 **Centrifuge** at  13200 rpm , for 1 minute . 1m
- 20 Making sure that the DNA pellet is not disturbed, **remove the ethanol** (this can be discarded). 2m
- 21 **Invert the (open) tube** onto clean absorbent paper; **leave for 10-15 minutes**. 15m  
*Towards the end of this time, set the heating bloc to  65 °C .*
- 22 Add  100 µL DNA Rehydration Solution **(Nuclease Free Water)** and **incubate** at  65 °C for 1 hour .  
*Mix the solution by gently tapping and shaking the tubes.*
- 23 Store the DNA at  2-8 °C .