



Aug 24, 2020

nCoV-2019 McGill RT Protocol, Lunascript

Forked from [Artic nCoV-2019 McGill modified Lunascript Reverse Transcriptase Nanopore sequencing protocol](#)Shu-Huang Chen¹, Sarah J Reiling¹, Josh Quick²¹McGill University; ²University of Birmingham**1** Works for me dx.doi.org/10.17504/protocols.io.bjgekjte

McGill Genome Centre

Sarah Reiling
McGill University

ABSTRACT

Artic nCoV-2019 McGill modified Lunascript Reverse Transcriptase sequencing protocol.

DOI

dx.doi.org/10.17504/protocols.io.bjgekjte

PROTOCOL CITATION

Shu-Huang Chen, Sarah J Reiling, Josh Quick 2020. nCoV-2019 McGill RT Protocol, Lunascript.

protocols.io<https://dx.doi.org/10.17504/protocols.io.bjgekjte>

FORK FROM

Forked from [Artic nCoV-2019 McGill modified Lunascript Reverse Transcriptase Nanopore sequencing protocol](#),
[Shu-Huang Chen](#)

KEYWORDS

Cov-19 Nanopore Sequencing, Cov-19 Illumina Sequencing

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 07, 2020

LAST MODIFIED

Aug 24, 2020

PROTOCOL INTEGER ID

40166

MATERIALS

NAME	CATALOG #	VENDOR
Water, nuclease free		
LunaScript RT SuperMix Kit	E3010L	New England Biolabs

cDNA preparation

1

Mix the following components in a 0.2 mL 8-strip tube;

Component	Volume
LunaScript RT SuperMix (5x)	4 µl
Nuclease-free water	5 µl
Template RNA	11 µl
Total	20 µl



Viral RNA input from a clinical sample should be between Ct 18-35. If Ct is between 12-15, then dilute the sample 100-fold in water, if between 15-18 then dilute 10-fold in water. This will reduce the likelihood of PCR-inhibition.



A mastermix should be made up in the **mastermix cabinet**. Tubes should be wiped down when entering and leaving the mastermix cabinet.

2

Gently mix by pipetting and pulse spin the tube to collect liquid at the bottom of the tube.

3

Incubate the reaction as follows:

⬆ 25 °C for ⌚ 00:02:00

⬆ 55 °C for ⌚ 00:20:00

⬆ 95 °C for ⌚ 00:01:00

Place on ice for ⌚ 00:01:00 or store cDNA at -20C