



Sep 26, 2022

Preparation of 0.5L of phosphate buffer (pH=6.0)

Alfonso Pérez Escudero¹, Alid Al-Asmar¹, gabrielmadirolas¹¹Research Centre on Animal Cognition (CRCA), Centre for Integrative Biology (CBI), Toulouse III University, CNRS, Toulouse, France

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.n2bvj8r7bgk5/v1

Perez-Escudero Lab



theapelab

ABSTRACT

Original source of the protocol: WormBook Methods

http://www.wormbook.org/chapters/www_strainmaintain/strainmaintain.html

This protocol is for making pH=6.0 phosphate buffer, which is used in the preparation of worm media, like NGM (Nematode Growth Medium).

DOI

dx.doi.org/10.17504/protocols.io.n2bvj8r7bgk5/v1

PROTOCOL CITATION

Alfonso Pérez Escudero, Alid Al-Asmar, gabrielmadirolas 2022. Preparation of 0.5L of phosphate buffer (pH=6.0). **protocols.io**
<https://protocols.io/view/preparation-of-0-5l-of-phosphate-buffer-ph-6-0-cgfttwjn>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 19, 2022

LAST MODIFIED

Sep 26, 2022

MATERIALS TEXT

Monopotassium phosphate:

 Potassium phosphate monobasic Sigma

Aldrich Catalog #795488-500G


Dipotassium phosphate:

 Potassium phosphate dibasic Sigma


Aldrich Catalog #P3786-500G

Filter: Rapid-Flow Nalgene 0.2µm aPES membrane, ref: 564-0020


Vaccum pump: Fisher Brand, ref: FB70155

- 1 Add  **54.2 g ±0.5** of monopotassium phosphate¹ (KH₂PO₄) to a clean 0.5L bottle.

¹Also referred to as potassium phosphate monobasic

- 2 Add  **17.8 g ±0.2** of dipotassium phosphate² (K₂HPO₄) to the bottle.

²Also referred to as potassium phosphate dibasic

- 3 Measure  **500 mL** of milliQ water in a measuring cylinder, and add it to the bottle.

- 4 Shake the bottle vigorously until all crystals are dissolved. This can take a few minutes.

- 5 If this is the first time you prepare this solution, you can ensure that the pH is around 6 using pH paper or a pH-meter. We usually have a pH between 5.8 and 6.

- 6 Filter-sterilize.

- 7 Store on shelf.
We usually use it without problem for at least a year after preparation.