



Jan 01, 2021

Chemically competent cells transformation

Elizabeth Fozo¹¹In-house protocols

1

Works for me

This protocol is published without a DOI.



Eadewunm

ABSTRACT

Chemically competent cells transformation protocol

PROTOCOL CITATION

Elizabeth Fozo 2021. Chemically competent cells transformation. **protocols.io**
<https://protocols.io/view/chemically-competent-cells-transformation-bq4xmyxn>

KEYWORDS

chemically competent cells, cell transformation, cells, protocol

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 01, 2021

LAST MODIFIED

Jan 01, 2021

PROTOCOL INTEGER ID

45943

GUIDELINES

1. Preparation of competent cells
2. Transformation

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Chemically competent cells transformation protocol

Preparation for competent cells

- 1 Inoculate a single colony into 5 mL appropriate media containing appropriate antibiotics

- 2 Incubate overnight using appropriate temperature and shaking conditions
- 3 Dilute overnight culture to an optical density of 0.01 in 10 mL of appropriate media containing appropriate antibiotics
- 4 Incubate at appropriate temperature and shaking conditions until optical density reaches 0.4 to 0.6. Do not allow to grow past this point.
- 5 Spin down cells at 3500 RPM at 4°C for 10 min
- 6 Resuspend pellet in 10 mL ice-cold 0.1M CaCl₂
- 7 Place cells on ice for 20 min (do not go under this time; can go longer)
- 8 Spin down cells as in step 5.
- 9 Resuspend pellet in 1 mL ice-cold 0.1M CaCl₂ + 13% glycerol (for a higher concentration of cells, use less)
- 10 Place on ice for 20 min
- 11 Use immediately or store at -80°C (0.25 mL/tube)

Transformation

- 12 Add 50 µL competent cells into a sterile Eppendorf tube
- 13 Add 1 µL plasmid DNA or ligation reaction (can vary this)
- 14 Place on ice for 30 min

- 15 Heat shock at 42°C for 1 min
- 16 Place on ice 2 min
- 17 Add 0.5 mL recovery media (SOC, LB, etc.)
- 18 Incubate for 1 hour at appropriate growth conditions.
- 19 Plate 200 µL on appropriate plate with appropriate antibiotics.
- 20 Incubate at 37°C overnight