



MAY 08, 2023



**DOI:**  
[dx.doi.org/10.17504/protocols.io.x54v9dm7pg3e/v1](https://dx.doi.org/10.17504/protocols.io.x54v9dm7pg3e/v1)

**Protocol Citation:** Kanae Nishii, Michael Möller, Robert G. Foster, Laura L. Forrest, Nathan Kelso, Sadie Barber, Caroline Howard, Michelle L. Hart 2023. High molecular weight plant DNA extraction for PacBio HiFi sequencing. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v9dm7pg3e/v1>

**MANUSCRIPT CITATION:** Kanae Nishii, Michael Möller, Robert G. Foster, Laura L. Forrest, Nathan Kelso, Sadie Barber, Caroline Howard, Michelle Hart (2023) A high quality, high molecular weight DNA extraction method for PacBio HiFi genome sequencing of recalcitrant plants. *Plant Methods*

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


# High molecular weight plant DNA extraction for PacBio HiFi sequencing

Robert G.  
 Kanae Nishii<sup>1,2</sup>, Michael Möller<sup>1</sup>, Foster<sup>3</sup>,  
 Laura L. Forrest<sup>1</sup>, Sadie  
 Nathan Kelso<sup>1</sup>, Barber<sup>1</sup>, Caroline Howard<sup>4</sup>,  
 Michelle L. Hart<sup>1</sup>

<sup>1</sup>Royal Botanic Garden Edinburgh; <sup>2</sup>Kanagawa University;

<sup>3</sup>Edinburgh Genomics, The University of Edinburgh;

<sup>4</sup>Wellcome Sanger Institute



Kanae Nishii

## ABSTRACT

High molecular weight and high quality DNA is mandatory for successful long-read sequencing. In addition, PacBio HiFi SMRTBell library construction requires omission of traditional DNA extraction chemicals such as guanidinium, chloroform and others. We developed a DNA extraction protocol working well for the recalcitrant plant *Streptocarpus*, and the extracted DNA was successfully used for PacBio HiFi sequencing.

## GUIDELINES

This is a DNA extraction protocol designed for *Streptocarpus* PacBio HiFi sequencing on the Sequel II system.

## MATERIALS

### Chemicals & Reagents

#### ■ Nuclear isolation base (NIB) buffer (final concentrations):

Tris-HCl, pH 8.0, 10 mM

EDTA, pH 8.0, 10 mM

Sucrose, 500 mM

KCl, 100 mM

#### \*To prepare 1L NIB buffer:

Tris-HCl, pH 8.0, 1M stock 10.0 ml

EDTA, pH 8.0, 500 mM stock 20.0 ml

Sucrose 171.2 g

KCl 7.5 g

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 26, 2022

**Last Modified:** May 08, 2023

**PROTOCOL integer ID:**  
70494

**Keywords:** PacBio HiFi, SMRTBell library, high molecular weight DNA extraction, plant genome, long-read sequencing

Add deionized water to make up 1 L.

Mix chemicals in the glass bottle and autoclave. Once opened, store at 4 °C.

■ Nuclear isolation (NI) buffer (final concentrations):

Tris-HCl, pH 8.0, 10 mM

EDTA, pH 8.0, 10 mM

Sucrose, 500 mM

KCl, 100 mM

Spermidine, 4 mM

Spermine, 1 mM

β-mercaptoethanol, 0.1 %

\*To prepare 400 ml NI buffer

NIB buffer	400.0 ml
------------	----------

Spermidine 4M stock	400.0 μl
---------------------	----------

Spermine 1M stock	400.0 μl
-------------------	----------

β-mercaptoethanol	400.0 μl
-------------------	----------

Mix chemicals on the day of DNA extraction, in glass bottle or 500 ml glass beaker, and leave on ice.

■ 10% Triton X-100/NIB:

\*To prepare, dilute 5 ml Triton X-100 with 45 ml NIB buffer in a 50 ml Falcon tube.

Heat in water bath to ca. 30 - 50 °C and once Triton X-100 is dissolved store at 4 °C.

■ Sorbitol buffer base solution (final concentrations):

Tris-HCl, pH 8.0, 100 mM

EDTA, pH 8.0, 5 mM

Sorbitol, 700 mM

\*To prepare 1L sorbitol buffer base solution

Tris-HCl, pH 8.0, 1 M stock	100.0 ml
-----------------------------	----------

EDTA, pH 8.0, 500 mM stock	10.0 ml
----------------------------	---------

Sorbitol	127.5 g
----------	---------

Add deionized water to make up 1 L.

Mix chemicals in glass bottle and autoclave. Once opened, store at 4 °C.

*NOTE:* Original sorbitol buffer base uses 350 mM sorbitol (Souza et al. 2012. Genet Mol Res 11: 756-64), but is here increased to 700 mM.

■ Sorbitol buffer (final concentrations):

Tris-HCl, pH 8.0, 100 mM

EDTA, pH 8.0, 5 mM

Sorbitol, 700 mM

PVP40, 1%

β-mercaptoethanol, 0.2%

\*To prepare 100 ml sorbitol buffer

Sorbitol buffer base solution	100.0 ml
PVP40	1.0 g (ca. half a dispensing spoon)
$\beta$ -mercaptoethanol	200.0 $\mu$ l

Prepare sorbitol buffer on day of DNA extraction. Mix chemicals in two of 50 ml Falcon tubes. Mix well.

- CTAB lysis buffer (final concentrations):

Tris-HCl, pH 8.0, 100 mM  
EDTA, pH 8.0, 20 mM  
NaCl, 1.4 M  
CTAB, 2%

\*To prepare 1L CTAB lysis buffer

Tris-HCl, pH 8.0, 1M stock	100.0 ml
EDTA, pH 8.0, 500 mM stock	40.0 ml
NaCl	81.8 g
CTAB	20.0 g

Add deionized water to make up 1 L.

Mix chemicals, stir well to dissolve CTAB and autoclave. Store at room temperature.

- 0.25N HCl:

\*Dilute 5N HCl with Milli-Q water in 50 ml Falcon tube. Store at room temperature.

- 4 M spermidine:

\*Dissolve spermidine in Milli-Q or Sigma water. Divide small aliquots into 2 ml tubes and store at -20 °C. Avoid repeated thawing and freezing.

- 1 M spermine:

\*Dissolve spermine in Milli-Q or Sigma water. Divide to a small amount into 2 ml tubes and store at -20 °C. Avoid repeated thawing and freezing.

- RNase A (e.g. RNase A 100 mg/ml, Qiagen):

\*It is recommended to incubate RNase A at 95 °C for 5 minutes in a thermocycler to deactivate DNase contamination. This is mandatory for old stocks.

- 3M NaOAc (pH 5.2):

\*Use of this chemical is optional when DNA needs to be concentrated by precipitation.

- Low (0.1  $\times$ ) TE buffer (final concentrations):

Tris-HCl, pH 8.0, 1.0 mM  
EDTA, pH 8.0, 0.1 mM

\*Prepare in Milli-Q water and autoclave before use.

- Tris-HCl, pH 8.0, 1 M stock solution
- EDTA, pH 8.0, 500 mM stock solution
- Liquid Nitrogen
- PVPP (Polyvinyl polypyrrolidone)
- Isopropanol
- 70% ethanol (molecular grade)
- Proteinase K Qiagen #19131
- Genomic-tip 100/G Qiagen #10243
- Buffer G2 Qiagen #1014636
- Buffer QBT Qiagen #19054
- Buffer QC Qiagen #19055
- Buffer QF Qiagen #19056
- pH indicator strips to check pH range 7.0 - 8.0

### **Plastics**

- Falcon tubes: 50 ml
- Eppendorf tubes: 1.5 ml, 2 ml
- Eppendorf LoBind tubes: 1.5 ml
- Pipetman tips: 1ml, 1 ml wide-bore, 200 µl, 20 µl
- Nylon mesh 100 µm pore size or Corning cell strainer 100 µm pore size #431752

### **Equipment**

- Fume bench
- Water bath
- Heat block
- Centrifuge for 50 ml Falcon tubes
- Centrifuge for 1.5 ml and 2.0 ml Eppendorf tubes
- Pipetman: 1 ml, 200 µl, 20 µl
- Liquid nitrogen container
- 500 ml glass beakers
- Pestle and mortar (20 – 30 sets)
- Spectrophotometer: Nanodrop (Thermo Fisher Scientific) or DeNovix DS-11 (DeNovix Inc.)
- Fluorometer: Qubit (Thermo Fisher Scientific) or DeNovix DS-11
- Fragment analyser: TapeStation (Agilent) or Femto Pulse (Agilent)



## SAFETY WARNINGS

- Buffers containing  $\beta$ -mercaptoethanol should be handled with great care and appropriate personal protective equipment (PPE; *e.g.*, gloves and lab coat), and work has to be carried out in the fume hood or on the fume bench.
- Liquid nitrogen may cause cold burns, frostbite, and eye damage and needs to be handled with great care with PPE such as cryogenic gloves and safety goggles.






## BEFORE START INSTRUCTIONS

- Prepare and check plastic and chemicals listed under section “Materials”.
- Ready-to-use nuclei isolation (NI) buffer and sorbitol buffer containing  $\beta$ -mercaptoethanol, spermidine, spermine, and PVP40 should be prepared on the day of DNA extraction.
- Leave NI buffer and Triton X-100/NIB on ice for precooling.
- Affix 100  $\mu$ m pore size nylon mesh on top of a 500 ml glass beaker with string or tape. Mesh can be replaced with Corning cell strainers (100  $\mu$ m pore size) and 50 ml Falcon tubes.
- Depending on plant species and if possible, it is recommended to cultivate > 100 g plant material to allow protocol optimization.

### Before starting DNA extraction


- 1 Prepare 400 ml NI buffer in a glass beaker and leave on ice.  
 On ice
- 2 Prepare 100 ml sorbitol buffer in two 50 ml Falcon tubes.  
 On ice
- 3 Affix nylon mesh on top of one empty 500 ml glass beaker with string or tape (or place Corning cell strainer on 50 ml Falcon tube).
- 4 Arrange liquid nitrogen and mortar and pestles on fume bench.

## Tissue grinding

- 5 Grind fresh leaf tissue in liquid nitrogen and mortar and pestle 3 times to a fine powder and add ground tissue sample to NI buffer prepared at step 1. Grind 1-2 grams leaf tissue at a time and in total approx. 30 grams in this protocol.  
 On ice
- 6 Filter sample-NI buffer mix through 100 µm pore size nylon mesh/beaker prepared at step 3. Keep all solutions on ice during filtering.  
 On ice
- 7 Divide filtrate equally to ten 50 ml Falcon tubes kept on ice.  
 On ice
- 8 Add 1/20<sup>th</sup> volume of 10% Triton X-100/NIB to tubes prepared at step 7. Gently mix by inverting tubes.  
 On ice
- 9 Centrifuge tubes at 2,000 × *g* for 10 minutes at 4 °C. 10m  
 2000 x g, 4°C,  
00:10:00
- 10 Discard supernatant gently by decantation, without disturbing or losing pellet.
- 11 Add 10 ml sorbitol buffer to each tube and mix gently.


- 12 Centrifuge tubes at  $3,000 \times g$  for 10 minutes at  $4^\circ\text{C}$ .

10m

 3000 x g,  $4^\circ\text{C}$ ,  
00:10:00

- 13 Discard supernatant by decantation (Optionally, repeat sorbitol buffer wash until supernatant is clear).

- 14 To remove sorbitol buffer completely, invert tubes on dry tissue briefly but take care not to lose pellet. Pellet of nuclei and small cells remaining in tubes can now be **frozen in liquid nitrogen** and stored at  $-80^\circ\text{C}$ .

  $-80^\circ\text{C}$

[SAFE STOP POINT for at least a few days]

## CTAB lysis

- 15 Add 3 ml CTAB lysis buffer directly to frozen pellet, a pinch of PVPP, and  $12\ \mu\text{l}$  RNase A to each tube and mix well by gently pipetting with wide-bore tips. Incubate at  $58^\circ\text{C}$  for 20 minutes.


  $58^\circ\text{C}$

- 16 Add  $60\ \mu\text{l}$  proteinase K to each tube. Incubate for more than 3 hours, but less than 5 hours, at  $58^\circ\text{C}$ . Occasionally shake tubes gently.

  $58^\circ\text{C}$


- 17 Centrifuge at  $4,400 \times g$  for 10 minutes at room temperature. Collect clear lysate to new 50 ml Falcon tube avoiding any cell debris.

10m

 4400 x g, Room temperature,  
00:10:00

- 18 To maximize lysate recovery, move remained debris/lysate to 2 ml tubes and centrifuge at 11,000 rpm for 5 minutes. Move clear lysate to same tube at step 17. In total approx. 30 ml lysate can be obtained.

5m

 11000 rpm, Room temperature,  
00:05:00


- 19 Adjust lysate with 0.25 N HCl to between pH 7.0 - 7.5. Check pH with pH indicator strips. Add 1 ml

or less 0.25N HCl at a time and check with pH paper each time.






- 20** Divide lysate to two 50 ml Falcon tubes. Add equal volume of Milli-Q water.

*NOTE:* Do not centrifuge tubes once water is added. Low salt condition tends to promote formation of a CTAB-DNA solidified complex.

## Qiagen Genomic-tip 100/G DNA extraction

- 21** Proceed with Qiagen Genomic-tip 100/G following the manufacturer's protocol. Set up six empty 50 ml Falcon tubes, labelled "QBT", "Sample", "QC1", "QC2", "QC3", "Final DNA". Set up three sets of each.
- 22** Set water bath to 50 °C and prewarm buffer QF.  
 50 °C
- 23** Place Genomic-tip 100/G column on 50 ml Falcon tube labelled "QBT". Equilibrate Genomic-tip 100/G column with 4 ml buffer QBT. Allow buffer to flow through column completely by gravity.
- 24** Move Genomic-tip 100/G to next Falcon tube labelled "Sample". Load one third of lysate (approx. 20 ml) obtained at step 20 to one Genomic-tip 100/G. Allow lysate to flow through column completely by gravity.
- 25** Move Genomic-tip 100/G to tube labelled "QC1". Load 7.5 ml buffer QC onto column. Allow buffer to flow through column completely by gravity.
- 26** Repeat buffer QC step two more times, on tubes labelled "QC2" and "QC3". In total, DNA in Genomic-tip 100/G column should be washed three times with buffer QC.



- 27** For final DNA elution, apply 5 ml QF buffer prewarmed to 50 °C to each Genomic-tip 100/G column.  
 50 °C
- 28** Divide eluted DNA in 1 ml aliquots to 2 ml tubes. Add 0.7 volume (0.7 ml) ice-cold isopropanol. Gently invert and mix and leave tubes at -20 °C overnight.  
 -20 °C  
[SAFE STOP POINT for at least a few days]
- 29** Centrifuge tubes at 11,000 rpm for 10 minutes. 10m  
 11000 rpm,  
00:10:00
- 30** Discard supernatant and add 1 ml 70% ethanol.
- 31** Centrifuge at 11,000 rpm for 10 minutes. 10m  
 11000 rpm,  
00:10:00
- 32** Discard supernatant and air-dry pellet by inverting tubes on clean tissue. Warming tubes at 37 °C for 10 minutes speeds up evaporation of ethanol, but do not over-dry.
- 33** Add 15-20 µl Low (0.1 ×) TE buffer.
- 34** Incubate tubes at 50 °C, and 300 rpm for 1 hour. Collect eluted DNA to 1.5 ml LoBind Eppendorf tube.  
 50 °C
- 35** Add 15-20 µl Low (0.1 ×) TE buffer for 2<sup>nd</sup> elution.

- 36** Incubate tubes at 50 °C, and 300 rpm for 1 hour. Collect DNA elute to 1.5 ml LoBind Eppendorf tube. Keeping 1<sup>st</sup> and 2<sup>nd</sup> DNA elution in separate tubes is recommended.

 50 °C

## DNA quality control (QC)

- 37** Proceed to DNA quality control. DNA quantification with Qubit and Nanodrop. For DNA quality, obtain A260/A280 and A260/A230 values with Nanodrop. Examine DNA fragment size distribution with TapeStation or FemtoPulse.