



NOV 23, 2023

OPEN ACCESS



Protocol Citation: Allan JW Lui 2023. Media changes and Passaging in 2- or 5-layer CellStacks. **protocols.io** <https://protocols.io/view/media-changes-and-passaging-in-2-or-5-layer-cellst-cjxhupj6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Dec 01, 2022

Last Modified: Nov 23, 2023


PROTOCOL integer ID:
73417

Media changes and Passaging in 2- or 5-layer CellStacks

Allan JW

Lui¹

¹Cancer Research UK Cambridge Institute, University of Cambridge



Allan JW Lui

ABSTRACT

Protocol for handling CellStacks

MATERIALS

Corning CellStacks (2-layer: Corning 3269, 5-layer: Corning 3313)
Complete cell culture media
Phosphate-buffered saline
Trypsin solution (usually 0.25%, with EDTA and phenol red; Gibco 25200)
Centrifuge
50ml centrifuge tubes

BEFORE START INSTRUCTIONS

Ensure shelves of incubator are truly level with a spirit level or app

Media Changes

1 Per flask, Prepare and pre-warm:

Reagent	2-layer	5-layer
Complete media for plating	220ml	550ml

2 Take flask out of incubator, and place in upright position. Wipe down cap areas with 70% ethanol if desired.

Holding the flask with the long side facing downwards, uncap the bottom cap, and pour out media in flask onto a walls of a large 10L tub to prevent splashing, shake flask to remove the final few drops.




IMPORTANT: Handle flasks for cell culture in a suitable biosafety cabinet. Image taken for demonstration of flask handling only.

3 Place CellStack on its flat position.

For 2-layer CellStacks, pipette in new warm media 50-60ml at a time until 220ml has been added
For 5-layer CellStacks, pour in 550ml of new warm media (often an entire bottle)

- 4 To evenly distribute media between layers: Flip flask towards you such that it rests on its longer side.



For 2-layer CellStacks, place a finger on the bottom ledge to ensure flask is level.

IMPORTANT: Handle flasks for cell culture in a suitable biosafety cabinet. Image taken for demonstration of flask handling only.

Then tilt into the upright position:



- 5 Return flask to incubator, lying flasks flat while minimising sudden or large movements to avoid media spilling between layers

Passaging


8m

- 6 Per flask, Prepare and pre-warm:

Reagent	2-layer	5-layer

Reagent	2-layer	5-layer
PBS	250ml	500ml
Trypsin	40ml	100ml
Complete media for trypsin inactivation	40ml x2	100ml x2
Complete media for plating	220ml	550ml


- 7 Pour out media in flask, wash out remaining media **twice** by:
1. Adding PBS
 2. Distributing evenly between layers to cover the entire culture surface area
 3. Pouring out PBS


- 8 Add Trypsin 3m
- Incubate flask for around  00:03:00 in incubator
- Tap flask to dislodge cells

- 9 Add 1 aliquot of media for trypsin inactivation and distribute evenly between layers.
Collect cell suspension into duran bottles (250ml for 2-layer, 500ml for 5-layer)
Wash out remaining cells with the other aliquot of inactivation media similarly.

Reduce clumping by pipetting up and down multiple times with a 10ml serological pipette against the bottom of the bottle

- 10 Count cells using a Vi-Cell and calculate volume of cells for seeding

- 11 If necessary, centrifuge cells for replating at  200 x g, 00:05:00 and remove supernatant. 5m

-  If plating into CellStacks:
1. Resuspend cells in media for plating
 2. Transfer suspension into flask by pouring or pipetting

- 12 Reseed cells at desired density.

If plating into CellStacks:

1. Calculate the required volume of cell suspension for plating
2. Remove the same volume of media from the pre-warmed aliquot of media for plating and replace with



the cell suspension

3. Transfer suspension into flask by pouring or pipetting

Pelleting cells

5m

13 Transfer desired number of cells (usually 5×10^7) into 50ml centrifuge tubes

14 Centrifuge cells at  200 x g, 00:05:00 and remove supernatant. 5m
Resuspend with  5 mL PBS , centrifuge and remove supernatant again.

15 Flick tube to loosen pellet and freeze on dry ice.