



intracellular calcium assay

PLOS One

Wakana Shoda¹, Naohiro Nomura¹, Fumiaki Ando¹, Hideaki Tagashira², Takahiro Iwamoto², Akihito Ohta¹, Kiyosih Isobe¹, Takayasu Mori¹, Koichiro Susa¹, Eisei Sohara¹, Tatemitsu Rai¹, Shinichi Uchida¹

¹Department of Nephrology, Graduate School of Medical and Dental Sciences, Tokyo Medical and Dental University;

²Department of pharmacology, Faculty of Medicine, Fukuoka University

1 Works for me dx.doi.org/10.17504/protocols.io.baihicb6

Wakana Shoda

ABSTRACT

This protocol was used for the experiment conducted in "Sodium-calcium exchanger 1 is the key molecule for urinary potassium excretion against acute hyperkalemia"

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235360>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Shoda W, Nomura N, Ando F, Tagashira H, Iwamoto T, Ohta A, Isobe K, Mori T, Susa K, Sohara E, Rai T, Uchida S (2020) Sodium-calcium exchanger 1 is the key molecule for urinary potassium excretion against acute hyperkalemia. PLoS ONE 15(6): e0235360. doi: [10.1371/journal.pone.0235360](https://doi.org/10.1371/journal.pone.0235360)

DOI

[dx.doi.org/10.17504/protocols.io.baihicb6](https://doi.org/10.17504/protocols.io.baihicb6)

PROTOCOL CITATION

Wakana Shoda, Naohiro Nomura, Fumiaki Ando, Hideaki Tagashira, Takahiro Iwamoto, Akihito Ohta, Kiyosih Isobe, Takayasu Mori, Koichiro Susa, Eisei Sohara, Tatemitsu Rai, Shinichi Uchida 2020. intracellular calcium assay . **protocols.io**
[dx.doi.org/10.17504/protocols.io.baihicb6](https://doi.org/10.17504/protocols.io.baihicb6)

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Shoda W, Nomura N, Ando F, Tagashira H, Iwamoto T, Ohta A, Isobe K, Mori T, Susa K, Sohara E, Rai T, Uchida S (2020) Sodium-calcium exchanger 1 is the key molecule for urinary potassium excretion against acute hyperkalemia. PLoS ONE 15(6): e0235360. doi: [10.1371/journal.pone.0235360](https://doi.org/10.1371/journal.pone.0235360)

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0235360>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 16, 2019

LAST MODIFIED

Jul 01, 2020

Calcium Kit - Fluo 4 (Dojindo, Kumamoto, Japan)

- Fluo 4-AM 50 µg × 10
- Dimethylsulfoxide 2 ml × 1
- Recording Medium (2X) 100 ml × 1
- Probenecid 1.3 ml × 1
- 5% Pluronic® F-127 2.5 ml × 1

35mm glass bottom dish (Iwaki, Shizuoka, Japan)

1 Creation of Loading Buffer: solution A+ solution B

■Solution A

5ml Recording Buffer

50µl Probenesid

80µl Pluronic F-127

approximate 5 ml pure water (up to 10ml)

■Solution B

50 µl DMSO

50µl Fluo4

Loading Buffer should be heated at 37 °C

2 Cultured cells in glass bottom dish should be washed with PBS, and replaced by 1.2 ml loading buffer. Then, Incubate at 37 °C for 1 hour. The drugs are added during this period.

3 Preparation for Recording Buffer while cells are incubated in loading buffer.

■Recording Buffer

5ml Recording Medium

50 µl Probenesid

approximate 5ml pure water (up to 10ml)

4 After 1 hour incubation of loading buffer, the cells are washed with PBS, then replaced by 1.2 ml of recording buffer.

5 Following the administration of stimuli (ex. 10mM KCl), the fluorescence intensities of Fluo-4 are quantified from five regions of interest using LSM 510 Meta confocal microscopy and the ZEN 2009 software (Carl Zeiss, Oberkochen, Germany).