

Oct 23, 2025

Immunofluorescence protocol of Pax7 and androgen receptor for frozen muscle sections with unmasking

DOI

dx.doi.org/10.17504/protocols.io.sujeeun

Hiroshi Sakai¹

¹Ehime University



Hiroshi Sakai

Ehime University

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.sujeeun>

Protocol Citation: Hiroshi Sakai 2025. Immunofluorescence protocol of Pax7 and androgen receptor for frozen muscle sections with unmasking. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.sujeeun>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: August 23, 2018

Last Modified: October 23, 2025

Protocol Integer ID: 14955

Keywords: Pax7, Immunofluorescence, androgen receptors, muscle, frozen tissues, androgen receptor for frozen muscle section, immunofluorescence protocol of pax7, androgen receptor, studying skeletal muscle regeneration, skeletal muscle regeneration, frozen muscle tissue, member of the nuclear receptor, nuclear receptor, frozen muscle section, muscle tissue, unmasking immunofluorescence protocol, simple protocol for pax7, immunofluorescence protocol, pax7, immunofluorescence



Abstract

Immunofluorescence protocol of Pax7 on muscle tissues is a critical step for studying skeletal muscle regeneration. Here, I describe the simple protocol for Pax7 for isopentane-frozen muscle tissues with unmasking. Androgen receptor, a member of the nuclear receptor superfamily, can be stained with this protocol to see the double positive cells.

Guidelines

Use **freshly cut sections** since Pax7 and AR on sections are not stable.

Materials

MATERIALS

- ✕ Anti-Androgen Receptor antibody [SP107] - N-terminal **Abcam Catalog #ab105225**
- ✕ Goat anti-Rabbit IgG (H L) Cross-Adsorbed Secondary Antibody Alexa Fluor 488 **Thermo Fisher Scientific Catalog #A-11008**
- ✕ ProLong™ Glass Antifade Mountant **Invitrogen - Thermo Fisher Catalog #P36980**
- ✕ Anti-Pax7 antibody **Developmental Studies Hybridoma Bank Catalog #pax7**
- ✕ Goat anti-Mouse IgG1 Cross-Adsorbed Secondary Antibody Alexa Fluor 568 **Thermo Fisher Scientific Catalog #A-21124**

Troubleshooting

Before start

REAGENT SETUP

Sodium citrate (10 mM, pH = 6)

- A solution (Keep at room temperature)
 - 2.1 g of Citric Acid Monohydrate ($\text{C}_6\text{H}_8\text{O}_7 \cdot \text{H}_2\text{O}$)
 - 100 ml of H_2O
- B solution (Keep at room temperature)
 - 14.7 g of Trisodium citrate dihydrate ($\text{C}_6\text{H}_5\text{Na}_3\text{O}_7 \cdot 2\text{H}_2\text{O}$)
 - 500 ml of H_2O
- Sodium citrate
 - 9 ml of A solution
 - 41 ml of B solution
 - 450 ml of H_2O



Fix

- 1 Fix the slides with 4% PFA/PBS for 5 min at room temperature.

Note

Do not air dry the samples before fixing. It may increase the background signal and lose androgen receptor signal.

- 2 Wash the slides in PBS three times for 5 min.

Unmasking

- 3 Place the slides into a slide chamber in PBS

- 4 Fill a beaker with 500 ml of Sodium citrate.

- 5 Place the beaker on the hotplate and turn it on full power.

Note

While waiting for the beaker to come to a boil, keep the sections in PBS.

- 6 Once boiling, transfer the slides from PBS to the beaker.

Note

Use care with hot solution. Use forceps.

- 7 Boil the slides for 10 min.

- 8 Remove the beaker from the hotplate and let the beaker cool for 20 min.

- 9 Remove the slide chamber from the beaker.



- 10 Wash slides in PBS three times for 5 min.

Blocking

- 11 Block the slides with 5% Goat serum/0.05% Triton X-100/PBS for 1h at room temperature. Prepare 100 μ L of 5% Goat serum/0.05% Triton X-100/PBS per slide.
- 12 Drain the slides for a few seconds (do not rinse) and wipe around the sections with tissue paper.

Primary antibodies

- 13 Dilute primary antibodies as follows in 5% Goat serum/0.05% Triton X-100/PBS. Prepare 100 μ L antibody solution per slide. Add primary antibody solution to the slides.

🧴 100 μ L 5% Goat serum/0.05% Triton X-100/PBS

🧴 5 μ L anti-Pax7 (mouse IgG1, 1/20)

🧴 0.25 μ L anti-androgen receptor (rabbit, 1/400)

- 14 Incubate overnight at 4°C.
- 15 Wash the slides in PBS three times for 5 min.

Secondary antibodies

- 16 Dilute secondary antibodies as follows in 5% Goat serum/0.05% Triton X-100/PBS. Prepare 100 μ L antibody solution per slide. Add secondary antibody solution to the slides.

🧴 100 μ L 5% Goat serum/0.05% Triton X-100/PBS

🧴 0.1 μ L anti-mouse IgG1 568 (1/1000)

🧴 0.1 μ L anti-rabbit IgG 488 (1/1000)

🧴 0.05 μ L DAPI (1/2000)

- 17 Incubate for 1h at room temperature protected from light.



18 Wash the slides in PBS three times for 5 min.

Mount

19 Mount using ProLong™ Glass Antifade Mountant and add coverslips.