

MAR 20, 2023

🌐 Plate Count Agar

Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo, The National Institute of Occupational H...

ABSTRACT

Plate Count Agar (PCA), is a microbiological growth medium to assess viable bacterial growth of a sample.

The total number of living aerobic bacteria can be determined using PCA which is a substrate for bacteria to grow on. The medium contains **casein** which provides nitrogen, carbon, amino acids, vitamins and minerals to aid in the growth of the organism. **Yeast extract** is the source for vitamins and **glucose** is the fermentable carbohydrate. **Agar** is the solidifying agent. This is a non-selective medium and the bacteria is counted as colony forming units per gram (CFU/g) in solid samples and (CFU/ml) in liquid samples.

GUIDELINES

Prepare enough for the necessary number of experiments to perform in order to prevent batch-to-batch variation between experiments. Input chemicals should be of a certain quality in order to prevent contaminants that can impact microbial growth and selection.

MATERIALS

Magnetic stirrer, autoclave, scale, tubes, flask

SAFETY WARNINGS



Be sure to wear appropriate PPE when working with antibiotics and chemicals. Take care when working with hot flasks and tubes.

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.eq2ly74delx9/v1

Protocol Citation: Andreas Sagen 2023. Plate Count Agar. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly74delx9/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
 We use this protocol and it's working

Created: Mar 08, 2023


Last Modified: Mar 20, 2023





PROTOCOL integer ID:
 78371





Keywords: Plate count agar, agar, plate counting, counting, CFU


500 mL Plate Count Agar


2h



- 1 In a bottle, add approx.  400 mL deionized water 5m




- 2 Measure and add:
 5 g Tryptone
 2.5 g Yeast extract
 1 g Glucose
 15 g Agar 15m

Materials:
 Tryptone **Merck Millipore (EMD Millipore) Catalog #T9410**
 Yeast Extract **Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875**
 Glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021**
 Agar **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296**

- 3 Adjust pH to  7.0 using sodium hydroxide 5m

- 4 Fill bottle to  500 mL total with deionized water 5m

- 5 Autoclave liquid at  121 °C for  00:45:00 45m

- 6 Aliquot liquid in  15 mL and  50 mL units and store refrigerated at  4 °C 15m