



MAR 30, 2023

🌐 NGS grade DNA isolation from plant

Anurag Daware¹

¹University of Hohenheim



Anurag Daware
University of Hohenheim

ABSTRACT

NGS grade DNA isolation from plant (maize) leaf tissue

MATERIALS

1. Leaf tissues from 2-week-old seedlings collected from *Arabidopsis thaliana*, tomato, soybean, sorghum, and maize.
2. BBs (Daisy, catalog number: 980060-446)
3. Freezer boxes (Fisher Scientific, catalog number: 03-395-465, 100/CTN)
4. Eppendorf tubes (Dot Scientific, catalog number: 609-GMT; RA2000-GMT)
5. Ethylenediaminetetraacetic acid (EDTA) (Fisher scientific, catalog number: BP2482100)
6. Sodium chloride (NaCl) (Fisher scientific, catalog number: BP358-1)
7. Sodium dodecyl sulfate (SDS) (Fisher scientific, catalog number: 28312)
8. Potassium acetate (Fisher scientific, catalog number: A16321.36)
9. Isopropanol (Fisher scientific, catalog number: 040983.M1)
10. 100% Ethanol (Fisher scientific, catalog number: T038181000)
11. DNA loading dye (NEB, catalog number: 10816015)
12. 100 bp DNA ladder (GOLDBIO catalog number: D003-500)
13. 2x PCR Master Mix (syd labs, catalog number: MB067-EQ2G-L)
14. Liquid nitrogen
15. Pipettes and tips

OPEN ACCESS

Protocol Citation: Anurag Daware 2023. NGS grade DNA isolation from plant. **protocols.io** <https://protocols.io/view/ngs-grade-dna-isolation-from-plant-cr5nv85e>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Mar 30, 2023

Last Modified: Mar 30, 2023

PROTOCOL integer ID:
79758

Steps

- 1 Pre-heat water baths (65°C and 37°C) before beginning the extraction
- 2 Prepare the extraction buffer as follows:

- 5 mL 1 M Tris, pH 8 (final conc. 100 mM)
- 2.5 mL 0.5 M EDTA (final conc. 25 mM)
- 2.5 mL 5 M NaCl (final conc. 250 mM)
- 0.5 g PVP-40 (final conc. 1% w/v)
- 40 mL Water

- 3 Grind 100 mg of desiccated leaf tissue into a fine powder using bead beater
- 4 Add 650 µl extraction buffer (EB) and 90 µl 10% SDS directly into each tube and incubate for 30 min at 65 °C while inverting the tubes 3-6 times.
- 5 Add 250 µl KOAc to each tube and mix well by inverting several times. A white precipitate should form immediately.
- 6 Spin for 5-15 minutes at 12,000 rpm at room temperature and then transfer 400-600 µl of the supernatant to a fresh tube.
- 7 It's recommended to spin again and move the supernatant to a fresh tube if for NGS purposes.
- 8 Add an equal volume of isopropanol to the tube containing the supernatant. It's recommended to prepare tubes filled with an equal amount of isopropanol in advance. At this point, samples can be left at -20 °C overnight.
- 9 Spin samples for 5 minutes at 10,000 rpm and pour off the supernatant.
- 10 Add 400-500 µl of 70-75% ethanol to each tube, then spin for 5 minutes at 7500 rpm at RT.

- 11** Directly pour off the supernatant and air-dry samples completely. Tubes can either be left with caps opened on a rack at the bench or in a fume hood. It is recommended to put them upside-down on paper towels. This is an important measure to make sure the residual ethanol is gone as it can directly inhibit downstream enzyme activities. It is not recommended to remove liquid residuals with a pipette, as it may increase the chance of cross-contamination or loss of DNA.
- 12** For genotyping, add 100-200 µl nuclease free water into each tube. There is no need to add RNase.
- 13** For NGS purposes, it is important to add RNase, followed by one additional precipitation step (from No. 5 to 12), and to adjust the volume of water for each sample following quantitation to ensure uniformity.