



FEB 07, 2024

OPEN ACCESS



External link:

<https://doi.org/10.1371/journal.pgen.1010802>

Protocol Citation: Amanda A. G. Ferreira, Bogdan Sieriebriennikov, Hunter Whitbeck 2024. HCR RNA-FISH protocol for the whole-mount brains of *Drosophila melanogaster*.


protocols.io

<https://protocols.io/view/hcr-rna-fish-protocol-for-the-whole-mount-brains-o-cx26xqhe>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
An alternative version of [dx.doi.org/10.17504/protocols.io.bzh5p386](https://doi.org/10.17504/protocols.io.bzh5p386) aiming to use less reagent.

HCR RNA-FISH protocol for the whole-mount brains of *Drosophila melanogaster*

 Forked from [HCR RNA-FISH protocol for the whole-mount brains of *Drosophila* and other insects](#)

Amanda A. G. Ferreira¹, Bogdan Sieriebriennikov², Hunter Whitbeck¹

¹New York University; ²NYU Grossman School of Medicine

Yen



Yen-Chung Chen
New York University

ABSTRACT

This is a protocol to perform RNA fluorescent in situ hybridization (RNA-FISH) using hybridization chain reaction (HCR) on whole-mount samples of the brains of the fly *Drosophila melanogaster* and other insects, e.g. the jumping ant *Harpegnathos saltator*. Probes and HCR reagents are purchased from [Molecular Instruments](#). This protocol is loosely based on the "[generic sample in solution](#)" protocol published by [Molecular Instruments](#). Our modifications include the description of fixation conditions, counterstaining by Hoechst, and altered washes. Additionally, we use larger concentrations of probes and hairpins following the protocol described by Younger, Herre et al. 2020. We have successfully employed this protocol to stain insect brains with up to 4 different probe sets simultaneously (hairpins conjugated with Alexa Fluor 488, 546, 496, and 647).

CITATION

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall (2020). Non-canonical odor coding ensures unbreakable mosquito attraction to humans. bioRxiv.

LINK

<https://doi.org/10.1101/2020.11.07.368720>

Created: Aug 01, 2023

Last Modified: Feb 07, 2024

PROTOCOL integer ID: 85822

Keywords: brain, drosophila, ant, hcr, in situ, larva, rna fish

GUIDELINES

WORKING PRACTICES:

Prepare all buffers using nuclease-free water. Use filter tips and nuclease-free tubes. If using spot plates, pre-clean them first with household bleach diluted 1:10 in water and then with 70% ethanol. Wear gloves and adhere to other practices aimed at minimizing RNA degradation in the sample. Working in a clean bench is not required if other RNase-free practice are followed.

PROBE DESIGN:

We select the target sequence or isoform of the gene of interest and let Molecular Instruments design the probes. For genes with multiple isoforms, either target the isoform that includes as many as possible constitutive exons and as few as possible alternatively spliced exons, or the isoform that has the highest RNA-seq coverage (assessed visually in IGV) if RNA-seq data are available. Aim for the highest number of probes in a set, ideally 40, although we have successfully performed experiments with probe sets containing <20 probes.

AMPLIFIER CHOICE:

We routinely perform multiplexed stainings with up to 4 different probe sets and a Hoechst counterstain. We use amplifiers conjugated with Alexa Fluor 488, 546, 594, and 647. We are able to detect clearly distinguishable signals with minimal bleed-through on our confocal microscope (Leica SP8). However, be aware that simultaneously using fluorophores with partially overlapping spectra (e.g. AF 546 and 594) requires setting narrower detection ranges, which reduces the amount of signal detected.

MATERIALS

REAGENTS TO PURCHASE:



Nuclease-Free Water (not DEPC-Treated) **Thermo Fisher Scientific Catalog #AM9937**



20X PBS (Phosphate Buffered Saline) pH 7.4 **growcells.com Catalog #MRGF-6396**



Schneider's Drosophila Medium **Thermo Fisher Catalog #21720024**



Paraformaldehyde, 16% (wt/vol) **Electron Microscopy Sciences Catalog #15710**

Safety information

Paraformaldehyde is toxic, consult the SDS sheet for proper handling instructions

Note

Avoid long-term storage of the paraformaldehyde solution after opening the ampoule



Triton X-100 **Sigma Aldrich Catalog #X100**



10% Tween 20 **BIO-RAD Catalog #1662404**



20X SSC **Quality Biological Catalog #351-003-131**



SlowFade™ Gold Antifade Mountant **Invitrogen - Thermo Fisher Catalog #S36936**

or any other antifade mountant



Hoechst 33258, Pentahydrate (bis-Benzimide), 100 mg **Thermo Fisher Catalog #H1398**

- dissolve in DMSO to 5 mg / mL, aliquot and store at -20 °C



Methanol **Fisher Scientific Catalog #A412-4**



HCR Probe Hybridization Buffer **Molecular Instruments**



HCR Probe Wash Buffer **Molecular Instruments**

Safety information

Hybridization and Wash Buffers contain formamide, consult the SDS sheet for proper handling instructions



HCR Amplification Buffer **Molecular Instruments**

Buffer to prepare:

Note

Prepare fresh using nuclease-free water, store at 4 °C if required after the 1st day of the protocol.

- 1% PBTx (1X PBS with 1% v/v Triton X-100 and 1mM glycine)
- 5X SSCT (5X SSC with 0.1% v/v Tween-20)
- 1X PBS

SAFETY WARNINGS



This protocol uses solutions of paraformaldehyde and formamide, which are highly toxic chemicals. Consult the SDS sheets of the reagents used in this protocol for proper handling instructions.

Day 1

1h 45m

1 Prepare all solutions

15m

2 Pre-heat an aliquot of Probe Hybridization Buffer if proceeding with hybridization on the same day (see step 2 case below)

2m

🧴 250 µL / sample (If preparing a new probe mixture)

🧴 100 µL / sample (If reusing probes)

🌡️ 37 °C

3 Dissect brains in cold Schneider's medium

10m

Note

Dissection can also be done in 1X Nuclease-Free PBS.

4 Fix brains in 800 µL of 4% PFA / Schneider's medium

20m

⌚ 00:20:00

🌡 Room temperature

🌀 24 rpm Nutator or 🌀 60 rpm Orbital shaker

Note

For fixation and all subsequent steps, samples can be placed either in Eppendorf tubes or in wells of a spot plate (e.g. Pyrex spot plate with 9 depressions, Catalog #CLS722085). Tubes are incubated on a nutator and plates are incubated on an orbital (horizontal) shaker.

5 Rinse 3x with 150 μ L of 1% PBTx

3m

🌡 Room temperature

6 Wash 3x 15 min with 150 μ L of 1% PBTx

45m

⌚ 00:15:00

⌚ 00:15:00

⌚ 00:15:00

🌡 Room temperature

🌀 24 rpm Nutator or 🌀 60 rpm Orbital shaker

7 Pre-hybridize samples by incubating them with 100 μ L of warm Probe Hybridization Buffer from step 2

10m



⌚ 00:10:00 can be extended to 30 min

🌡 37 °C

Note

Hybridization buffer is viscous, and the brains might be stuck on the wall. Be careful not to remove them or let them dry out by accident.

8 In the meantime, prepare a 8nM probe solution by adding 0.4 pmol of each probe mixture (e.g. 0.4 μ L of 1 μ M stock) to the warm Probe Hybridization Buffer for the total volume of 100 μ L.

🧴 0.4 μ L 1 μ M probe per sample

🧴 100 μ L Warm probe hybridization buffer per sample

9 Remove the pre-hybridization solution from the sample and add the probe solution from step 9

3m

Note

Take extra care while removing the liquid. Hybridization Buffer is viscous and brains may float.

10 Incubate samples with the probes



🕒 Overnight We usually do ~24 h, but incubation can be extended to ~48 h. Minimum recommended is 12 h

🌡️ 37 °C

🔄 24 rpm Nutator or 🔄 60 rpm Orbital shaker

Day 2

1h 20m

11 Equilibrate an aliquot of Amplification Buffer to room temperature

2m

🌡️ Room temperature

🧴 60 µL per amplifier per sample (if preparing new) + 🧴 110 µL per sample

Note

If using 3 different amplifiers for 2 samples, $60 * 3 * 2 + 110 * 2 = 580\mu\text{L}$ is to be equilibrated.

12 Pre-heat an aliquot of Probe Wash Buffer

2m

🌡️ 37 °C

🧴 450 µL per sample

13 Remove excess probes by washing 4 x 10 min with 100µL pre-heated Probe Wash Buffer from step 12

40m

🕒 00:10:00 Wash 100µL pre-heated probe wash buffer

🕒 00:10:00 Wash 100µL pre-heated probe wash buffer

🕒 00:10:00 Wash 100µL pre-heated probe wash buffer

🕒 00:10:00 Wash 100µL pre-heated probe wash buffer

🌡️ 37 °C

🔄 24 rpm Nutator or 🔄 60 rpm Orbital shaker

14 Wash samples 2 × 5min with 1mL of 5X SSCT at room temperature.

10m

🕒 00:05:00 5x SSCT

🕒 00:05:00 5x SSCT

15 Equilibrate each sample with 100 µL of Amplification Buffer

10m

🕒 00:10:00 can be extended to 30 min

🌡️ Room temperature

🔄 24 rpm Nutator or 🔄 60 rpm Orbital shaker

16 During sample equilibration, warm up the thermo cycler to 95°C.

17 Prepare **separately** each hairpin (h1 and h2) of each amplifier. Each sample require 2 µL of **each** 3 µM stock of hairpin to be mixed in a final volume of 100µL/sample amplification buffer. Aliquot the required amount of hairpin to **individual** PCR tubes. For example, if amplifiers B1 and B2 are being used, prepare 4 PCR tubes that contain B1-h1, B1-h2, B2-h1, and B2-h2, respectively.

30m

Incubate the tubes in a thermocycler at 95 °C for 90 sec. Immediately take them out of the machine (while it is still at 95 °C), place them in a rack and incubate them at room temperature **in the dark** for at least 30 min.

🕒 00:30:00 or longer

🌡️ Room temperature

18 Briefly spin down the hairpin solutions, add 100 µL/sample amplification buffer, and use **the same** buffer to resuspend all hairpins used.

🌡️ Room temperature

19 Remove 100 µL liquid from the samples and add the Amplification Buffer + hairpins from step 17

2m

Note

Take extra care while removing the liquid. Amplification Buffer is viscous and the brain may float.

20 All the next steps are light sensitive and must be done in the dark!



21 Incubate in **the dark**.



🕒 Overnight

🌡 Room temperature

🌀 24 rpm Nutator or 🌀 60 rpm Orbital shaker

Day 3

3h 15m

22 Prepare 1.2 ml/sample 5X SSCT

23 Save hairpin mixture and keep at -20°C

2m

24 Wash 2 x 5 min with 150 μ L of 5X SSCT

10m

🕒 00:05:00

🕒 00:05:00

🌡 Room temperature

🌀 24 rpm Nutator or 🌀 60 rpm Orbital shaker

25 Wash 2 x 15 min with 150 μ L of 5X SSCT

30m

🕒 00:15:00

🕒 00:15:00

🌡 Room temperature

🌀 24 rpm Nutator or 🌀 60 rpm Orbital shaker

26 If not using DAPI, skip this step and the next and proceed to step 28.

2h



If using DAPI: incubate samples with 150 μ L of SSCT + DAPI (1:1000 of the 5 mg/mL stock - final concentration 10 μ g/mL) for 1 h

🕒 02:00:00

- 🌡 Room temperature

🔄 24 rpm Nutator or 🔄 60 rpm Orbital shaker
- 27 Wash 2 x 10 min with 150 µL of 5X SSCT 20m

🕒 00:10:00

🕒 00:10:00

🌡 Room temperature

🔄 24 rpm Nutator or 🔄 60 rpm Orbital shaker
- 28 Rinse with 2X Nuclease-Free PBS to remove the detergent 2m
- 29 Remove the PBS and add Slowfade (or another antifade mountant) 2m
- 30 If necessary, finish dissecting to remove extra tissue. Mount brains on slides. 10m
- 31 Proceed with imaging

