



FEB 02, 2024

Clusterin cellular uptake assay

Patricia Yuste Checa¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

ABSTRACT

This protocol details how to efficiently monitor Clusterin and Clusterin/Substrate uptake in different cell types, like HEK293T, iNeurons and iMicroglia.

ATTACHMENTS

[975-2532.docx](#)

GUIDELINES

- To study Clusterin-A488 uptake in the presence of substrate, incubate [M] 1 micromolar (μM) Clusterin-A488 with the corresponding amount of substrate, e.g. denatured Luciferase or A β aggregates, in PBS for 00:20:00 at 37 °C or 42 °C (denatured Luciferase) in a total volume of 30 μL or 40 μL in the case of HEK293T or iNeurons, respectively. After the incubation time dilute the mix 1/10 in media and add it to the cells resulting in a final concentration of [M] 0.1 micromolar (μM) Clusterin-A488 (5 undetermined).
- To monitor substrate uptake e.g., denatured luciferase or A β aggregates, in the presence of Clusterin, mix [M] 0.2 micromolar (μM) of denatured Luciferase-pHrodo or [M] 0.5 micromolar (μM) of A β -pHrodo aggregates with the corresponding amount of unlabeled Clusterin in a total volume of 30 μL or 40 μL in the case of HEK293T or iNeurons, respectively. Dilute the mix 1/10 in media. pHrodo Red dye is pH sensitive dye which fluoresces brightly only in acidic environments and therefore can be used to specifically monitor phagocytosis and endocytosis, but substrates labeled with A488 can be also used.
- The indicated Clusterin-A488 or substrate concentrations and incubation times for each cell line are tentative. These parameters should be experimentally tested to be in the linear range of the assay.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.14egn3k2yl5d/v1

Protocol Citation: Patricia Yuste Checa, F Ulrich Hartl 2024. Clusterin cellular uptake assay. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn3k2yl5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 25, 2024

PROTOCOL integer ID: 94599

Keywords: ASAPCRN

Funders Acknowledgement:

Aligning Science Across
Parkinson's
Grant ID: ASAP-000282

ACCUTASE™ 100 mL STEMCELL Technologies Inc. Catalog #7920

Trypan Blue Solution 0.4% Thermo Fisher Scientific Catalog #15250061

Clusterin cellular uptake assay - HEK293T cells

25m

1 Plate 100,000 HEK293T cells per well in a 24-well plate.

2 On the next day, add 5 undetermined of Clusterin-A488 together with 300 µL of fresh DMEM (without fetal bovine serum, 1.5 µg in 300 µL medium) to the cells and place the cells back in the incubator.



3 After 04:00:00 incubation, place the plate on ice to stop endocytosis.



4h

4 Wash the cells gently with cold 1x PBS.



5 Add 100 µL TrypL Express Enzyme (Gibco). Incubate for few minutes On ice .






6 Collect the cells with  400 μL of cold medium and transfer them to an Eppendorf tube placed  On ice .

7 Centrifuge at  1000 x g, 4°C, 00:05:00 .

5m



8 Discard the supernatant and fix the cells by resuspending the cell pellet with  200 μL 4% Paraformaldehyde (PFA) in 1x PBS. Incubate for  00:10:00 at  Room temperature .

10m



9 Centrifuge at  1000 x g, 4°C, 00:05:00 .

5m



10 Wash the cell pellet with 1x PBS.



11 Centrifuge at  1000 x g, 4°C, 00:05:00 .

5m



12 Resuspend the cell pellets with  160 μL 1x PBS  7.2 and store at  4 °C until analyzed.

iNeurons

1h

13



Add 5 undetermined of Clusterin-A488 to 250,000 iNeurons cultured in a well of a 12-well plate (add 2 μg Clusterin-A488 to 200 μL of fresh medium and add the mix to the well with cells containing 200 μL conditioned medium) and place the cells back in the incubator.

14

After 01:00:00 incubation place the plate On ice to stop endocytosis.

1h

15



Wash the cells gently with cold 1x PBS.

16



Add 100 μL Accutase. Incubate for 5-10 minutes On ice.

17

Collect the cells with 400 μL of cold medium and transfer them to Eppendorf low binding tubes placed On ice.

18


















Centrifuge at 1000 x g, 4°C for 00:05:00 (swing-bucked centrifuge preferred).

5m







Note



The use of low binding tubes and swing-bucked centrifuge significantly reduces cell loss.


- 19 Discard the supernatant and fix the cells by resuspending the cell pellet with  200 μ L 4% PFA in 1x PB  10m
 Incubate for  00:10:00 at  Room temperature .
- 20 Centrifuge at  1000 x g, 4°C, 00:05:00 (swing-bucked centrifuge preferred).  5m

- 21 Wash the cell pellet with 1x PBS.

- 22 Centrifuge at  1000 x g, 4°C, 00:05:00 (swing-bucked centrifuge preferred).  5m

- 23 Resuspend the cell pellets with  160 μ L 1x PBS  7.2 and store at  4 °C until analyzed.

iMicroglia

1h 15m

- 24 Dispense 150,000 iMicroglia cells per well with  300 μ L medium into a Geltrex-coated 24-well plate.
- 25 Add 5  5 undetermined of Clusterin-A488 ( 1.5 μ g in  300 μ L medium) and place the cells back in the incubator.



26 After  00:30:00 incubation place the plate on ice to stop endocytosis, collect the cells and transfer them to Eppendorf low binding tubes placed  On ice. 30m

27 Centrifuge at  1000 x g, 4°C, 00:05:00 (swing-bucked centrifuge preferred). 5m






Note


The use of low binding tubes and swing-bucked centrifuge significantly reduces cell loss.

28 Wash the cell pellet with 1x PBS.



29 Discard the supernatant and fix the cells by resuspending the cell pellet with  200 µL 4% PFA in 1x PBS 10m
Incubate for  00:10:00 at  Room temperature.




30 Centrifuge at  1000 x g, 4°C, 00:05:00 (swing-bucked centrifuge preferred). 5m



31 Wash the cell pellet with 1x PBS.



32 Centrifuge at  1000 x g, 4°C, 00:05:00 (swing-bucked centrifuge preferred). 5m



33 Resuspend the cell pellets with  160 µL 1x PBS  7.2 and store at  4 °C until analyzed.




Uptake quantification

1h 15m

34



Quantify Clusterin or substrate uptake by measuring A488 or pHrodo Red intensity inside the cells by flow cytometry. If A488 is used, add  50 μ L of Trypan blue solution 0.4% (refer materials section) right before measuring to quench the 488 fluorescence outside the cells.

Note

An Attune NxT flow cytometer (Thermo Fisher Scientific) can be used with the following settings:

- Alexa485: Excitation 488 nm - Emission 550/30.
- pHrodo Red: Excitation 561 nm - Emission 585/16.