



JAN 22, 2024

🌐 S-Complete Medium

📁 In 1 collection

Adrien Assie¹, Buck Samuel¹

¹Baylor College of Medicine

Samuel Lab



Adrien Assie

Baylor College of Medicine

OPEN  ACCESS



External link:

http://www.wormbook.org/chapters/www_strainmaintain/strainmaintain.html

Protocol Citation: Adrien Assie, Buck Samuel 2024. S-Complete Medium. **protocols.io**

<https://protocols.io/view/s-complete-medium-zzaf72e>

MANUSCRIPT CITATION:

Lewis, J.A. and Fleming, J.T. (1995). In: Methods in cell biology, Vol. 48, H.F. Epstein and D.C. Shakes, eds. (San Diego: Academic Press), p. 3.

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

DISCLAIMER

We are not the author of this protocol. It was first described in Lewis and Fleming (1995). Protocol inspired from the one available on the Wormbook website.

ABSTRACT


Large quantities of *C. elegans* can be grown in a liquid medium. Liquid cultures of *C. elegans* are usually grown on S Medium using concentrated *E. coli* OP50 as a food source. This is the S-Complete medium recipe.

Created: Apr 10, 2019


Last Modified: Jan 22, 2024

PROTOCOL integer ID: 22274

Keywords: C. elegans, Liquid culture, S medium

1  977 mL of S-Basal medium (Link below)

Protocol



NAME
S-Basal Medium

CREATED BY
Adrien Assie

PREVIEW


1.1 Start with  700 mL water

1.2  5.9 g NaCl ( 100 millimolar (mM))


1.3  50 mL of  1 Molarity (M) Potassium Phosphate Buffer, pH 6.0

OR


 1 g K₂ HPO₄ and  6 g KH₂PO₄

1.4 Adjust water to  1000 mL

1.5 Autoclave

2  1 mL of 5 mg/ mL cholesterol (dissolved in Et-OH)

3  10 mL of  1 Molarity (M) Potassium Citrate Buffer, pH 6.0

4  10 mL of Trace metals solution (sterile)

Protocol



NAME



Trace Metals Solution

CREATED BY

Adrien Assie

PREVIEW


4.1 Start with  700 mL water

4.2  10 mL of *500 mM EDTA, pH 8.0 or  1.86 g Na₂EDTA.2H₂O (5 mM)

4.3  0.69 g $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ (2.5 mM)

4.4  0.20 g $\text{MnCl}_2 \cdot 4\text{H}_2\text{O}$ (1 mM)

4.5  0.29 g $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$ (1 mM)

4.6  0.016 g CuSO_4 (0.1 mM)

4.7 Adjust to  1000 mL water

4.8 Autoclave and store in the dark

5  3 mL of 1 M CaCl_2 (sterile)

6  3 mL of 1 M MgSO_4 (sterile)