



NOV 13, 2023

## 🌐 iPSCs Maintenance and Banking

Mahmoud ElAchwah<sup>1</sup>, sol.diazdeleon<sup>1</sup>, pangzh<sup>1</sup>

<sup>1</sup>Rutgers University

MINND



Mahmoud ElAchwah  
Rutgers University

### ABSTRACT

This protocols offers a thorough description of the maintenance and banking of induced pluripotent stem cells.

### MATERIALS

#### Matrigel:

Matrigel catalog number: **354234**

Matrigel stock concentration: **9-11 µg/µL**  
**110 µg/mL**

Matrigel Stock aliquot volume: **400 µL**

Matrigel LoT number: **13823002**

Matrigel working concentration: **90-**

#### Culture Antibiotic:

Primocin catalog number: **ant-pm-1**

Primocin stock concentration: **50 mg/mL**  
concentration: **100 µg/mL**

Primocin stock aliquot volume: **1 mL**

Primocin LoT number: **NA**

Primocin working

#### CEPT:

Chroman catalog number: **HY-15392**  
concentration: **50 nM**

Chroman stock concentration: **0.5 mM**

Chroman working

Emricasan catalog number: **S7775**  
concentration: **5000 nM**

Emricasan stock concentration: **50 mM**

Emricasan working

Polyamine supplement catalog number: **P8483**

Polyamine supplement

OPEN ACCESS



DOI:  
[dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1)

**Protocol Citation:** Mahmou d ElAchwah, sol.diazdeleon, pangzh 2023. iPSCs Maintenance and Banking. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewov1qd5ogr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 15, 2023

**Last Modified:** Nov 13, 2023

**PROTOCOL integer ID:** 87845

working concentration: **NA**  
Polyamine supplement stock concentration: **NA**

Trans-ISRIB catalog number: **5284**  
concentration: **700 nm**  
Trans-IRIB stock concentration: **7 mM**

Trans-IRIB working

#### Culture Media:

MEM catalog number: **51200-038**  
**2451155**

MEM LoT number:

Accutase catalog number: **AT-104**  
**2B2527A**

Acutase LoT number:

mTesr+ basal medium catalog number: **100-0274**  
**AH29845535**  
mTesr+ 10x supplement: **100-0275**  
mTesr Plus kit: **100-0276**

mTesr+ LoT number:

mTESR+ supplemented media



- Add 100ml of mTESR 10x supplement into 400ml of mTESR+ media
- Filter and aliquot 45ml into conical tubes, keep at -20°C.









#### mFreSR Media:

mFreSR media catalog number: **05855**  
mFreSR media stock concentration: **NA**  
concentration: **NA**  
mFreSR media Stock aliquot volume: **50 mL**

mFreSR media LoT number: **NA**  
mFreSR media working






## Diluting Matrigel

- 1 Thaw Matrigel stock [M] 9-11 µg/µL stored at  -80 °C on ice
- 2 Make  400 µL aliquots of Matrigel to avoid multiple freeze-thaw cycles using previously cooled tips and tubes

- 3 Store  400  $\mu\text{L}$  Matrigel aliquots at   $-20\text{ }^{\circ}\text{C}$  until use
- 4 To prepare matrigel working solution: thaw  400  $\mu\text{L}$  Matrigel aliquots  On ice or at   $4\text{ }^{\circ}\text{C}$
- 5 Dilute  400  $\mu\text{L}$  Matrigel aliquot in  40 mL of MEM (Working Concentration 90-110  $\mu\text{g}/\text{mL}$ )
- 6 Mix well before use. Store at   $4\text{ }^{\circ}\text{C}$  protected from the light













## Matrigel Coating

1h



- 7 Coat wells with  1 mL of Matrigel per well of 6 well plates or per 35 mm plate
- 8 Leave coated plates for  02:00:00 to  04:00:00 in the incubator at   $37\text{ }^{\circ}\text{C}$
- 9 Store diluted matrigel at   $4\text{ }^{\circ}\text{C}$

6h

## Medium Preparation For Passaging


- 10 Prepare mTsr+ by adding  100 mL of mTsr+ 10x supplement into  400 mL of mTESR+ media and filter it.
- 11 Add  1 mL of  50 Molarity (m) Primocin (final concentration 100ug/mL)
- 12 Take out  3 mL of MEM per well of 6 well plates or per 35 mm plate in a conical tube ( 1 mL /well to wash old medium and  2 mL /well to wash out acutase) and leave at room temperature to warm up
- 13 Take out  1 mL of acutase per well of 6 well plate or per 35mm plate in  15 mL conical tube to warm up at RT
- 14 Take out  3 mL of mTESR+ per well of 6 well plate or 35 mm plate ( 2 mL /well for plating and  1 mL /well for resuspension) to warm up at RT



## Washing iPSCs

- 15 Put iPSC plate out of the incubator at  Room temperature
- 16 Aspirate old medium
- 17 Add  1 mL of MEM per well of 6 well plate or per 35 mm plate, rock the plate to wash

**Detaching Adherent Cells**

4m


19 Add  1 mL of acutase per well of 6 well plate or per 35 mm plate

20 Place plate in the incubator at  37 °C for  00:05:00

5m

**Centrifuging Cells**

21 Get plate out of incubator and check that cells are detached, avoid pippeting too much.

Add  2 mL of MEM

22 Transfer detached cells with MEM ( 3 mL total volume) into 15 mL conical tube

23 Place conical tubes in centrifuge and spin at  1000 rpm for  00:05:00 at  Room temperature

5m

**Preparing Seeding Medium**

24 Add  1 µL of CET per  1 mL of mTESR+ (1:1000 ratio)

25 Add 1  $\mu$ L of P per 1 mL of mTesr+ (1:1000 ratio)

26 Mix mTESR+ with CETP

27 Aspirate matrigel from coated plates

28 Add 2 mL of mTesr+ containing CETP to each well

## Collecting Cells

29 Aspirate supernatant from 15 mL conical tube containing cells (should have pellet)

30 Resuspend pellet with 1 mL of mTESR+ very gently to avoid dissociating iPSCs into single cells

31 Add 50  $\mu$ L of cell solution per well in 6 well plate or per 35 mm plate (on top of the mTESR+ with CETP)

32 Gently mix the plate to evenly spread the cells


## Incubation

5m



- 33 Ensure the presence of cells in each well by viewing the plate under the microscope (should see floating cells)
- 34 Incubate plate in the incubator at  37 °C  Overnight

5m

## Changing iPSC medium

- 35 Aspirate old medium the next day to wash the CETP
- 36 Add  2 mL of fresh mTESR+ to each well of the 6 well plate or per 35mm plate (no wash with MEM first day after splitting)

## Wash and iPSC Medium Change



- 37 Wash CETP from each well/plate with  1 mL of MEM
- 38 Aspirate wash
- 39 Add  2 mL of fresh mTESR+ to each well of the 6 well plate or per 35 mm plate
- No need to keep washing plates after the first wash, you can directly change medium

## iPSC Maintenance



- 40 Keep feeding cells everyday until they are ready for splitting again (cells are usually ready for another splitting in 2-3 days, when their confluency reach 70%)

Check cells under microscope to estimate confluency

## iPSC Freezing

- 41 **For freezing purposes**, follow the same protocol shown above for cell maintenance for the following sections: Washing iPSCs, Detaching Afdherent cells, centirfuging cells (steps 10-23)
- 42 Instead of collecting cells have freezing vials set and barcoded with iPSC line number, passage number, gene mutation and date
- 43 Resuspend pellet with  1 mL of mFreSR very gently
- 44 Add  500 µL of cell solution to freezing vials

## iPSC Banking

- 45 Store freezing vials at  -80 °C in a cryogenic freezing container to prevent ice crystals from forming within the cells in the freezing vials
- 46 After 24 hours, move freezing containers with vials from  -80 °C to a nitrogen tank