

JAN 23, 2024

🌐 CompDuplex: Accurate detection of somatic mutations by duplex-seq with comprehensive genome coverage

Muchun Niu¹, Chenghang Chuck Zong¹

¹Department of Molecular and Human Genetics, Baylor College of Medicine



Muchun Niu

DISCLAIMER

Patent application covering CompDuplex chemistry has been filed by Baylor College of Medicine.

ABSTRACT

Somatic mutations continuously accumulate in the human genome, posing vulnerabilities towards aging and increased risk of various diseases. However, accurate detection of somatic mutations at the whole genome scale is still challenging. By tagging and independently sequencing the two complementary strands of DNA, the recent development of duplex-sequencing methods has greatly improved the detection accuracy, however, the limited genome coverage and the compromised compatibility with existing sequencing platforms have constrained the broad applications of these methods.

To overcome these technical challenges, here we developed a duplex sequencing method with comprehensive genome coverage, which we refer to as CompDuplex-seq. The streamlined chemistry of CompDuplex assay allows efficient generation of libraries readily compatible with standard Illumina 2x150 paired-end sequencing. In addition, we validated the accuracy of somatic mutation calling and comprehensive genome coverage of CompDuplex by profiling a single-cell expanded clone. To summarize, CompDuplex chemistry supports genome-wide coverage while maintaining high accuracy, which we believe will facilitate the whole genome characterization of somatic mosaicism in various biological systems.

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.kxygx3x4og8j/v1

Protocol Citation: Muchun Niu, Chenghang Chuck Zong 2024. CompDuplex: Accurate detection of somatic mutations by duplex-seq with comprehensive genome coverage. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.kxygx3x4og8j/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 17, 2024



















Last Modified: Jan 23, 2024

PROTOCOL integer ID: 93708

Keywords: Duplex-seq, Somatic mutation, Accurate mutation calling

Funders Acknowledgement: SMaHT Network (NIH) Grant ID: UG3NS132132

PROTOCOL MATERIALS


-  NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X** In [2 steps](#)
-  0.5M EDTA **Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML** Step 8
-  T4 Polynucleotide Kinase Reaction Buffer **New England Biolabs Catalog #B0201S**
- Step 9
-  1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025** In [2 steps](#)
-  iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**
- Step 11
-  Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516** Step 3
-  Water **Invitrogen - Thermo Fisher Catalog #46-2224** Step 5
-  BtgZI - 500 units **New England Biolabs Catalog #R0703L** Step 7
-  10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML**
- Step 8
-  5M NaCl **Invitrogen - Thermo Fisher Catalog #AM9760G** Step 1
-  Tn5 transposase, unloaded **Diagenode Catalog #C01070010-20** Step 3
-  Tagmentation Buffer 1 **Illumina, Inc. Catalog #20015171** Step 5
-  rCutSmart Buffer **New England Biolabs Catalog #B6004S** Step 7
-  StickTogether DNA Ligase Buffer **New England Biolabs Catalog #B0535S** Step 10
-  Ampure XP beads **Beckman Catalog #A63881** In [5 steps](#)
-  Protease **Qiagen Catalog #19155** Step 8
-  T7 DNA Ligase - 750,000 units **New England Biolabs Catalog #M0318L** Step 10
-  Water **Invitrogen - Thermo Fisher Catalog #46-2224** Step 1

Custom Tn5 transposase assembly

1h

1

2m

Prepare  500 μ L of **2X Annealing Buffer**.

-  40 μ L  1M Tris-HCl, pH 8.0 **Invitrogen - Thermo Fisher Catalog #15568025**
-  10 μ L  5M NaCl **Invitrogen - Thermo Fisher Catalog #AM9760G**
-  450 μ L  Water **Invitrogen - Thermo Fisher Catalog #46-2224**

2 Transposon annealing.

27m

Prepare the following mix in a PCR tube.

🧴 20 µL	2X Annealing Buffer
🧴 10 µL	[M] 200 micromolar (µM) ME_REV
🧴 10 µL	[M] 200 micromolar (µM) T_BtgZ1_ME

Oligonucleotides sequences

ME_REV: /5Phos/CTGTCTCTTATACACATCT

T_BtgZ1_ME: ATGTGTGGAGCGATG AGATGTGTATAAGAGACAG

Mix up the reaction, spin down, and perform the following reactions on a thermal cycler.

🔥 95 °C 5 min
🔥 65 °C 5 min, Ramp rate 0.1°C/s
🔥 4 °C Hold, Ramp rate 0.1°C/s

This results in a [M] 50 micromolar (µM) transposon mix .

3 Custom Tn5 assembly.

31m

Prepare the following mix in a 1.5 mL tube.







🧴 20 µL	🧴 Tn5 transposase, unloaded Diagenode Catalog #C01070010-20
🧴 20 µL	[M] 50 micromolar (µM) transposon mix


Gently pipet to mix, and incubate at 🔥 23 °C for 30 min.

Add 🧴 20 µL 🧴 Glycerol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G5516** to the tube, mix well, aliquot, label as “**Tn5_BtgZ1**”, and store at 🔥 -20 °C .

CompDuplex Procedure






4 1. Re-purification of genomic DNA

Dilute  20 ng  Sample to  10 µL nucleus-free water in a PCR tube. Add  5 µL  Ampure XP beads **Beckman Catalog #A63881** (0.5X) to the tube, mix well, and incubate at  Room temperature for 5 min.

Put the tube on the magnetic stand. When the aqueous phase becomes transparent, remove the supernatant, wash twice with  120 µL **[M] 80 % volume Ethanol**. Air-dry the beads.


5 2. On-bead custom Tn5 tagmentation


Prepare the following **tagmentation mix**.

 2 µL  Tagmentation Buffer 1 **Illumina, Inc. Catalog #20015171**
 0.48 µL 1:30 diluted **Tn5_BtgZ1**
 7.52 µL  Water **Invitrogen - Thermo Fisher Catalog #46-2224**


Resuspend the air-dried beads with  10 µL **tagmentation mix**.

Incubate on a thermal cycler using the following program:  20 °C for 5 min,  55 °C for 15 min.

Quench the reaction by adding  3 µL of 0.2M EDTA.


Incubate on a thermal cycler at  50 °C for 30 min.



6 3. Gap filling

Add  3 µL of 0.2M MgCl₂ to the tube.

Add  16 µL  NEBNext Ultrall Q5 Master Mix **New England Biolabs Catalog #M0544X** to the tube.

Incubate on a thermal cycler at  65 °C for 3 min.

Quench the reaction by adding  2 µL of 0.5M EDTA.

Perform double size selection with 0.48X/0.75X  Ampure XP beads **Beckman Catalog #A63881**, elute in  7.5 µL water.

7 4. Restriction enzyme digestion

Prepare the following restriction enzyme digestion mix (Total  15 µL).

7.5 μL Gap filling product

1.5 μL rCutSmart Buffer New England Biolabs Catalog #B6004S

0.7 μL BtgZI - 500 units New England Biolabs Catalog #R0703L

5.3 μL Water

Incubate on a thermal cycler with the following program: 25 °C for 10 min, 37 °C for 3 h, 10 °C for 12 h.

8 5. Restriction enzyme release

Prepare the following **restriction enzyme release mix** (Total 500 μL).

15 μL 1M Tris-HCl, pH 8.0 Invitrogen - Thermo Fisher Catalog #15568025

12 μL 0.5M EDTA Merck MilliporeSigma (Sigma-Aldrich) Catalog #E7889-100ML

10 μL 1M KCl

5 μL 10% Triton X-100 Merck MilliporeSigma (Sigma-Aldrich) Catalog #93443-100ML

18 μL 10% NP40

420 μL Water

20 μL 20 mg/mL Protease Qiagen Catalog #19155

Add 15 μL **restriction enzyme release mix** to the sample.

Incubate on a thermal cycler with the following program: 50 °C for 45 min, 25 °C for 2 h.

Purify with 0.85X Ampure XP beads Beckman Catalog #A63881, elute in 10 μL water.

9 6. Y-shape ligation adapter annealing

Prepare the following 15 micromolar (μM) T1T2_TCTT ligation adapter.

1 μL T4 Polynucleotide Kinase Reaction Buffer New England Biolabs Catalog #B0201S

10.5 μL 100 micromolar (μM) T1_BtgZ1_ME

7.5 µL [IM] 100 micromolar (µM) **Rev_T2_BtgZ1**

31 µL Water.

Oligonucleotides sequences

T1_BtgZ1_ME: TCGTCGGCAGCGTC AGATGTGTAT

Rev_T2_BtgZ1: /5Phos/ TCTT ATACACATCT CCGAGCCCACGAGAC

Mix up the reaction, spin down, and perform the following reactions on a thermal cyclor.

65 °C 5 min

20 °C Hold, Ramp rate 0.1°C/s

10 7. Y-shape ligation adapter ligation

Prepare the following ligation mix (Total 50 µL).

10 µL Sample from Step 8

25 µL  StickTogether DNA Ligase Buffer **New England Biolabs Catalog #B0535S**

3 µL [IM] 15 micromolar (µM) T1T2_TCTT ligation adapter

3 µL  T7 DNA Ligase - 750,000 units **New England Biolabs Catalog #M0318L**

9 µL Water

Incubate at 25 °C for 30 min.

Purify with 0.8X  Ampure XP beads **Beckman Catalog #A63881** , elute in 51 µL water.

11 8. Quantification of library complexity

Dilute 1 µL of ligation product with 9 µL water. Use the 1:10 diluted product to quantify library complexity using qPCR. Any Illumina Nextera libraries with known concentration can be used as a standard. An example qPCR procedure is shown below.

5 µL  iTaq Universal SYBR Green Supermix **Bio-Rad Laboratories Catalog #172-5112**

0.25 µL [IM] 10 micromolar (µM) **T1ME**

0.25 µL [IM] 10 micromolar (µM) **T2ME**

3.5 µL water

1 µL template

Oligonucleotides sequences

T1ME: TCGTCGGCAGCGTC AGATGTGTATAAGAGACAG

T2ME: GTCTCGTGGGCTCGG AGATGTGTATAAGAGACAG

94 °C 2 min

30 cycles of

94 °C 20s

68 °C 20s

72 °C 1 min, fluorescence scanning

Melting curve

12 9. Library amplification

We recommend to amplify 14 cycles for a library with 20 million DNA fragments.

25 µL NEBNext UltraII Q5 Master Mix **New England Biolabs Catalog #M0544X**

2.5 µL 10 micromolar (µM) **Illumina Nextera N5XX index primer**

2.5 µL 10 micromolar (µM) **Illumina Nextera N7XX index primer**

20 µL ligation product and water

Oligonucleotides sequences

Illumina Nextera N5XX index primer: AATGATACGGCGACCACCGAGATCTACAC NNNNNNNN
TCGTCGGCAGCGTC

Illumina Nextera N7XX index primer: CAAGCAGAAGACGGCATACGAGAT NNNNNNNN
AGATGTGTATAAGAGACAG

PCR cycles:

98 °C 30s

14 cycles of

98 °C 10s

63 °C 30s

 72 °C 1 min

 72 °C 3 min

Purify with 0.8X  Ampure XP beads **Beckman Catalog #A63881**, elute in  20 µL water.