

Aug 18, 2022

# Growth Curve Stress Test (Instructor Protocol)

Brian Teague<sup>1</sup><sup>1</sup>University of Wisconsin - Stout

1 Works for me

Share

This protocol is published without a DOI.

Yeast ORFans CURE

Brian Teague  
University of Wisconsin - Stout

## ABSTRACT

This is the instructor protocol for



Growth Curve Stress Test  
by Brian Teague,  
University of Wisconsin - Stout

PREVIEW

RUN

## PROTOCOL CITATION

Brian Teague 2022. Growth Curve Stress Test (Instructor Protocol).  
**protocols.io**  
<https://protocols.io/view/growth-curve-stress-test-instructor-protocol-cffitjke>



## KEYWORDS

saccharomyces, growth, phenotype

## LICENSE

\_\_\_\_\_ This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## IMAGE ATTRIBUTION

Chalik1, CC BY 3.0 , via Wikimedia Commons

## CREATED

Aug 17, 2022

LAST MODIFIED

Aug 18, 2022

PROTOCOL INTEGER ID

68810

MATERIALS TEXT

### Equipment

- Shaking incubator (or preferably a roller drum in an incubator)
- Spectrophotometer for measuring OD600

### Materials

- [Yeast Nitrogen Base \(YNB\) Sunrise](#)
- [Science Catalog #1500-100](#) Step 1
- [CSM-Leu-Ura Powder Sunrise](#)
- [Science Catalog #1038-010](#) Step 1
- [Ammonium sulfate Sigma](#)
- [Aldrich Catalog #A4418](#) Step 1
- [Uracil Yeast Culture Grade Sunrise](#)
- [Science Catalog #1906-010](#) Step 1
- [L-Leucine Yeast Culture Grade Sunrise](#)
- [Science Catalog #1980-010](#) Step 1
- [α-D-Glucose Sigma](#)
- [Aldrich Catalog #158968](#) Step 3
- Test tubes for growing yeast cultures

SAFETY WARNINGS

Yeast and yeast media are not hazardous.

The yeast we are using are genetically modified. Make sure to inactivate them (by bleach or autoclaving, for example) before disposing of them, following your institution's guidelines.



Prepare 2x synthetic media 30m

- 1 Put **200 mL** of deionized water in a 250 ml bottle. Then add:

- [Yeast Nitrogen Base \(YNB\) Sunrise](#)
- **0.855 g** [Science Catalog #1500-100](#)
- [CSM-Leu-Ura Powder Sunrise](#)
- **0.335 g** [Science Catalog #1038-010](#)



 **Ammonium sulfate Sigma**

■  **2.5 g Aldrich Catalog #A4418**

■  **5 mL** 100X uracil solution (  **2 mg/mL**

 **Uracil Yeast Culture Grade Sunrise**

**Science Catalog #1906-010** )

■  **5 mL** 100X leucine solution (  **12 mg/mL**

 **L-Leucine Yeast Culture Grade Sunrise**

**Science Catalog #1980-010** )

2 Autoclave  **121 °C** for  **00:30:00**

30m

3 Cool, then add  **50 mL**  **40 Mass / % volume** gluco

 **α-D-Glucose Sigma**




**Aldrich Catalog #158968**

using aseptic technique.

#### Prepare the yeast cultures

4 The afternoon or evening before lab, prepare a test tube with 5 ml 2xSDC media for each yeast strain you will be testing.

5 Pick a colony from each strain's plate and resuspend it in the media.

6 Grow  **Overnight** at  **30 °C** , either on a roller drum (preferred) or on an orbital shaker at  **180 rpm** .

7 The next morning, several hours before class, measure the OD600 of each culture.

8 Prepare another test tube with 5 ml of 2xSDC media for each yeast strain. Subculture the overnights to an OD600 of 0.1. Place the tubes back in the incubator or shaker.

## 9 Optimally, grow the cultures for at least 3 hours before lab.

The experiments will be more reproducible if the cells always start the experiment in logarithmic growth phase!

### Instructor Tips & Common Student Errors

#### 10 Instructor Tips

- This assumes that, once you or your students have identified yeast colonies with successful knockouts, you have struck them out on plates
- Log-phase growth is the most reproducible growth mode for yeast. This is why the experiment should start in log phase and why doubling time should be computed in log phase.
- Growing the cells in double-strength media means that they can be diluted 1:1 and still be in "complete" media for the stress experiment.
- Yeast cells are heavy and settle rapidly toward the bottom of a test-tube. Make sure you vortex them before measuring OD600; vortex again before subculturing; and vortex again before putting them in an assay plate.
- Computing doubling time is not optional! Because we're looking at an exponential phenomenon, a small difference in starting concentration is rapidly amplified and can "look" on a graph like a strain is growing slower.
- I thought everyone graduated highschool with a basic ability to use a spreadsheet? I am often surprised by the number of students that have trouble with "make a graph of this data."

#### 11 Common student errors

- Not computing doubling time to compare growth, but just drawing conclusions from eyeballing a graph.
- Not vortexing the yeast cells before putting them in the assay plate.
- General data manipulation and interpretation difficulty