

•



Aug 23, 2020

© ELISA for quantification of human immunoglobulin D (IgD) in serum or plasma.

Angel A Justiz-Vaillant¹, Belkis Ferrer-Cosme²

¹University of the West Indies St. Augustine; ²"Saturnino Lora Torres' Provincial Teaching Clinical Surgical Hospital. Cuba

In Development dx.doi.org/10.17504/protocols.io.bj7wkrpe

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant University of the West Indies St. Augustine

ABSTRACT

IgD is a monomer with a molecular weight of 184 Kd. IgD is present in a meager amount in the serum (0.03 mg/mL) and has an unknown function against pathogens. It is regarded as a BCR. IgD may play an essential role in antigen-triggered lymphocyte differentiation. [1]

DOI

dx.doi.org/10.17504/protocols.io.bj7wkrpe

PROTOCOL CITATION

https://dx.doi.org/10.17504/protocols.io.bj7wkrpe

LICENSE

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 23, 2020

LAST MODIFIED

Aug 23, 2020

PROTOCOL INTEGER ID

40918

- 1 An anti-human IgD coating antibody is adsorbed onto the microwells by incubation overnight at 4°C with carbonatebicarbonate buffer.
- 2 Add 50 μl of human serum or plasma. Human IgD present in the serum or plasma binds to antibodies adsorbed into the microwells.
- 3 The microplate is blocked with 3% non-fat milk-PBS buffer and later wash to remove unbound proteins.
- 4 Fifty (50) μl of biotin-conjugated anti-IgD antibody is added. The optimal dilution must be investigated.

protocols.io
1
08/23/2020

Citation: Angel A Justiz-Vaillant, Belkis Ferrer-Cosme (08/23/2020). ELISA for quantification of human immunoglobulin D (IgD) in serum or plasma.. https://dx.doi.org/10.17504/protocols.io.bj7wkrpe

5	The microplate is rewashed with PBS-Tween 20 buffer, pH 7.4.
6	One hundred μ I of streptavidin-HRP conjugate is added and it binds to the biotin-conjugated anti-IgD antibody. The optimal dilution of this conjugate must be investigated.
7	The plate is washed following incubation to remove the unbound Streptavidin-HRP.
8	Add 100 μl of 3,3',5,5'- tetramethylbenzidine (TMB; Sigma-Aldrich) into each well.
9	Incubate the microwells in the dark for 20 min.
10	A colored product is formed in proportion to the quantity of human IgD present in the sample or standard.
11	The reaction is terminated by addition of 100 $\mu l3MH2SO4$ and the absorbance is measured at 450 nm.
12	A standard curve is made from 7 human IgD standard dilutions and the human IgD sample concentration is determined.
13	For better results place the microplate on a microplate shaker in every incubation.