



Sep 06, 2022


🌐 Harvesting and irradiation of mouse embryonic fibroblasts (MEFs) for hPSC cultures

 In 1 collection

Hanqin Li¹, Oriol Busquets², Steven Poser², Dirk Hockemeyer¹, Frank Soldner²

¹University of California, Berkeley; ²Albert Einstein College of Medicine

1 *Works for me*

 Share

dx.doi.org/10.17504/protocols.io.b4n7qvhn



Devin E Snyder

ABSTRACT

This protocol describes the process of harvesting and irradiating mouse embryonic fibroblasts (MEFs) to use as feeder cells for human pluripotent stem cell (hPSC) culture.

General notes

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.
2. MEFs were obtained as described in *Manipulating the Mouse Embryo: A Laboratory Manual*, Third Edition (ISBN: 0879695919)

Andras Nagy, Marina Gertsenstein, Kristina Vintersten, & Richard Behringer. *Manipulating the Mouse Embryo: A Laboratory Manual*, 3rd ed.. Cold Spring Harbor Laboratory Press.

DOI

dx.doi.org/10.17504/protocols.io.b4n7qvhn

PROTOCOL CITATION

Hanqin Li, Oriol Busquets, Steven Poser, Dirk Hockemeyer, Frank Soldner 2022.
Harvesting and irradiation of mouse embryonic fibroblasts (MEFs) for hPSC
cultures. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4n7qvhn>



FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's
Grant ID: ASAP-000486

COLLECTIONS




**Maintenance and inactivation of mouse embryonic fibroblasts (MEFs) as
feeder cells for human pluripotent stem cell culture**

KEYWORDS

ASAPCRN

LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

57791


PARENT PROTOCOLS

Part of collection

[Maintenance and inactivation of mouse embryonic fibroblasts \(MEFs\) as feeder cells for human pluripotent stem cell culture](#)

MATERIALS TEXT

Item	Vendor	Catalog #
DMEM	Corning	10-013-CV
FB Essence	Avantor	10803-034
FBS	Gibco	10437028
200mM L-Glutamine	Sigma	G8540
Penicillin & Streptomycin	Gibco	15140-122
MEM Non-Essential Amino Acids	Gibco	11140-050
Gelatin powder	Sigma	G2625
0.25% Trypsin/EDTA (Trypsin)	Gibco	25200-056
DPBS w/o Ca & Mg (DPBS)	Corning	MT21031CV
50ml centrifuge tubes	Corning	1495949A
10cm petri dish	Fisher	08757100D
10ml serological pipet	Corning	7200574
15cm tissue culture dish	Corning	0877224

- 1 Wash the plates once with 10 ml DPBS for each 15-cm plate.
- 2 Add 3 ml Trypsin and incubate for in 37°C; 5% CO₂ for  00:10:00 10m
- 3 Add 9 ml MEF medium to neutralize the Trypsin. Collect cell suspension in a 50 ml conical tube.

3.1 MEF medium

A	B
DMEM	435 ml
FB Essence/FBS*	75 ml
200mM L-Glutamine	5 ml
Penicillin & Streptomycin (100x)	5 ml
MEM Non-Essential Amino Acids	5 ml

*We have successfully used either FB Essence or FBS and have not observed an obvious difference. Final volume: 500ml

- 4 Rinse the plate with 9 ml of new MEF medium to collect the remaining cells.
- 5 If there are cell clumps that not fully dissociated, let the clumps settle by gravity for 3 min then transfer the single cell suspension to a new conical tube. Cells in clumps are usually not uniformly irradiated and that may lead to proliferating feeders.
- 6 Parafilm the conical tubes and bring cells to irradiation.
- 7 Irradiate the cells in conical tubes at 3500 cGy on Precision X-RAD 320 irradiator.