



AUG 09, 2023

Induction of starvation stress

Louise Uoselis¹

¹Lazarou Lab, WEHI



Louise Uoselis
WEHI

ABSTRACT

Induction of starvation stress using EBSS in HeLa cells

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.5jyl8pq7g2w/v1

Protocol Citation: Louise Uoselis 2023. Induction of starvation stress.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.5jyl8pq7g2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 09, 2023



Last Modified: Aug 09, 2023

Day 1

- 1 Seed cells, aiming for a confluency of 90 – 100% at the time of treatment the next day.

Day 2

9h

- 2 Feed cells for  01:00:00 in an appropriate volume of standard growth media. 1h
- 3 Remove the media, and wash each well 3 times with an ample amount of room temperature PBS (eg. Wash with 3 mL of PBS per well of a 6 well plate).
- 4 After the third wash, add an appropriate volume of Earle's Balanced Salt Solution (EBSS) to each well/plate.
- 5 After  08:00:00 of incubation in EBSS, harvest all samples for analysis. 8h