

JUN 20, 2023

OPEN ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.5qpvo3n9bv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvo3n9bv4o/v1)

**Protocol Citation:** Laura.SabioRodriguez 2023. Lipid (Oil Red O) Staining .  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpvo3n9bv4o/v1>

**MANUSCRIPT CITATION:**  
 Lipid (Oil Red O) Staining Kit.  
 SIGMA-ALDRICH, Catalog  
 Number: MAK194

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
 We use this protocol and it's working

**Created:** Jun 20, 2023

**Last Modified:** Jun 20, 2023


**PROTOCOL integer ID:**  
 83716

## 🌐 Lipid (Oil Red O) Staining

Laura.SabioRodriguez<sup>1</sup>

<sup>1</sup>CeMi

Centre for the Cellular Microenvironment [CeMi]



Laura.SabioRodriguez

### ABSTRACT

Lipid droplets (LDs) are dynamic, ubiquitously present lipid-storage organelles, predominantly present in the adipocytes. Triglycerides, neutral lipids, and cholesterol esters stored in LDs are the largest sources of energy. The presence of excess LDs in adipocytes results in obesity and obesity-linked pathologies such as dyslipidemia and diabetes type.

### ATTACHMENTS

[Lipid Oil Red O Staining kit.pdf](#)

### SAFETY WARNINGS



#### Hazard Classifications

Acute Tox. 3 Inhalation - Acute Tox. 4 Oral - Carc. 1B - Eye Irrit. 2 - Muta.  
 2 - Skin Irrit. 2 - Skin Sens. 1 - STOT SE 3

**Keywords:** Lipid droplets, Staining, Adipocytes, Oil Red O

## BEFORE START INSTRUCTIONS

Kit: **Lipid (Oil Red O) Staining Kit. SIGMA-ALDRICH, Catalog Number: MAK194**

Briefly centrifuge vials before opening. Use ultrapure water for the preparation of reagents.


Reconstitute the Oil Red O Stock solution with  20 mL of



Isopropanol (IPA) 100% Contributed by users

. Mix well and leave undisturbed

for  00:20:00 to make the Stock Solution, which is stable for 1 year.

The Oil Red O Working Solution is obtained by adding 3 parts of Oil Red O Stock Solution to 2 parts of water. Mix well and leave undisturbed for  00:10:00 and filter through Whatman No. 1 filter paper. The Working Solution is stable for 2h and must be prepared 15 min before use.


## Procedure

30m

### 1 All components:

Use  100  $\mu$ L / well for a 96 well plate


Use  2 mL / well for a 6 well plate

Use  6 mL / well for a 100 mm culture dish


### 2 Fixation

30m


Remove the medium and gently wash the cells twice with

 1X PBS (Phosphate-buffered saline)

. Add

 10 % formalin Contributed by users


to the

cells and incubate for  00:30:00 to 1 h

*Note: Do not add formalin directly onto the cells. Pipette onto the wall and mix gently rotating.*

### 3

Discard the formalin and wash the cells twice using

 Distilled Water Contributed by users

5m


Add 60%  Isopropanol Contributed by users to the cells and incubate for  00:05:00 .

### 4


Discard 60% isopropanol and cover the cells evenly with **Oil Red O Working Solution**. Rotate the

10m

plate or dish, and incubate for

 00:10:00 to 20 min

5 Discard the Oil Red O Solution and wash the cells 2-5 times with


 Distilled Water Contributed by users


until no excess stain is seen.

6 Add  Hematoxylin Contributed by users to the cells and incubate for  00:01:00 .

1m

Discard hematoxylin and wash the cells 2-5 times with

 Distilled Water Contributed by users

7 Cover the cells with  Distilled Water Contributed by users and view under microscope. Lipid

droplets appear red and nuclei appear blue.

*Note: Keep the cells covered with water to avoid drying.*