



Version 1

Dec 21, 2020

# qPCR: Bacterial SSU rRNA 338F-516P-805R V.1

Roey Angel<sup>1</sup>, Eva Petrova<sup>1</sup>, Ana Lara<sup>1</sup><sup>1</sup>Soil and Water Research Infrastructure

1 Works for me dx.doi.org/10.17504/protocols.io.bqwymxfw

SoWa RI Anaerobic and Molecular Microbiology (public)  
Tech. support email: [eva.petrova@bc.cas.cz](mailto:eva.petrova@bc.cas.cz)Roey Angel  
Soil and Water Research Infrastructure

## ABSTRACT

Universal 16S rRNA probe-based-qPCR assay for bacteria.  
The primers and probe are taken from [Yu et al. \(2005\)](#).

Yu Y, Lee C, Kim J, Hwang S (2005). Group-specific primer and probe sets to detect methanogenic communities using quantitative real-time polymerase chain reaction. Biotechnology and bioengineering. <http://dx.doi.org/10.1002/bit.20347>

## THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Yu, Y., Lee, C., Kim, J., and Hwang, S. (2005). Group-specific primer and probe sets to detect methanogenic communities using quantitative real-time polymerase chain reaction. Biotechnol Bioeng 89, 670–679. doi:10.1002/bit.20347.

## DOI

[dx.doi.org/10.17504/protocols.io.bqwymxfw](https://dx.doi.org/10.17504/protocols.io.bqwymxfw)

## PROTOCOL CITATION

Roey Angel, Eva Petrova, Ana Lara 2020. qPCR: Bacterial SSU rRNA 338F-516P-805R. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bqwymxfw>  
Version created by Roey Angel

## MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Yu, Y., Lee, C., Kim, J., and Hwang, S. (2005). Group-specific primer and probe sets to detect methanogenic communities using quantitative real-time polymerase chain reaction. Biotechnol Bioeng 89, 670–679. doi:10.1002/bit.20347.

## KEYWORDS

qPCR, dual-labelled probe, 16S rRNA gene, bacteria

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Dec 20, 2020

LAST MODIFIED

Dec 21, 2020

PROTOCOL INTEGER ID

45752

MATERIALS TEXT

MATERIALS

 [iQ™ SYBR® Green Supermix](#) **BioRad**

**Sciences Catalog #1708880**

 [TaqMan™ Fast Advanced Master Mix](#) **Thermo Fisher**

**Scientific Catalog #4444556**

Step 2

ABSTRACT

Universal 16S rRNA probe-based-qPCR assay for bacteria.

The primers and probe are taken from [Yu et al. \(2005\)](#).

Yu Y, Lee C, Kim J, Hwang S (2005). Group-specific primer and probe sets to detect methanogenic communities using quantitative real-time polymerase chain reaction. Biotechnology and bioengineering. <http://dx.doi.org/10.1002/bit.20347>

## Primers and probe

1

Name	Type	Sequence	Target region <sup>1</sup>
BAC338F	Forward	ACT CCT ACG GGA GGC AG	338-354
BAC516P <sup>2</sup>	Probe	TGC CAG CAG CCG CGG TAA TA	516-536
BAC805R	Reverse	GAC TAC CAG GGT ATC TAA TC	785-805

1. Relative to *E. coli* SSU rRNA gene

2. The probe must be dual-labelled either with 5'-6-FAM, 3'-BHQ1 or any other valid combination

## qPCR mixture

2

A	B	C	D
Reagent	Final concentration	1 tube (20 µl)	plate (20 µl x 100)
PCR H <sub>2</sub> O		2.2	220
2x TaqMan Fast Advanced Master mix	1x	10	1000
BSA (20 µg/µl)	0.4mg/ml	0.4	40
338F	0.5 µM	1.0	100
<b>805R</b>	0.5 µM	1.0	100
<b>516P</b>	0.2 µM	0.4	40
Template		5	5 x 100







 **TaqMan™ Fast Advanced Master Mix** **Thermo Fisher**

**Scientific Catalog #4444556**

Thermocycler programme 1h 30m

1h 30m

3

1.  **95 °C** for  **00:05:00**
2. x 40 {
  - 2.1  **95 °C** for  **00:00:30**
  - 2.2  **62 °C** for  **00:01:00** take snapshot
- }