



Mar 12, 2021

# Making MES Buffers for Protein EDAC Particle Coupling

Kenneth Schackart<sup>1</sup>, Kattika Kaarj<sup>1</sup><sup>1</sup>University of Arizona

1 Works for me dx.doi.org/10.17504/protocols.io.bmndk5a6

Yoon Lab

 Kenneth Schackart  
University of Arizona

SUBMIT TO PLOS ONE

## ABSTRACT

This protocol describes how to make MES buffers for covalently coupling proteins to carboxylated polystyrene particles via EDAC. Adapted from Bangs Laboratories Inc. TechNote 205 Covalent Coupling.

Three buffers are made from 2-(*N*-morpholino)ethanesulfonic (MES) acid with varying pH. The three buffers are referred to as :

- Activation buffer (acidic for activating the carboxyl group on particles for EDAC reaction)
- Coupling buffer (alkaline for coupling protein to o-acylisourea on particle surface created by EDAC reaction)
- Storage buffer (for storing particles while preventing degradation and self aggregation)

## ATTACHMENTS

[bangs\\_laboratories\\_technote205\\_covalent\\_coupling.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.bmndk5a6](https://dx.doi.org/10.17504/protocols.io.bmndk5a6)

## PROTOCOL CITATION

Kenneth Schackart, Kattika Kaarj 2021. Making MES Buffers for Protein EDAC Particle Coupling .  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bmndk5a6>

## KEYWORDS

MES buffer, EDAC, particle conjugation, protein conjugation, carboxylated microparticles, carboxyl modified microparticles, polystyrene particles, carbodiimide

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Sep 23, 2020

## LAST MODIFIED

Mar 12, 2021

## PROTOCOL INTEGER ID

42405

## GUIDELINES

When adjusting pH, be sure to mix well and check pH continuously or frequently (if using pH strips).

## MATERIALS TEXT

### MATERIALS

☒ [MES, free acid, monohydrate](#) **Bio Basic**

**Inc. Catalog #MB0341.SIZE.100g**

☒ [BSA Sigma](#)

**Aldrich Catalog ##A8806**

☒ [Hydrochloric Acid](#) **Fisher**

**Scientific Catalog #A144S**

☒ [Sodium hydroxide](#) **Sigma-**

**aldrich Catalog #306576**

### Reagents:

- MES (2-(*N*-morpholino)ethanesulfonic) acid
- Deionized water (DW)
- Sodium hydroxide (NaOH) solution, **[M]1 Molarity (M)**
- Hydrochloric acid (HCl), **[M]1 Molarity (M)**
- Bovine serum albumin (BSA)

### Materials:

- pH meter or strips

## BEFORE STARTING

Ensure that pH meter is calibrated.

### Make MES Buffer Bulk

- 1 Combine:
  - **12.8 g** MES free acid.
  - **600 mL** DW.
  - Swirl until MES is dissolved.
- 2 Adjust to ~6 pH. If pH is low add NaOH, if pH is high add HCl.

Typically the solution is about **pH3**. In this case, adding about **16 mL** of NaOH **[M]1 Molarity (M)** should get you close. If the desired activation buffer pH is lower than **pH6**, adjust at this step to that pH.

### Separate and Adjust

- 3 Separate the bulk buffer into the three buffers by transferring one third of the bulk to each of three glass bottles (~**200 mL** each).
- 4 Adjust pH to the desired level for each buffer. NaOH is used to raise pH, HCl is used to lower pH. Recommended ranges from Tech Note 205 are shown below.
  - Activation buffer: **pH4.5** - **pH7.5**
  - Coupling buffer: **pH7.2** - **pH8.5**

- Storage buffer: pH7.0

To reach ~ pH8 Coupling buffer, about 13 mL 1 Molarity (M) NaOH should be added.

To reach ~ pH7 Storage buffer, about 10 mL 1 Molarity (M) NaOH should be added.

#### Finish Storage Buffer

- 5 Add bovine serum albumin (BSA) to make 0.1 % w/v to Storage buffer (not Activation or Coupling buffer).

- For 200 mL , add 200 mg BSA.

Swirl to dissolve BSA.

BSA coats the particle surface, preventing nonspecific interaction (e.g. hydrophobic), also reducing self aggregation. The concentration of BSA can be varied from 0.05 % w/v - 0.1 % w/v . Additionally, surfactants may be used (e.g. Tween-80). Other blocking agents may be used in place of BSA.

#### Store Buffers

- 6 Store all buffers covered at 20 °C .